# SciComm-Bulletin

ZÜRICH, SCHWEIZ – 21.08.2009

## TiUnite®

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## Referenzen

- Allegrini S, Jr., Allegrini MR, Yoshimoto M, Konig B, Jr., Mai R, Fanghanel J, Gedrange T. Soft tissue integration in the neck area of titanium implants – an animal trial. *J Physiol Pharmacol* 2008;59 Suppl 5:117-32.
- Bahat O. Technique for placement of oxidized titanium Implants in compromised maxillary bone: prospective study of 290 implants in 126 consecutive patients followed for a minimum of 3 Years after loading. *Int J Oral Maxillofac Implants* 2009; 24: 325-334
- Calandriello R, Tomatis M. Immediate occlusal loading of single lower molars: five-year clinical follow-up from a prospective multi-center study using TiUnite implants. *Europero* 6, 2009, Stockholm, abstract #403.
- Friberg B, Jemt T. Rehabilitation of edentulous mandibles by means of five TiUnite implants after one-stage surgery: a 1-year retrospective study of 90 patients. *Clin Implant Dent Rel Res* 2008;10(1):47-54
- Friberg B, Jemt T. Clinical experience of TiUnite implants: a 5-year cross-sectional retrospective follow-up study. *Clin Implant Dent Relat Res* 2009, accepted for publication.
- Glauser R, Portmann M, Ruhstaller P, Lundgren AK, Hämmerle C, Gottlow J. Stability measurements of immediately loaded machined and oxidized implants in the posterior maxilla: a comparative clinical study using resonance frequency analysis. *Appl Osseointegration Res* 2001;2:27-9.
- Glauser R, Zembic A, Ruhstaller P, Windisch S. Five-year results of implants with an oxidized surface placed predominantly in soft quality bone and subjected to immediate occlusal loading. *J Prosthet Dent* 2007 Jun;97(6 Suppl):S59-68. Erratum in *J Prosthet Dent* 2008;99(3):167. AO 2009, San Diego, oral presentation with seven-year follow-up data.
- Hall J, Miranda-Burgos P, Sennerby L. Stimulation of directed bone growth at oxidized titanium implants by macroscopic grooves: an in vivo study. *Clin Implant Dent Relat Res* 2005;7 Suppl 1:76-82.
- Ivanoff CJ, Widmark G, Johansson C, Wennerberg A. Histologic evaluation of bone response to oxidized and turned titanium micro-implants in human jawbone. *Int J Oral Maxillofac Implants* 2003;18(3):341-8.
- Kielbassa AM, Martinez-de Fuentes R, Goldstein M, Arnhart C, Barlattani A, Jackowski J, Knauf M, Lorenzoni M, Matorana C, Mericske-Stern R, Rompen E, Sanz M. Randomized controlled trial comparing a variable-thread novel tapered and a standard tapered implant: interim one-year results. *J Prosthet Dent* 2009;101:293-305.
- Mura P. Immediate and early loading of Replace Select Tapered TiUnite implants placed in post-extractive sockets and healed sites. Retrospective analysis of the 5-year clinical outcome. EAO 2009, Monaco, abstract accepted.
- Nickenig HJ, Wichmann M, Schlegel KA, Nkenke E, Eitner S. Radiographic evaluation of marginal bone levels adjacent to parallel-screw cylinder machined-neck implants and rough-surfaced micro-threaded implants using digitized panoramic radiographs. *Clin Oral Implants Res* 2009, epub March 2009.
- Poulos N, Rodriguez N, Lee J, Rueggeberg F, Schüpbach P, Hall J, Susin C, Wikesjö U ME. Evaluation of a novel calcium phosphate coated titanium porous oxide implant surface: a study in rabbits, *Clin Implant Dent Relat Res* 2009, submitted for publication.
- Schüpbach P, Glauser R. The defense architecture of the human peri-implant mucosa: a histological study. *J Prosthet Dent* 2007;97:15-25.
- Xiropaidis AV, Qahash M, Lim WH, Shanaman RH, Rohrer MD, Wikesjö U ME, Hall J. Bone-implant contact at calcium phosphate-coated and porous titanium oxide (TiUnite)-modified oral implants. *Clin Oral Implants Res* 2005;16:532-9
- Zechner W, Tangl S, Fürst G, Tepper G, Thams U, Mailath G, Watzek G. Osseous healing characteristics of three different implant types. *Clin Oral Implants Res* 2003;14(2):150-7.