

Clinical case

NobelActive® – Dr. Eric Rompen

January, 2012





Dr. Eric Rompen

Professor and Department
Head of Periodontology/ Dental
Surgery, University of Liège

Belgium

„An incredibly effective implant that provides excellent stability even in compromised sites; because of its unique features to preserve the biology of marginal hard and soft tissues, NobelActive has become my No.1 implant in the esthetic area.“



Anterior restoration

Male patient
20 years old

Patient congenitally missing
two upper lateral incisors,
labial profile with concavities
and lack of bone volume.

Implants:

NobelActive NP 3.5x13mm

FDI POS:

Teeth 12, 22


Initial clinical situation



Diagnosis

CT scan confirmed sufficient bone volume (labial/palatal width) for placement of narrow diameter implants due to slight subsidence of facial bone table.

Implant placement



Minimally invasive surgical approach with a semi-lunar primary incision of the ridge and two very short release incisions. Placement of two NobelActive NP 3.5x13mm implants.

Soft tissue contour development



Intra-oral resizing of Procera Esthetic Abutments Zirconia using profuse irrigation.

Provisional restoration



Immediate function with provisional crowns. Picture shows tissue response to provisional restoration six weeks after implant placement.




Precise adaptation of the provisional restoration

Final restoration



Placement of NobelProcera Crowns four months after surgery. The soft tissue and radiographic results after 24 months show excellent soft tissue contours and stable bone levels.

24 months follow-up



X-rays 24 months after implantation showing stable bone levels.