
GoldAdapt™ Engaging GoldAdapt™ Non-Engaging

INSTRUCTIONS FOR USE


73506 C West 1201 Printed in Sweden. Date of issue 2011-06-28.
All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare.
Product images in this folder are not necessarily to scale.

GoldAdapt™ Engaging

GoldAdapt™ Non-Engaging

1

GoldAdapt™ screw chart

	Abutment screw	Lab screw
Internal tri-channel connection NP	36818	31170
Internal tri-channel connection RP, WP, 6.0	29475	29293
External hex connection NP	29282	31168
External hex connection RP	29283	29290
External hex connection WP	29284	31169
Internal conical connection 3.0	36776	36805
Internal conical connection NP	36917	31168
Internal conical connection RP	36918	29290

ENGLISH

Important: Please read.

GoldAdapt™ Engaging, GoldAdapt™ Non-Engaging
Internal conical connection for NobelActive™ and NobelReplace® CC.
Internal tri-channel connection for the following implant systems: NobelReplace®, Replace Select™ and NobelSpeedy™ Replace.
External hex connection for Brånemark System® and NobelSpeedy™ Groovy.

Disclaimer of liability: These prosthetic components are part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. Nobel Biocare assumes no liability whatsoever for damage arising thereof.

Description: A premanufactured endosseous dental implant abutment directly connected to the endosseous dental implant intended for use as an aid in prosthetic rehabilitation.

Indications: GoldAdapt™ Engaging are intended for single tooth screw-retained and cement-retained crowns as well as for multiple unit cement-retained implant bridges.

GoldAdapt™ Engaging Conical Connection 3.0 are intended to be used together with NobelActive™ 3.0 implant to replace a single lateral incisor in the maxilla or central or lateral incisor in the mandible. The screw-retained solution may be used when the screw access hole is located through the occlusal surface or through the cingulum. GoldAdapt™ Non-Engaging are intended for screw-retained multiple teeth fixed prostheses. This screw retained solution is mainly indicated when the screw access hole are located through the occlusal surface or through the cingulum as well as for limited interocclusal space. Indicated for implants with less than 40° overall divergences to allow path of insertion.

Contraindications: It is contraindicated to use GoldAdapt™ Engaging Conical Connection 3.0 in other positions than for lateral incisors in the maxilla or central and/or lateral incisors in the mandible. GoldAdapt™ Engaging Conical Connection 3.0 is not to be used for multiple unit restorations.

Instruction for clinician: We strongly recommend that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge. For more info visit www.nobelbiocare.com. Always work with an experienced colleague the first time you employ a new treatment method. Nobel Biocare has a global network of mentors for this purpose.

Procedural Precautions: Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, avoiding excessive transverse loading forces, particularly in immediate loading cases. If the prosthesis metal substructure is made of gold alloy, this should have a high gold content. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Sterility: The GoldAdapt™ Abutments co-packed with abutment screw are delivered non-sterile. Steam sterilize abutments and the prosthetic kit, if applicable for 5 min at 135°C/274°F.
Note: Use of non-sterile abutments may lead to infection of tissues or infectious diseases.

Handling procedures

Clinical Procedure

- Place the impression coping implant level onto the implant and take an implant level impression (**A+B**).
 - Connect the healing abutment or temporary restoration.
- #### Laboratory Procedure
- Assemble the impression coping and implant replica and position into impression.
 - Fabricate a working model with removable gingival material.
 - Attach the GoldAdapt™ into implant replica and secure with lab screw.
 - Connect abutment and reduce the plastic sleeve to appropriate height (**C**) and wax-up a framework.
 - Fabricate the final abutment or framework using standard C&B techniques (**D**).
Note: Do not sandblast the seating surfaces.

Clinical procedure

- For single units:** Connect the customized abutment. It is recommended to verify the final abutment seating using radiographic imaging.
 - For multiple units:** Connect the customized abutments or implant bridge/bar. It is recommended to verify the final implant restoration seating using radiographic imaging.
 - Tighten the customized abutment(s) or implant bridge/bar abutment, except for Conical Connection 3.0, to **35 Ncm** using Unigrip™ Screwdriver and Manual Torque Wrench prosthetic (**E:1**). For Conical Connection 3.0 tighten abutment to **15 Ncm** using screwdriver and wrench described above (**E:2**).
Caution: For Conical Connection 3.0; Never exceed **15 Ncm** prosthetic tightening torque for the abutment screw. Over tightening of abutment may lead to a screw fracture.
 - Close screw access hole.
 - For Abutments: Cement final restoration if applicable.
Note: Do not use temporary cement when cementing ceramic crowns and bridges, due to increased risk of micro fractures.
- GoldAdapt™ screw chart (1).

Alloy Recommendations:

GoldAdapt™: Melting range: 1400–1490°C/2550–2720°F.
Coefficient of thermal expansion: 12 µm/m*°K.

Recommended casting alloys:

Conventional gold alloys: High gold content (min 75% Au + Pt metal) alloys, standard ISO 1562 type 4. Soldering in the range of 800–890°C/1472–1634°F.
Ceramic bonding alloys: High gold content (min 75% Au) alloys, standard ISO/DIS 9693, NiOM type A. Soldering in the range of 800–890°C/1472–1634°F.

For additional information on restorative and dental laboratory procedures please consult treatment guidelines available at www.nobelbiocare.com or request latest printed version from a Nobel Biocare representative.

Materials: GoldAdapt™: Gold alloy 60% Au, 19% Pt, 20% Pd, 1% Ir. Plastic: POM.
Abutment screws: Titanium alloy 90% Ti 6% Al 4%V.

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Please note that not all products may have been licenced in accordance with Canadian law.

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