NobelProcera™ Implant Bar Overdenture

INSTRUCTIONS FOR USE

73545 West 1109 Printed in Sweden. Date of issue 2011-07-21.
All rights reserved. Dolder® is a trademark of Cendres+Métaux SA. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale.

Implant Platform Screw Compatibility Table

Manufacturer	Implant Platform	Screw Art #
Nobel Biocare	Brånemark System® NP	29282
Nobel Biocare	Brånemark System® RP	29283
Nobel Biocare	Brånemark System® WP	29284
Nobel Biocare	NobelReplace® NP	36818
Nobel Biocare	NobelReplace® RP	29475
Nobel Biocare	NobelReplace® WP	29475
Nobel Biocare	NobelReplace® 6.0	29475
Nobel Biocare	NobelActive™, NobelReplace® CC, NP	36917
Nobel Biocare	NobelActive™, NobelReplace® CC, RP	36918
Nobel Biocare	Multi-unit NP/RP	29285
Nobel Biocare	Multi-unit WP	29286

NobelProcera™ and Procera® Abutment (clinical) Screw Tightening Torque

Brånemark System®, NobelReplace®, NobelActive™	35Ncm
--	-------

ENGLISH

Important: Please read.

NobelProcera™ Implant Bar Overdenture

Indications. The NobelProcera™ Implant Bar Overdenture is indicated for use as an overdenture bar that attaches to implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function.

Clinical and Laboratory input: The technician discusses the type of bar/special abutment and relevant retentive elements with the restoring dentist and patient.

Impression Implant or Abutment Level Open Tray Technique: A custom-made open impression tray is essential for an accurate impression and for access to the impression coping guide pins. A wax lid can cover the opening.

Handling procedures

Clinical Procedure

Obtaining impression

- Use a custom made rigid impression tray with occlusal screw access holes following an open tray impression protocol.
- Syringe an elastomeric impression material around the impression copings intraorally and take an impression.
- After complete setting of the impression material unscrew the impression coping screws and remove the impression tray.
- In addition, take preliminary interocclusal records.
- Inspect the impression for discrepancies and send to the laboratory.

Note: For NobelActive™ or NobelReplace® Conical Connection, ensure that Impression Coping Bridge Open Tray Conical Connection is used.

Laboratory Procedure

Fabricating definitive cast, registration records.

- In the laboratory attach implant or abutment analogues to the impression copings and pour the impression with low expansion dental stone.
- Use an acrylic resin base (optional: integrate provisional abutment cylinders) and fabricate a wax-rim for interocclusal records to properly orient the models in the articulator.

Clinical Procedure

Intraoral jaw relation records

- Use the wax-rims and a facebow for jaw relation records.

Laboratory Procedure

Fabricating set-up

- Make a diagnostic tooth set up on the definitive cast.

Clinical Procedure

Intraoral try-in of set-up

- Try-in the diagnostic tooth set-up to verify functional and esthetic parameters.

Laboratory Procedure

Design bar/abutments

- Scan the definitive cast and the respective implant/abutment locations, and the wax set-up and virtually design the anticipated framework in the NobelProcera™ software.
- Send scan data file to Nobel Biocare production facility.
- Upon return check for precision of fit on the definitive cast.

Bar Design Parameters

The basic design parameters for each bar type are given in the table below.

Bar Type	Bar Shape	Max. Canti-lever	Max. Span	Min. Diagonal section	Min. Bar Link Height
Dolder® Resilient Regular	Fixed shape	20mm	40mm	N/A	N/A
Dolder® Resilient Small	Fixed shape	13,5mm	35mm	N/A	N/A
Dolder® Rigid Regular	Fixed shape	20mm	40mm	N/A	N/A
Dolder® Rigid Small	Fixed shape	13,5mm	35mm	N/A	N/A
Hader	Fixed shape	20mm	40mm	N/A	N/A
Milled	Free form	30mm	45mm	2,9mm	2,5mm
Montreal and Montreal Lingual Metallic	Free form	40mm	45mm	2,9mm	2,5mm
Round	Fixed shape	10mm	25mm	N/A	N/A
Paris	Free form	22mm	40mm	2,9mm	2,5mm
Wrap Around	Free form	25mm	50mm	N/A	2,5mm

In addition to the basic design parameters in the table above, the following general design parameters should be followed:

- When there is a cantilever with clip, distal cylinders should have a minimum 1mm wall thickness.
- When using beveled distal cylinders, the minimum thickness point between the inner screw head surface and the external cylinder surface should not be less than 0,25mm.
- When using micro and macro fixed shape extension links over 10mm, a 1,5mm radius reinforcement is required between the cylinder and link.
- When screwed attachments are used there should be a minimum 1mm wall thickness between the external thread shape and the outside buccal and lingual wall of the bar.

Clinical Procedure

Intraoral verification of framework fit

- Control precision of fit of delivered overdenture bar intraorally. A passive and precise fit is essential to the success of the implants. A forced fit can result in stress concentration on one or more implants, potentially leading to implant failure. Problematic screw loosening or bone loss may result from non-passive bar fit. If the bar does not fit passively, a new bar must be fabricated.

Laboratory Procedure

Fabricate overdenture

- Manufacture overdenture or modify an existing denture according to standard protocol.
- Stability of the restoration can further be reinforced by a non-precious metal framework integrated into the overdenture.
- Integrate the housings of the selected retentive elements in the overdenture.
- Remove the restoration from the definitive cast after complete polymerization and finish the tissue-bearing areas for maximum support of the surrounding tissues and optimal hygiene around the implants.

Clinical Procedure

Delivery

- The Implant Bar Overdenture should be cleaned and sterilized. Sterilize by using steam sterilization for 5 minutes at 135°C.
- Implant level – Place the overdenture bar and tighten to **35Ncm** using the Manual Torque Wrench Prosthetic and Unigrip™ Screwdriver Machine.
- Abutment level – Verify the tightness of the abutments to **35Ncm** – For angled abutments to **15Ncm**. Connect the restoration to the abutments with prosthetic screws. Tighten the prosthetic screws to **15Ncm** using the Manual Torque Wrench Prosthetic and Unigrip™ Screwdriver Machine.
- Only Nobel Biocare provided screws should be used with the overdenture bar. Refer to the tables to identify the proper screw to be used for each implant or abutment platform.
- Deliver the denture. Adjust the retention as necessary. Verify the occlusal relationship; confirm resiliency and possible hinge axis movement.

Distributed in North and South America by Nobel Biocare USA LLC, Yorba Linda CA

Manufacturer: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Please note that not all products may have been licenced in accordance with Canadian law.

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Distribué en Amérique du nord et du sud par Nobel Biocare USA LLC, Yorba Linda CA
Fabricant: Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Suède.