


# Easy Abutment™

## INSTRUCTIONS FOR USE

22679 A ML West 1001 Printed in Sweden. Date of issue 2009-02-03. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Product images in this folder are not necessarily to scale.

### G Easy Abutment™ screw chart

	0.5–1 mm Abutment	1.5–2 mm Abutment
Repl NP	29839	29840
Repl RP, WP, 6.0	29841	28786
Brånemark NP	30663	30664
Brånemark RP	30725	30726
Brånemark WP	30727	30728

## ENGLISH

**Important: Please read.**

### Easy Abutment™

Internal connection for the following implant systems: NobelReplace™, Replace® Select and NobelSpeedy™ Replace.

External connection for the following implant systems: Brånemark System™, NobelSpeedy™ Groovy.

**Disclaimer of liability:** These prosthetic components is part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and surgical application of this product are under the control of the clinician, they are his/her responsibility. We assume no liability whatsoever for damage arising thereof.

**Description:** A premanufactured endosseous dental implant abutment directly connected to the endosseous dental implant intended for use as an aid in prosthetic rehabilitation.

**Indications:** Indication range from single tooth to multiple teeth restorations using cement-retained supra-constructions.

**Instruction for clinician:** We strongly recommend that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge. For more info visit [www.nobelbiocare.com](http://www.nobelbiocare.com). Always work with an experienced colleague the first time you employ a new treatment method. Nobel Biocare has a global network of mentors for this purpose.

**Procedural Precautions:** Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment. It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments, adjusting occlusion to the opposing jaw, avoiding excessive transverse loading forces, particularly in immediate loading cases. If the prosthesis metal substructure is made of gold alloy, this should have a high gold content. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

**Sterility:** The abutments co-packed with abutment screw, impression coping and healing cap are delivered sterile and are for single use only prior to the labelled expiration date. If re-sterilization of abutments is required use steam sterilization at 135°C/274°F, for 5 minutes. For steam sterilization of Prosthetic kit, sterilize at 135°C/274°F, for 5 minutes.

### Handling procedures

Modifications of abutments could be performed using a carborundum disk and carbide bur. Use profuse irrigation if adjustments are made intra-orally.

**Note:** *Occlusal reduction of the Easy Abutment™ must not exceed 1 mm, if so the abutment screw will be compromised. Never reduce the margin of the abutment as this will compromise seating when using Easy Abutment™ impression coping.*

### Clinical Procedure

- Select appropriate abutment.
- Connect the abutment (**Illustration A**). It is recommended to verify the final abutment seating using radiographic imaging.
- Tighten the abutment to **35Ncm** using Screwdriver Unigrip™ and Manual Torque Wrench prosthetic (**B**).
- Seal the screw access hole.
- Take an impression of abutment. If abutment margins have not been modified the Impression Cap could be used (**C, D and E**). A "snap" indicates that the impression cap is in place. Inject the impression material slowly into the hole on top of the impression cap and take an impression.
- Fabricate a temporary restoration and cement on abutment.  
**Caution:** *Do not use Temporary Coping Easy Abutment™ with polyurethane cements. The cement will not cure.*

### Laboratory Procedures

- Produce a working model with removable gingival material.
- Fabricate a crown or bridge with Procera® or with conventional casting technique using the plastic copings as burn-out patterns. A ceramic cap is also available for a single crown restoration.
- Complete the restoration with ceramic.

### Clinical Procedures

- Remove temporary restoration if applicable.
- Use the Screwdriver Unigrip™ and Manual Torque Wrench prosthetic to verify the tightening of the abutment to **35Ncm**.
- Cement the final crown or bridge using conventional procedures.  
**Note:** *Do not use temporary cement when cementing ceramic crowns and bridges, due to increased risk of microfractures.*

Easy Abutment™ screw chart (**table G**).

For additional information on restorative and dental laboratory procedures please consult treatment guidelines available at [www.nobelbiocare.com](http://www.nobelbiocare.com) or request latest printed version from a Nobel Biocare representative.

**Materials:** Abutment and abutment screws: Titanium alloy 90% Ti 6% Al 4%V. Healing caps: Polysulphone. Impression coping Easy Abutment™: Copolyester. Temp/plastic coping: Polycarbonate. Ceramic cap in aluminum oxide (for Easy Abutment™).

**Manufacturer:** Nobel Biocare AB, Box 5190, SE-402 26 Göteborg, Sweden. Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. [www.nobelbiocare.com](http://www.nobelbiocare.com).

**Caution:** The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

**Argentina**  
Nobel Biocare Argentina  
Phone: +54 11 4825 9696  
Cust. support:  
toll free 0800 800 66235

**Brazil**  
Nobel Biocare Brazil  
Phone: +55 11 5102 7000  
Cust. support: 0800 169 996

**Canada**  
Nobel Biocare Canada  
Phone: +1 905 762 3500  
Cust. support:  
+1 800 939 9394

**China**  
Nobel Biocare China  
Phone: +86 21 5206 6655  
Cust. support:  
+86 21 5206 0974

**Colombia**  
Nobel Biocare Colombia  
Phone: +57 1 820 9410  
Cust. support:  
toll free +1 800 012 1255

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