Snappy™ Abutment

INSTRUCTIONS FOR USE

73502 F West 1312 Printed in Sweden. Date of issue 2011-06-28.
All rights reserved. Nobel Biocare, the Nobel Biocare logo and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare.
Product images in this folder are not necessarily to scale.

A


B

C

D

E

F

1

Snappy™ Abutment screw chart

	Abutment screw	Lab screw
Internal tri-channel connection NP	36818	31170
Internal tri-channel connection RP, WP, 6.0	29475	29293
External hex connection NP	29282	31168
External hex connection RP	29283	29290
External hex connection WP	29284	31169
Internal conical connection NP	36917	31168
Internal conical connection RP	36918	29290

ENGLISH

Important: Please read.

Snappy™ Abutment

Internal conical connection for NobelActive® and NobelReplace® CC.
Internal tri-channel connection for the following implant systems: NobelReplace®, Replace Select™ and NobelSpeedy® Replace.
External hex connection for Bränemark System® and NobelSpeedy® Groovy.

Disclaimer of liability: These prosthetic components are part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Nobel Biocare. Use of products made by third parties in conjunction with Nobel Biocare prosthetic components will void any warranty or other obligation, express or implied, of Nobel Biocare. The user of Nobel Biocare prosthetic components has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Nobel Biocare disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Nobel Biocare products. The clinician is also obliged to study the latest developments in regard to the Nobel Biocare products and their applications regularly. In cases of doubt, the clinician has to contact Nobel Biocare. Since the processing and surgical application of this product are under the control of the clinician, he is his/her responsibility. Nobel Biocare assumes no liability whatsoever for damage arising thereof.

Description: A premanufactured endosseous dental implant abutment directly connected to the endosseous dental implant intended for use as an aid in prosthetic rehabilitation.

Indications: Indication range from single tooth to multiple teeth restorations using cement-retained supra-constructs.

Instruction for clinician: We strongly recommend that clinicians, new as well as experienced implant users, always go through special training before undertaking a new treatment method. Nobel Biocare offers a wide range of courses for various levels of knowledge. For more info visit www.nobelbiocare.com. Always work with an experienced colleague the first time you employ a new treatment method. Nobel Biocare has a global network of mentors for this purpose.

Procedural Precautions: Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

It is especially important to achieve proper stress distribution through passive adaptation and fitting of the bridge to the implant abutments; adjusting occlusion to the opposing jaw; avoiding excessive transverse loading forces, particularly in immediate loading cases. If the prosthesis metal substructure is made of gold alloy, this should have a high gold content. Because of the small size of prosthetic components, care must be taken that they are not swallowed or aspirated by the patient.

Sterility: The abutments co-packed with abutment screw, impression coping, healing cap and temporary coping (eng) are delivered sterile and are for single use only prior to the labelled expiration date. If re-sterilization of abutments is required use steam sterilization at 135°C/274°F, for 5 minutes. For steam sterilization of Prosthetic kit, sterilize at 135°C/274°F, for 5 minutes.

Note: Use of non-sterile abutments may lead to infection of tissues or infectious diseases.

Handling procedures

Modifications of abutments could be performed using copious water irrigation: Use high-speed dental and a fine diamond drill.

Note: Occlusal reduction of the Snappy™ Abutment should not be performed when planning to use Snappy™ impression coping as retention may be compromised.

Clinical Procedure

- Select appropriate abutment and check occlusal clearance.
- Connect the abutment (A). It is recommended to verify the final abutment seating using radiographic imaging.
- Tighten the abutment to **35 Ncm** using Unigrip™ Screwdriver and Manual Torque Wrench prosthetic (B).
- Press the impression coping onto the abutment. A "snap" indicates that the impression coping is in place (C).
- Take an impression (D).
- Provisionalize using plastic temporary coping or a healing cap (E).

Caution: Do not use plastic temporary coping with polyurethane cements. The cement will not cure.

Laboratory Procedure

- Produce a working model with removable gingival material.
- Caution:** When fabricating the model, use only Snappy™ Abutment 4.0 Abutment Replica in corresponding 4.0 Impression Coping and Snappy™ Abutment 5.5 Abutment Replica in corresponding 5.5 Impression Coping. Verify correct fit before casting the model.
- Fabricate a crown or bridge with NobelProcera® or with conventional casting technique using the plastic copings as burn-out patterns.
- Complete the restoration with ceramic.

Clinical Procedure

- Remove temporary restoration if applicable.
- Use the Unigrip™ Screwdriver and Manual Torque Wrench prosthetic to verify the tightening of the abutment to **35 Ncm**.
- Cement the final restoration using conventional procedures after sealing of access hole (F).

Note: Do not use temporary cement when cementing ceramic crowns and bridges, due to increased risk of microfractures.

Snappy™ Abutment screw chart (1).

For additional information on restorative and dental laboratory procedures please consult treatment guidelines available at www.nobelbiocare.com or request latest printed version from a Nobel Biocare representative.

Materials: Abutment and abutment screws: Titanium alloy 90% Ti 6% Al 4%V. Healing caps: Polysulphone. Impression coping Snappy™ Abutment: Nylon. Temp/plastic coping: Polycarbonate.

Manufacturer: Nobel Biocare AB, Box 5190, 402 26
Västra Hamngatan 1, 411 17 Göteborg, Sweden.
Phone: +46 31 81 88 00. Fax: +46 31 16 31 52. www.nobelbiocare.com

Caution: The caution text "Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist" is shown on labels with "Rx Only".

Please note that not all products may have been licenced in accordance with Canadian law.

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Fabricante: Nobel Biocare AB, Box 5190, 402 26