
# NobelClinician™ Viewer

Installation guide

# For Windows®.

Please follow the instructions carefully.

## Installation of NobelClinician Viewer software

1. First of all verify whether you are working on a Windows® 32-bit or 64-bit operating system. Both versions have specific installers.
2. Check that your computer has access to the Internet.
3. Make sure that your computer can view .pdf files. (Double-click on a .pdf file to test.) Adobe® Reader® is an example of software typically used to read .pdf files.
4. Log on to the computer with a user account that has Administrator rights.
5. Download the install executable from <http://www.nobelbiocare.com/nobelclinicianviewer>. Run the downloaded install executable.
6. A welcome message to the NobelClinician Viewer install wizard appears. Click “Next”.
7. The NobelClinician End-User License Agreement is displayed. Read it carefully. You can also print this information. To accept the agreement, click on “I accept the terms in the license agreement”, and then click “Next”. Otherwise, click “Cancel” and the installation procedure will abort.
8. Choose the “destination folder” for the NobelClinician Viewer software. This is the folder where the software files will reside. Then, click “Next”.
9. Click on “Install” to start the actual installation.
10. Select the preferred language to read the instructions for use. Then, click “Finish”.
11. The software starts.
12. Enter registration data such as name, e-mail and address to register the software.
13. “NobelClinician Viewer” icon is now available on your desktop. Double-click on it to start the software.

NobelClinician Viewer can also be started by clicking the Windows® Start button. Browse the list of all programs that appears and click on Nobel Biocare. Then click on the NobelClinician Viewer icon to start the software. The instructions for use are also available here in a printable format.

At all times the License Agreement is available in the “About” dialog within the NobelClinician Viewer software.

It is advised to use the software while being connected to the Internet for receiving software updates.

## Technical Support

If you encounter any problems during the installation, please contact your local Customer Support. (See page 20).

Windows® is a registered trademark of Microsoft Corporation, Adobe® and Reader® of Adobe Systems Incorporated, Apple® and Mac® of Apple Incorporated in the United States and/or other countries.

# For Mac® OS X.

Please follow the instructions carefully.

## Installation of NobelClinician Viewer software

1. Check that your computer has access to the Internet.
2. Log on to your computer with a user account that has Administrator rights.
3. Download the install executable from <http://www.nobelbiocare.com/nobelclinicianviewer>. Run the downloaded install executable.
4. Double-click the Apple® Disk Image file (.dmg) to mount it. The mounted NobelClinician Viewer disk will appear on your desktop.
5. Double-click this file. An installation dialog pops up, showing different options:
  - Install NobelClinician Viewer
  - Instructions for Use (IFU)
  - License Agreement
6. To install the NobelClinician Viewer software double-click the first option.
7. Enter the username and password of your operating system user account, if requested.
8. Click “Continue” to go to the End-User License Agreement. Read it carefully. You can also print this information. To accept the agreement, click on “I accept the terms in the license agreement”, and then click “Next”. Otherwise, click “Cancel” and the installation procedure will abort.
9. Click “Continue”. A warning is displayed, asking for agreement with the End-user License Agreement.
10. Click “Agree” to continue with the installation. When clicking “Disagree” the installation procedure is aborted.
11. Next you can keep the default location for the installation or you can change it to the location of your choice. This is the disk where the software files will reside. Then click “Continue”.
12. Click “Install” to start the actual installation.
13. When the installation is finished, the software is started and a registration wizard is shown.
14. Fill out the details and click “Register”.
15. Now a NobelClinician Viewer icon is available in Applications. Click on it to start the software.

At all times the License Agreement is available in the “About” dialog within the NobelClinician Viewer software.

It is advised to use the software while being connected to the Internet for receiving software updates etc.

## Technical Support

If you encounter any problems during the installation, please contact your local Customer Support. (See page 20).

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## Installation der NobelClinician Viewer Software

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