

# Regenerative solutions How-to reference guide



## Step-by-step + clinical case

# Horizontal ridge augmentation



Treatment protocol  
as performed by  
**Naheed Mohamed**  
DMD, MSD, Dip. Perio,  
Dip. ABP, FRCD(c)

### Product range

**Grafting material:** creos™ allo.gain, creos™ xenogain

**Membranes:** creos™ xenoprotect

**Sutures:** RESORBA® GLYCOLON™


**Fixation systems:** Pro-fix™ Tenting kit and screws, Master-Pin-Control

## Step-by-step protocol – Dr. Mohamed

**Initial clinical situation:** Atrophic alveolar crests with horizontal bone loss and vertical dimension unchanged – diagnosed with CBCT scan



**1**  
Create an initial incision using a new 15C blade on the horizontal and two vertical incisions 1-2 teeth away from augmentation area.




**2**  
Open the flap to reveal the bony ridge. Assess the amount of bone width that needs to be created for an esthetic outcome.



**3**  
From the planning, place NobelActive® TiUltra implants if there is enough bone to place the implants in the ideal 3D position.



**4**  
Decorticate the bone to stimulate blood flow into the site.




**5**  
Place a resorbable collagen membrane such as creos™ xenoprotect or Cytoplast™ RTM collagen on the lingual side and fixate with pins until the membrane is secure.



**6**  
Mix autogenous bone and mineralized allogenic bone or xenogenic bone and hydrate with saline. Once the bone is ready, place around the defect site and fixate the membrane buccally with pins until it is secure.



**7**  
Using a new 15C blade, release the periosteal fibers to prepare for tension-free primary closure.



**8**  
Suture using Cytoplast™ PTFE sutures on the horizontal and RESORBA® GLYCOLON™ on the verticals. The PTFE will remain closed to allow for healing. RESORBA® GLYCOLON™ sutures will absorb, leaving the verticals with an esthetic result.

# Horizontal ridge augmentation – continued

## Clinical case – Dr. Mohamed



### Initial clinical situation: Perio compromised centrals and laterals

Soft tissue management comes first. Tooth with root showing. Patient presented with periodontally compromised centrals and right lateral.


#### Phase 1

Removed the teeth and let the soft tissue heal. After the soft tissue has healed, move into the surgery to evaluate the vertical and horizontal bone deficiencies.


#### Phase 2



1 Tissue healed after phase 1.



2 Two vertical incisions to open the defect. Deficiency observed.



3 Use Pro-fix™ tenting screw to plan the height of the bone that you want to fill. Overfill with bone substitutes by about 15–20% more than the desired ridge volume.



4 Use a creos™ xenoprotect collagen membrane or Cytoplast™ RTM membrane and pin on the palatal side.



5 Fill the space with creos™ xenogain and creos™ allo.gain bone substitutes.



6 Stretch the collagen membrane over and secure with pins from the Master-Pin-Control kit to mimic Dr. Urban's sausage technique\*.



7 Primary, tension-free closure is achieved utilizing a horizontal mattress and single interrupted stitch.



8 Healing after about 7 months. Observe the dense, stable bone and the patient is now ready for implant placement.



9 Final bone volume achieved. Patient is now ready for implant placement.

\* Efficacy of the Sausage Technique in Rebuilding the Crestal Buccal Bone Thickness: A Retrospective Analysis – [bit.ly/4kYKoJR](https://bit.ly/4kYKoJR)

### Nobel Biocare regenerative products used in this case




creos™  
allo.gain



creos™  
xenogain



creos™  
xenoprotect



RESORBA®  
GLYCOLON™



Pro-fix™  
Tenting kit



Pro-fix™  
Tenting screws



Master-  
Pin-Control

MKT-6035 Rev 01 (06/25, (GMT 97257)) © Nobel Biocare USA, LLC, 2025. All rights reserved. Distributed by: Nobel Biocare. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. Nobel Biocare does not take any liability for any injury or damage to any person or property arising from the use of this clinical case. This clinical case is not intended to recommend any measures, techniques, procedures or products, or give advice, and is not a substitute for medical training or your own clinical judgment as a healthcare professional. Viewers should never disregard professional medical advice or delay seeking medical treatment because of something they have seen in this clinical case. Full procedure is not shown. Certain sequences have been cut. creos™ allo.gain is manufactured by Solvita, formerly Community Tissue Services (CTS), 2900 College Drive, Kettering, OH 45420, USA, and distributed by Nobel Biocare. creos™ xenogain is manufactured for/by Nibec Co. Ltd., Chungcheongbuk-do, 27816 Korea, and distributed by Nobel Biocare. creos™ xenoprotect is manufactured for/by Matricel GmbH, Kaiserstrasse 100, 52134 Herzogenrath, Germany, and distributed by Nobel Biocare. Pro-fix™ is manufactured by Osteogenics Biomedical, Inc., 4620 71st Street, Bldg 78-79, Lubbock, TX 79424, USA, and distributed by Nobel Biocare. RESORBA® GLYCOLON™ Absorbable Sutures are manufactured by Resorba® Repair and Regenerate an Advanced Medical Solutions Group plc company, Am Flachmoor 16, D-90475 Nuremberg, Germany, and distributed by Nobel Biocare. The Master-Pin-Control Bone Management® system is manufactured by Hager & Meisinger GmbH, Hansemanstr. 10, 41468 Neuss, Germany, and distributed by Nobel Biocare.