



Nobel Biocare N1™

Guided SURGERY

Guided Pilot Drill

The osteotomy is created using either the OsseoDirector or Guided Pilot Drill, and OsseoShapers. OsseoShapers are threaded devices that are inserted and removed at low speed without irrigation. They replace the usual drills used for the creation of the osteotomy.

1

The Guided Pilot Drill can be used as alternative or before using the OsseoDirector.

💧 ≤2000 rpm



Caution The OsseoDirector, Guided Pilot Drill and Twist Step Drill extend up to 1 mm longer than the implant length when seated. Allow for this additional length when drilling near vital anatomical structures.

Warning The OsseoDirector, Guided Pilot Drill, OsseoShapers and Twist Step Drills are sharp instruments. Handle with care to prevent injury.

2

Seat the surgical template.



For detailed instructions on NobelGuide components refer to Nobel Biocare IFU2001 and IFU2009.

3

Engage the Guided Pilot Drill with the handpiece and drill the full planned depth.

The Guided Pilot Drill must proceed at high speed, maximum 2000 rpm and under constant and profuse irrigation by sterile saline at room temperature.

💧 ≤2000 rpm



4

Remove the surgical template and proceed with the OsseoDirector to widen the osteotomy or directly to the OsseoShaper 1.



Order online

Order our complete range of implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

Order by phone

Call our customer service team or contact your sales representative.

Lifetime warranty

The warranty covers all Nobel Biocare implants including prefabricated prosthetic components.

nobelbiocare.com/warranty



Manufacturer

Nobel Biocare AB Box 5190, 402 26
Västra Hamngatan 1
411 17 Göteborg Sweden
www.nobelbiocare.com

Distributed in USA by

Nobel Biocare USA, LLC
22715 Savi Ranch Parkway
Yorba Linda, CA, 92887 USA

Caution Federal law restricts this device to sale by or on the order of a licensed physician or dentist

Note In order to improve readability, Nobel Biocare does not use [™] or [®] in the running text. By doing so, however, Nobel Biocare does not waive any right to the trademark or registered mark and nothing herein shall be construed to the contrary.

Disclaimer Some products may not be regulatory cleared, released or licensed for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability.



nobelbiocare.com/n1



GMT 79630 GB 2204 © Nobel Biocare Services AG, 2022. All rights reserved. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.ons.