

Anterior restoration using advanced tissue grafting for improved esthetics

Juan Zufía
Spain



Juan Zufía

Patient

32, male

Surgical solution

NobelReplace CC PMC 3.5 x
13 mm Slim Healing Abutment
CC NP 7 mm Temporary
Abutment Engaging CC NP

Restorative solution

NobelProcera ASC Abutment
with zirconia crown

Surgery date

July 22, 2014

Tooth position

Upper-right central incisor

Total treatment time

Six months

"Currently, one of the most important advances in implant dentistry is the improvement in treatment of gingival tissue around the implant.

The grafting of connective tissue in the cervical area of the implant improves esthetics and long-term marginal sealing. In this sense, the slim abutment increases the chances of success, allowing the simultaneous placement and stabilization of the gingival graft.

Slim abutments maximize space for grafting, reduce treatment time and improve the final result."

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Initial Situation: 32-year-old patient. Central incisor was lost through trauma two years prior.
Soft tissue is lacking both horizontally and vertically.
Overall Health: Very good. No medication. Non-smoker.



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



The NobelClinician Software is used for accurate implant treatment planning. Because of the narrow width of the ridge and presence of the nasopalatine nerve, placement of a NobelReplace CC PMC 3.5 x 13 mm is the most appropriate solution for rehabilitation of this region.



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Preview of surgical template in the NobelClinician Software allows for the template to be visualized before sending order for production.


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



The flap is raised using a microblade.

Partial thickness flap procedure is performed and the bone is not exposed.


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Surgical template is placed in the correct position after raising the flap. It fits perfectly in the patient's mouth. No instability of the surgical template could be detected.

The Guided Twist Drill $\text{Ø } 2.0$ mm is used to prepare the site for the final depth.


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Preparation of the osteotomy is carried out with \varnothing 3.5 mm guided drill.


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



A NobelReplace CC PMC 3.5 x 13 mm implant is placed into the osteotomy with using Implant Driver CC NP for Slim Abutment.

The marker on the driver indicates the correct depth of the implant.



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



A connective tissue graft from the palate is used to improve the soft tissue environment around the implant. The patient had no tuberosity at all.

The connective tissue graft is prepared so that it covers the implant from the buccal side to the palatal side.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



The prepared palatal connective tissue graft is then placed on the Slim Healing Abutment.


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



The Slim Healing Abutment is connected to the implant.

The graft covers the edentulous space completely.



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



A monofilament suture is used to displace the tissue coronally.

The suture must generate as little tension as possible on the surrounding soft tissue to avoid ischemia of the flap and to allow the revascularization of the graft.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Situation two weeks after the surgery.

After three months of healing the tissue has grown vertically and horizontally.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Occlusal view three months after the surgery, the soft tissue contour developed correctly.



The Slim Healing Abutment is removed to be replaced by a larger-diameter healing cap.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



The healing cap will work as a transition between the use of Slim Healing Abutment and the provisional restoration.

Provisional restoration is placed after three months. Screw retained Temporary Abutment Engaging CC NP is used. A screw-retained final restoration is to be placed at a later date.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Six months after surgery, a NobelProcera ASC Abutment and a zirconia crown are placed as the final restoration .

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



To ensure symmetry, a composite restoration is performed on the left central incisor.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Final clinical picture six months after the surgery.

Case courtesy of Juan Zufía



GMT 72420 GB 2011 © Nobel Biocare Services AG, 2020. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://www.nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. Nobel Biocare does not take any liability for any injury or damage to any person or property arising from the use of this clinical case. This clinical case is not intended to recommend any measures, techniques, procedures or products, or give advice, and is not a substitute for medical training or your own clinical judgement as a healthcare professional. Viewers should never disregard professional medical advice or delay seeking medical treatment because of something they have seen in this clinical case. Full procedure is not shown. Certain sequences have been cut.