

Nobel Biocare N1™ concept: biologically friendly innovation reduces cell death



Clinical relevance

Preserves the vital bone

and reduces cell damage, demonstrated to promote osseointegration.

Simplified osteotomy protocol

with most implants placed after only two drilling steps.

Successful implant placement

with predictable insertion torque even in sites with poor bone quality and quantity.




An Osteotomy Tool That Preserves Bone Viability: Evaluation in Preclinical and Clinical Settings. Bahat O, Yin X, Holst S, Zabalegui I, Berroeta E, Pérez J, Wöhrle P, Sörgel N, Brunski J, Helms JA. *Journal of Clinical Medicine*. 2022;11(9):2536.


The results: pre-clinical

- Less bone damage and cell death compared to conventional protocol.
- Osteoconductive bone coagulum retained in OsseoShaper osteotomy.

OsseoShaper




Conventional drill



Bone coagulum visibly retained in osteotomy with OsseoShaper (a). Conventional osteotomy is devoid of cells (b).

The results: clinical

- 14 of the 15 (93.3%) osteotomies were prepared in only two drilling steps.
- All implants accurately followed the osteotomies.
- Predictable mean final implant insertion torque of 46.9 ± 10.3 Ncm, despite poor bone quality and quantity at all sites.



Study design

Animal model and retrospective clinical case series

Histomorphometric analyses - 36 rodents




7 patients

Challenging bone conditions

15 implants

Delayed loading

Sample clinical case from the study



Images courtesy of Dr. Oded Bahat, USA

A 38-year-old female patient required rehabilitation of a missing maxillary central incisor. Intraoral exam and CBCT showed a pronounced buccal concavity (a, b) and soft bone (b). The osteotomy was created with the OsseoDirector (c) followed by OsseoShaper (d) at low speed (50 rpm) and without irrigation. Note the bone chips within the OsseoShaper threads upon its withdrawal from the osteotomy (e). A Nobel Biocare N1 concept implant was placed with the final insertion torque of 45 Ncm (f, g). Clinical view at 3 years (h).

Authors' conclusion

This innovation in implant site preparation preserves vital anatomical and cellular structures while simplifying surgical protocols.

More to explore



Read on
PubMed



FOR - Safety and
effectiveness of
the N1 Concept



LinkedIn
Dr. Bahat



nobelbiocare.com/N1

GMT 81062 © Nobel Biocare Services AG, 2022. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://www.nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions

**Nobel
Biocare™**