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
Quellenangaben

1. Milleret V, Lienemann PS, Gasser A, Bauer S, Ehrbar M, Wennerberg A. Rational design and in vitro characterization of novel dental implant and abutment surfaces for balancing clinical and biological needs. Clin Implant Dent Relat Res. 2019 Mar;21 Suppl 1:15-24.
2. Robust papilla regeneration and marginal bone gain: 2-Year results from a prospective study. Fabbri G, Ban G, Lim HC, Noh K. Clin Oral Impl Res;34(S27):159-160.
3. All-on-4 concept using TiUltra surface implants and Multi-unit Xeal abutments: Report on the outcomes between 6 to 36 months. Ferro A, De Araújo Nobre M. Clin Oral Impl Res;34(S27):137-138.
4. A Retrospective Observational Study Assessing the Clinical Outcomes of a Novel Implant System with Low-Speed Site Preparation Protocol and Tri-Oval Implant Geometry. Fabbri G, Staas T, Urban I. J Clin Med. 2022 Aug 18;11(16):4859.
5. Navigation guided socket-shield technique for implant and pontic sites in the esthetic zone: A proof-of-concept 1-year prospective study with immediate implant placement and loading. Pozzi A, Arcuri L, Kan J, Londono J. J Esthet Restor Dent. 2022 Jan;34(1):203-214.
6. Dynamic navigation guided surgery and prosthetics for immediate loading of complete-arch restoration. Pozzi A, Hansson L, Carosi P, Arcuri L. J Esthet Restor Dent. 2021 Jan;33(1):224-236.
7. Evaluation of gradually anodized surface on four implant systems. A prospective, thousand-patient real-world study. Eshraghi et al. Poster presented at: AO 2024 annual meeting; March 7-9 2024; Charlotte, NC, USA.
8. Dynamic Navigation System for Immediate Implant Placement in the Maxillary Aesthetic Region. Battista E, Gasparro R, Cacciola M, Sammartino G, Marenzi G. Applied Sciences. 2022;12(11):5510.
9. Donker VJJ, Raghoobar GM, Vissink A, Meijer HJA. Digital Workflow for Immediate Implant Placement and Chairside Provisionalization in the Esthetic Zone. Case Rep Dent. 2022 Apr 1;2022:5114332.
10. Susin C, Finger Stadler A, Muszkopf ML, de Sousa Rabelo M, Ramos UD, Fiorini T. Safety and efficacy of a novel, gradually anodized dental implant surface: A study in Yucatan mini pigs. Clin Implant Dent Relat Res. 2019 Mar;21 Suppl 1:44-54.
11. Schupbach P, Glauser R, Bauer S. Al2O3 Particles on Titanium Dental Implant Systems following Sandblasting and Acid-Etching Process. Int J Biomater. 2019 Jun 2;2019:6318429.
12. Schupbach P., Glauser R., 2019, "Preventing aluminium oxide contamination of dental implant surfaces" Research OUTREACH (111). Verfügbar unter: <https://researchoutreach.org/articles/preventing-aluminium-oxide-contamination-of-dental-implant-surfaces> (Zugriff am 28.05.2024).
13. Fischer NG, Aparicio C. Junctional epithelium and hemidesmosomes: Tape and rivets for solving the "percutaneous device dilemma" in dental and other permanent implants. Bioact Mater. 2022 Mar 19;18:178-198. doi: 10.1016/j.bioactmat.2022.03.019.
14. Fabbri G, Sorrentino R. A Biologically Driven Concept to Design the Emergence Profile Around Dental Implants: Surgical and Prosthetic Considerations to Optimize Hard and Soft Tissue Integration. Int J Periodontics Restorative Dent. 2021 Nov-Dec;41(6):913-921.
15. Nossowitz M, Teale M, Mathes S, Venturato A, Gasser A. Evaluation of anodized surfaces designed for improved soft tissue integration. Foundation for Oral Rehabilitation (FOR) 2019, pp. 1-7. www.for.org/scientificarticle/001
16. Sasada Y, Cochran DL. Implant-Abutment Connections: A Review of Biologic Consequences and Peri-implantitis Implications. Int J Oral Maxillofac Implants. 2017 Nov/Dec;32(6):1296-1307.
17. Urban IA, Tattan M, Ravid A, Saleh MH, Tavelli L, Avila-Ortiz G. Simultaneous Alveolar Ridge Augmentation and Periodontal Regenerative Therapy Leveraging Recombinant Human Platelet-Derived Growth Factor-BB (rhPDGF-BB): A Case Report. Int J Periodontics Restorative Dent. 2022 September/October;42(5):577-585.

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89339 DE 2410 © Nobel Biocare Services AG, 2024. Alle Rechte vorbehalten. Vertrieb durch: Nobel Biocare.