





























Prosthetic components

<p>Plastic Temporary Abutment Plastic Temporary Coping Healing Abutment (incl. Healing Abutment Anatomical PEEK and Slim Healing Abutment) On1 (IOS) Healing Cap Healing Cap Multi-unit</p>						<p>By hand</p>
<p>Slim Temporary Abutment Prosthetic Screws for restorations on Multi-unit Abutment level Titanium Temporary Coping Multi-unit</p>				 <p>15 Ncm</p>		
<p>Ball Abutment 17° Multi-unit Abutment 30° Multi-unit Abutment</p>						
<p>Abutments for NobelActive 3.0</p> <p> Max 15 Ncm</p>						
<p>All other temporary and final abutments for Nobel Biocare implant systems (incl. Temporary Abutment Anatomical PEEK)</p>						 <p>35 Ncm</p>
<p>On1 temporary and final components</p>						

Note: Prosthetic components for other implant systems may require different torque values. Always consult the respective instructions for use.