

Step-by-step digitally planned esthetic restoration

Igor Ristić and
Nikola Vasić
Serbia

Igor Ristić and
Nikola Vasilic

Patient

Female, 43

Clinical situation

loose frontal incisor, bit dislocated over time

Diagnosis

vertical root fracture with cast metal post
and PFM crown on upper left central incisor

Surgical solution

NobelActive® NP

Restorative solution

Temporary: Temporary Snap Abutment
and TempShell

Final: NobelProcera® ASC Abutment with
e.max Ceram veneering and thin feldspatic
veneers on 12,11, 22 FDI

Surgery date

Dec 10 / 2016

Final prosthetic delivery

Jun 21 / 2017

Total treatment time

6 months

Patient with vertically fractured root of
an upper central incisor underwent digital
planning and immediate narrow platform
implant placement in the fresh extraction
socket. Connective tissue graft was made
to achieve optimal result.

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Visible dislocation of upper left central incisor, creating esthetic and functional problems

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Intraoral pre-treatment situation



Intraoral pre-treatment situation - lateral view

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



(CB)CT image of the full arch



Bimaxillary intraoral scan with digitally removed upper-left central incisor


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



SmartFusion of a (CB)CT scan and intraoral scans in NobelClinician® software

SmartSetup – a digital wax-up proposed by the software of upper left central incisor in NobelClinician® software

Initial clinical situation



Treatment planning

Surgical procedure

Restorative procedure

Outcome

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic



Analysis of sagittal aspect of tooth, soft tissues and bone



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Digital position of implant in planning phase, Temporary Snap Abutment included


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Digital design of surgical template in NobelClinician® software

Printed NobelGuide® template ready for surgical procedure



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Digitally designed TempShell in DTX Studio™ design software

TempShell milled in local dental laboratory, made from PMMA

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Intrasulcular incision



Forceps adapted to the crown for efficient luxation of the tooth

Initial clinical situation

Treatment planning


Surgical procedure

Restorative procedure

Outcome



Atraumatic removal of the remaining root structures



Cleaning up of the alveolar socket, removing all granulation tissues

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Try-in of the surgical template intraorally – occlusal view

Try-in of the surgical template intraorally – frontal view

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Guided Twist Drill (10+) 7-13 mm, Ø 2 mm

Initial drilling using Guided Twist Drill*

* As this was a case with immediate implant placement in an extraction site, the Screw Tap was not used, due to the size of the alveolar socket


Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



NobelActive® NP 15mm

Implant insertion using the NobelGuide® surgical template

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Final implant positioning

Verification of implant position according to the surgical template

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Harvesting of the connective tissue graft from tuberosity



Connective tissue graft

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Try-in of the TempShell – vestibular view



Try-in of the TempShell – vestibular view



Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Extraoral picture of the Temporary Snap Abutment

Intraoral pictures after the Temporary Snap Abutment try-in


Initial clinical situation

Treatment planning


Surgical procedure

Restorative procedure

Outcome



Intraoral position of Temporary Snap Abutment is as previously planned



Planning position of Temporary Snap Abutment

**Initial clinical
situation**

**Treatment
planning**

**Surgical
procedure**

**Restorative
procedure**

Outcome



Shortening of the Temporary Snap Abutment

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Loading of light cure material for TempShell relining

Positioning of the TempShell

Light curing and setting in place of TempShell on the Temporary Snap Abutment

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Temporary Snap Abutment in TempShell after polymerization in right position



Extraoral final contouring of emergence profile on temporary crown

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Finishing and polishing of emergence profile



Removing wings for final placement of TempShell

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Finalizing Temporary Snap Abutment for fixation



Drilling the screw access hole with Apical Drill

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Widening the screw access hole with a diamond burr



Placing of temporary crown

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Torque value is 15 Ncm for Temporary Snap Abutment

Immediate postoperative situation with temporary crown

Initial clinical situation



Soft tissue final re-shaping prior to impression procedure

Treatment planning



Minimal preparation of surrounding teeth

Surgical procedure

Restorative procedure



Individualized impression coping prior to impression procedure

Outcome

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

Courtesy of Dr. Igor Ristić and Dr. Nikola Vasilic

NobelProcera® ASC Abutment and veneers on model

Restorations on model, palatal view

Individualized NobelProcera® ASC Abutment in zirconia veneered with e.max Ceram

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Final torque of Individualized NobelProcera® ASC Abutment is 35 Ncm

Implant restoration immediately after tightening

Rest of the restorations bonded on teeth. Screw access hole filled with Teflon™ and light cured composite

Initial clinical situation

Treatment planning


Surgical procedure

Restorative procedure

Outcome



Postoperative intraoral X-ray



Superimposition of the digital planning as previously done in NobelClinician® software

Initial clinical situation



Mock-up for finalization of the treatment. Veneers on upper central right and lateral incisors and upper left lateral incisor and canine are planned

Treatment planning



Intraoral view of the mock-up

Surgical procedure

Restorative procedure

Outcome



Lateral mock-up view

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



Intraoral view of finished NobelProcera® ASC Abutment



Perfect blend of restorative materials with natural tissues, immediate postoperative photo

Initial clinical situation

Treatment planning

Surgical procedure

Restorative procedure

Outcome



X-ray with implant and restoration after final prosthetic delivery



Final smile of the satisfied patient

Case courtesy of Igor Ristić and Nikola Vasilic

