

NobelPearl™

Abutment Screw Retrieval Kit




Procedure for screwed implant abutment connections after screw failure



1. Place the Rescue Shell XT Regular or Wide on the implant and make sure it is seated securely. A light-cured composite (e.g. primopattern LC) can be used to fix the Rescue Shell in connection with the adjacent teeth.
2. Using the green handpiece, intermittently drill the centre of the screw fragment using the Rescue Pilot Drill (Komet item No. H1SE.205.010) at a speed of 2,000 rpm, using as little pressure as possible. The pilot hole should be 2-3 mm in depth. Remove the Rescue Shell, then suction and clean the hole. Note: Use a dental dam if suctioning the particles is not sufficient.
3. Attach the left-hand drill (Rescue Drill «left», Komet item No. H1L.204.014) to the ratchet adapter. Do not use a handpiece! Place the drill vertically on the pilot hole and push it slightly into the pilot hole. Turn the adapter by hand to the **left** (counter-clockwise) while maintaining pressure on the drill. The drill should engage with the screw fragment so that it can be removed with minimal effort. (If the fragment comes loose from the drill after being rotated, it can be removed with an appropriate gripper.

If the drill fails to engage with the fragment or the fragment cannot be loosened, repeat the procedure while making sure that the pilot hole is centric and sufficiently deep. The pilot hole can also be drilled without the Rescue Shell if it is considered to be an obstacle and you have a clear view of the fragment.*

4. Clean the interior of the implant with the cleaning brush and handpiece. Make sure that there are no black residues in the implant.
Caution: Clean the inside of the implant with a brush (e.g. periodontal brush, pressurised air...)
5. Please contact our specialist consultants if you would like to review the reason for the screw failure.



If the fragment cannot be removed using the left-hand drill, the pilot drill can be used to drill out/mill the entire fragment as a last resort. Make sure that adequate cooling is provided.

www.nobelbiocare.com/pearl

GMT 78588 GB 2202 © Nobel Biocare Services AG, 2022. All rights reserved. Distributed by: Nobel Biocare. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/ trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. NobelPearl™ is manufactured for/by Dentalpoint AG, Bodenackerstrasse 5, CH 8957 Spreitenbach, Switzerland and distributed by Nobel Biocare Services AG.