


The logo for X-GUIDE, featuring a stylized 'X' with a green and yellow gradient, followed by the word 'GUIDE' in a bold, dark blue sans-serif font.The Nobel Biocare logo, consisting of a large, stylized 'N' in a dark blue color, followed by the words 'Nobel Biocare' in a dark blue sans-serif font, with a trademark symbol.

# Perform edentulous cases with X-Guide®

with the EDX system

Navigate a wide range of cases including terminal dentition, full-arch and dual full-arch

Secure and stable attachment

Patient tracker remains out of the way during surgery

Special bone-screw design – flapped or flapless

Different arm designs based on handedness of the user


# Elevate your X-Guide game

## Full-arch workflow with X-Mark™

1. Scan and plan
2. Attach patient tracker using edentulous hardware
3. Virtually register the patient with X-Mark
4. Perform navigated surgery

### EDX-s series

## Selected components from our experts



### 2x EDX Screws and EDX Nuts

XNVP010921 (8mm)

Check our online store to see all 10 variations of arms and 6 different screw sizes.



[nobelbiocare.com/x-guide](https://nobelbiocare.com/x-guide)



MKT-5498 Rev 00 (11/23, GMT 87803) © Nobel Biocare USA, LLC, 2023. All rights reserved. Distributed by: Nobel Biocare. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a licensed dentist or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. X-Guide is manufactured by X-Nav Technologies, LLC and distributed by Nobel Biocare.