

# OsseoSet™ 300

SIMPLIFY  
the treatment  
process

Product overview



# OsseoSet 300

Simplifies the treatment process.

## Easy to use

Color touchscreen with user-friendly menu guidance and a glass surface for easy cleaning.

## For the dental assistant

The pump design makes inserting the coolant hose quick and secure.

## Clear, bright and wide field of view

Excellent lighting conditions facilitate accurate treatment results.

## Secure placement

The automatic torque control provides accurate torque delivery. It can be adjusted from 5 to 80 Ncm.

## Ideal for multiple users

Drill protocols can be customized for multiple surgeons in the practice.

## Automatic thread cutting function

Supports you when placing implants in dense bone. It minimizes unnecessary compression of the bone.



# OsseoSet 300 product line

## OsseoSet 300

87289	OsseoSet 300 1015, NW, 19LC, WS-75, 120V <sup>1</sup>
87290	OsseoSet 300 1015, N2, 19LC, WS-75, 120V <sup>2</sup>



## Accessories for OsseoSet 300

NB30033001	Surgical Contra-angle WS-75 L, 20:1
NB30056000	Surgical Handpiece S-11, straight 1:1
NB30058001	Surgical Handpiece S-11 L, straight 1:1
NB30095000	Surgical Handpiece Zygoma SZ-75, 20:1
NB30281000	Motor with Cable 1.8 m EM-19 LC
NB04363600	Irrigation Tubing Set 2.2 m, for OsseoSet 300 (6/pkg)
NB02610500	Cannula Internal Cooling
NB07962790	Transportation Case
NB30285001 or NB30285000 <sup>3</sup>	Foot control S-N2 CAN-Bus
NB30264002 or NB30264000 <sup>3</sup>	Foot control S-NW CAN-Bus (wireless)



<sup>1</sup> Includes OsseoSet 300 drill unit, motor with cable 1.8 m EM-19 LC, surgical contra-angle WS-75 L, 20:1 and wireless foot control.

<sup>2</sup> Includes OsseoSet 300 drill unit, motor with cable 1.8 m EM-19 LC, surgical contra-angle WS-75 L, 20:1 and foot control.

<sup>3</sup> Depending on market.

---

## ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[store.nobelbiocare.com/us/en/](https://store.nobelbiocare.com/us/en/)

---

## ORDER BY PHONE

Contact your sales representative, or call our customer service team.

**800 322 5001**

---

## LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

[nobelbiocare.com/warranty](https://nobelbiocare.com/warranty)

[nobelbiocare.com/drillunits](https://nobelbiocare.com/drillunits)



MKT-5404 Rev 00 (07/23, GMT 85936). © Nobel Biocare USA, LLC, 2023. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed dentist or physician. See Instructions for Use for full prescribing information, including indications, contraindications, warnings and precautions.