

NobelProcera® Scan and Design Services Implant Bar Overdenture

Platform Master List and Screw Identification Chart

Manufacturer	Platform	Clinical screw reference number from legal manufacturer
Nobel Biocare		
NobelActive®, NobelParallel™ CC and NobelReplace® Conical Connection	Internal conical connection NP	37891*
	Internal conical connection RP/WP	37892*
NobelReplace®, Replace Select™, NobelSpeedy® Replace, NobelReplace® Platform Shift	Internal tri-channel connection NP	36818*
	Internal tri-channel connection RP	29475*
	Internal tri-channel connection WP	29475*
	Internal tri-channel connection 6.0	29475*
Brånemark System® and NobelSpeedy® Groovy	External hex connection NP	29282*
	External hex connection RP	29283*
	External hex connection WP	29284*
	Multi-Unit RP	29285*
	Multi-Unit WP	29286*
	Steri-oss HL 3.25	2380*
	Steri-oss HL 3.8/4.5	2380*
	Replace Hex 3.5	29283*
	Replace Hex 4.3	29283*
	Steri-oss HL/Replace Hex 5.0	29283*
	Steri-oss HL/Replace Hex 6.0	29283*
	Internal Hexagon	32983*
	Nobel Novum	27266*
AB Dental Devices	Narrow Platform	P3-al
Alpha-Bio Tec	Internal Hex 4.5	STLAS
	TCT System	6090
Anthogyr	Anthofit™ HE System Conical Abutment 4.1	OIT1002
	Anthofit™ OI System 3.5	OICA005
	Anthofit™ OI System Conical Abutment 3.5	OIT1002
	Axiom™ REG/PX Platform 4.0	OPTS160
	Axiom™ REG/PX Straight Conical Abutment	OPSC140

* Clinical screws included with Implant Bar Overdenture

Manufacturer	Platform	Clinical screw reference number <small>from legal manufacturer</small>
BEGO Semados®	S/RI-Line 3.25/3.75/4.1/4.5	53999
	S/RI-Line 5.5	53999
BioHorizons™	Abutment for Screw 3.5/4.0	PXABS
	Abutment for Screw 5.0	222-110
	Conical internal connection 3.5	PXAS
	Conical internal connection 4.5	PXAS
	Conical internal connection 5.7	PXAS
	External implant system 3.5	130-300
	External implant system 4.0	140-300
	External implant system 5.0	140-300
	External implant system 6.0	140-300
	Single Stage Restoration 3.5	PXAS
	Single Stage Restoration 4.5	PXAS
	Single Stage Restoration 5.7	PXAS
Biomet 3i®	Certain® Internal 3.4	IUNIHT
	Certain® Internal 4.1	ILRGHT and ILRGHG
	Certain® Internal 5.0	ILRGHT and ILRGHG
	Certain® Internal 6.0	ILRGHT and ILRGHG
	Conical Abutment 3.4	GSH30
	Conical Abutment 4.1/5.0	GSH30
	Conical Abutment 6.0	GSH30
	External Hex 3.4	UNIHT
	External Hex 4.1	UNIHT
	External Hex 5.0	UNIHT
	External Hex 6.0	UNIHT
	IOL® Abutment	GSH30
	Low Profile abutments	LPCTSH
	Standard abutment	GSH30
Biotech Dental Implants	B.I.S. Conic 3.5/3.9	FMC39452R
	B.I.S. Conic 4.0/4.4	FMC39452R
	B.I.S. Conic 5/5.4	FMC39452R
	B.I.S. Conic Base Analog 4.4	E007
BTI™	Externa Platform 4.1	TTUH
	Externa Platform 4.1 (abutment level)	TTUH
	Interna Platform 4.1	INTTUH
	Tiny™ Platform 3.5	TTTH
Camlog®	Camlog® Bar Abutment 3.3	J4005.1602
	Camlog® Bar Abutment 3.8	J4005.1602
	Camlog® Bar Abutment 6.0	J4005.1602
	Camlog® Implant System 3.3	J4005.1601
	Camlog® Implant System 3.8	J4005.1601
	Camlog® Implant System 4.3	J4005.1601
	Camlog® Implant System 5.0	J4005.2001
	Camlog® Implant System 6.0	J4005.2001
	CONELOG® Implant System 3.8	C4006.1601

Manufacturer	Platform	Clinical screw reference number <small>from legal manufacturer</small>
Dentsply® (Astra Tech Implant System™)	20° UniAbutment 3.5/4.0	22435
	45° UniAbutment 3.5/4.0	22435
	Astra Aqua	24449
	Astra Lilac	24209
Dentsply® (XiVE®)	FRIADENT 3.0	46-4304
	FRIADENT 3.4	46-4305
	FRIADENT 3.8	46-4305
	FRIADENT 4.5	46-4305
	FRIADENT 5.5	46-4305
	Friadent MP 3.8	25-4200
	Friadent MP 4.5	25-4200
	Friadent MP 5.5	25-4200
	Xive® TG Platform 3.4-4.5	25-4200
Implant Direct (Attachment International)	3.25	45-300052
	3.3	48-400052
Implant Direct (Sybron)	Endopore® Internal Hex 4.1	05B-2RS
	Entegra™ External Hex 3.5	05-2RS
	Entegra™ External Hex 4.1	05B-2RS
	Entegra™ External Hex 5.0	05B-2RS
	PITT-EASY® Paracentric Line 3.75	pitt easy screw 90270
	PITT-EASY® Paracentric Line 4.0	pitt easy screw 90270
	PITT-EASY® Paracentric Line 4.9	pitt easy screw 90270
	SybronPRO™ Series 3.5 Hex	5X-184
	SybronPRO™ Series 4.1 Hex	5X-184
LifeCore Dental (Keystone)	Genesis System 3.8	G71127
	Genesis System 4.5	G71127
	Genesis System 5.5	G71127
	PrimaConnex® System 3.5	45060
	PrimaConnex® System 4.1	45060
	PrimaConnex® System 5.0/7.0/8.0	45060
	PrimaConnex® System Multi-unit 3.5/4.1	45164
	PrimaConnex® System Multi-unit WD	45164
	Renova™ Implant System RDI 4.5mm	43116
	Renova™ Implant System SDI 3.75mm	43116
	Restore™ Implant System WD	R9203-50-48
	Restore™ System RD 3.75/4.0	R9202-40-48
	Restore™ System SD 3.3	R9203-34-48
	Stage-1 system RDS	S2408K
	MIS®	Internal Hex
Multi-Unit system		MM-S0330
Standard Platform		EA-GH480
Wide Interier Hex		MD-S0219
Neoss®	Access Abutment	31221
	Neoss Implant System 3.25	31150

Manufacturer	Platform	Clinical screw reference number <small>from legal manufacturer</small>
Osstem™	TS & GS System mini 4.0	GSASM
	TS & GS System standard 4.5	GSASR
	US System Esthetic Low	MTS200
	US System Regular platform	ASR200
SIC Implants	Sic Platform 3.3	936658
	Sic Platform 4.2	936658
Simpler Implants	Implant Analog 3.25	S17013
	Implant Analog 4.0	S17014
Southern Implants	External Hex 4.0	TSH3
	External Hex 5.0	TSH3
Straumann®	BoneLevel NC (3.3)	36017*
	BoneLevel RC (4.1/4.8)	36018*
	Multi-base NC 3.5	025.0900-04
	Multi-base RC 4.5	025.0900-04
	Multi-base RC 6.5	025.0900-04
	Octagon (Standard / Standard Plus) 4.8	31984*
	Octagon (Standard / Standard Plus) 6.5	31984*
	Ortho System 4.2	048.350
	SynOcta® RN	048.350
SynOcta® WN	048.350	
Thommen Medical	SPI®RETAIN DIRECT System	4.03.506
	SPI®VARIOmulti System 4.0/4.5/5.0	4.03.507
	VARIOeco system 3.5	4.03.500
	VARIOeco system 4.0	4.03.501
	VARIOeco system 4.5	4.03.501
	VARIOeco system 5.0	4.03.501
	VARIOeco system 6.0	4.03.501
Titan implant	Internal Hex 4.5 platform	C-ITS
	Internal Hex 5.7 platform	C-ITS
Zimmer® Dental	AdVent® 4.5	AVGCS
	Calcitek 3.25	calcitek 3.25 screw
	Calcitek 4.0	calcitek 4.00 screw
	Screw-Vent®	3.5 MHLAS
	Screw-Vent®	4.5 MHLAS
	Screw-Vent®	5.7 MHLAS
	Specta-Cone	SCTS
	Spline® Implant System 3.25	1537
	Spline® Implant System 3.75/4.0	1537
	Spline® Implant System 5.0	1537
	SwissPlus® Implant System 3.8	MHLAS
	SwissPlus® Implant System 4.8	GPCAS
	Tapered abutment 4.5	SCTS

* Clinical screws included with Implant Bar Overdenture



Refer to Nobel Biocare Implant
Bar Overdenture Instruction of
Use for more information.



[nobelbiocare.com](https://www.nobelbiocare.com)



GMT 80300 GB 2208 © Nobel Biocare Services AG, 2022. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://www.nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.