


# PureSet™

Product overview



## Contents

---

### Freehand surgery

Nobel Biocare N1™ PureSet 87293

NobelActive® PureSet 87294

NobelParallel™ CC PureSet 87295

NobelReplace® CC PureSet 87296

NobelSpeedy® Groovy PureSet 87297

Brånemark System® PureSet 87298

Replace Select™ TC PureSet 87299

### Guided surgery

NobelActive® Guided PureSet 87305

NobelParallel™ CC Guided PureSet 87306

NobelReplace® CC Guided PureSet 87307

### Prosthetics

Prosthetic PureSet Basic 87301

Prosthetic PureSet Nobel Biocare N1™ 87302

Prosthetic PureSet 87353

### Accessories

PureSet Tray accessories

# PureSet™

PureSet sets for freehand and guided surgery are available for NobelActive, NobelParallel CC, NobelReplace CC and Nobel Biocare N1.

The PureSet sets NobelSpeedy Groovy, Brånemark System and Replace Select TC are available for freehand surgery.

For even greater efficiency during freehand surgery, two of Nobel Biocare's leading implant systems – NobelActive and NobelParallel CC – can both be stored in a single tray, requiring fewer instruments, thanks to their similar drilling protocols.

# Ease, speed and hygiene, all-in-one with PureSet

## STRAIGHTFORWARD AND INTUITIVE LAYOUTS

Designed to reflect surgical workflows and make the components easy to identify

## SECURE RETAINMENT OF INSTRUMENTS


Surgical instruments are contained and secured using special metal springs

## UNCOMPROMISED HYGIENE

PureSet Trays are made of stainless steel<sup>1</sup> and designed for automated cleaning with few manual steps


## SIGNIFICANT TIME SAVINGS

Significant time savings as tray and instruments remain together during cleaning and sterilization




<sup>1</sup> PureSet Trays are made of stainless steel except for the PureSet plates and grommets.

# Freehand SURGERY



# Nobel Biocare N1™ system 87293



# Nobel Biocare N1™ system

## 87293

INCLUDED IN THE SET

PUR0400 Nobel Biocare N1™ PureSet™ Tray

301074 Radiographic Template Nobel Biocare N1™

301075 Wall Chart Nobel Biocare N1™ PureSet

### Implant Driver Nobel Biocare N1™ TCC NP Short

28 mm 300905



### Implant Driver Nobel Biocare N1™ TCC NP Long

37 mm 300906



### Implant Driver Nobel Biocare N1™ TCC RP Short

28 mm 300907



### Implant Driver Nobel Biocare N1™ TCC RP Long

37 mm 300908



### Direction Indicator Nobel Biocare N1™

300918




### OsseoShaper Extension Nobel Biocare N1™

300923



### Depth Probe Nobel Biocare N1™

300924



### Manual Torque Wrench Surgical Nobel Biocare N1™

301083



### Screwdriver Machine Omnigrip Mini

28 mm 300853



### Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



### Handle for Machine Instruments

29161



# Nobel Biocare N1™ system

87293

**Screwdriver Machine Nobel Biocare N1 Base**  
28 mm 300903

**OsseoShaper Extension NB N1**  
300923

**Direction Indicator NB N1**  
300918

**Depth Probe NB N1**  
300924


**Screwdriver Machine Omnigrip Mini**  
28 mm 300853



Handle for machine instruments  
29161



Manual Torque Wrench Surgical NB N1  
301083



**Implant Driver NB N1 TCC NP Short**  
28 mm 300905

**Implant Driver NB N1 TCC NP Long**  
37 mm 300906

**Implant Driver NB N1 TCC RP Short**  
28 mm 300907

**Implant Driver NB N1 TCC RP Short**  
37 mm 300908

# Nobel Biocare N1™ system

## 87293

NOT INCLUDED IN THE SET  
(available for separate purchase)

### Implant site preparation instruments

301401	OsseoDirector Nobel Biocare N1™ Ø 1.8–2.4 x 8–18 mm
301077	OsseoDirector Nobel Biocare N1™ Ø 1.8–2.4 x 8–14 mm
301081	Guided Pilot Drill Nobel Biocare N1™ Ø 2.0 mm
301078	Ø 3.5 Twist Step Drill Ø 2.5/3.4 x 10–14 mm
301422	Ø 4.0 Twist Step Drill Ø 3.3/3.8 x 8–18 mm
301079	Ø 4.0 Twist Step Drill Ø 3.3/3.8 x 8–14 mm



### Abutment Retrieval Tool Nobel Biocare N1™ TCC

NP RP

300921 300922



### Torque Wrench for Prosthetic Devices

301082 Manual Torque Wrench Prosthetic Nobel Biocare N1™



### Screwdrivers

300852	Screwdriver Machine Omnigrip™ Mini 20 mm
300854	Screwdriver Machine Omnigrip™ Mini 36 mm
300855	Screwdriver Manual Omnigrip™ Mini 20 mm
300856	Screwdriver Manual Omnigrip™ Mini 36 mm
300904	Screwdriver Machine Nobel Biocare N1™ Base 36 mm
301135	Screwdriver Manual Nobel Biocare N1™ Base 28 mm
301136	Screwdriver Manual Nobel Biocare N1™ Base 36 mm
29158	Screwdriver Machine Multi-unit 21 mm
29156	Screwdriver Manual Multi-unit 25 mm





# Nobel Biocare N1™ system

## 87293

### OsseoShaper 2 Nobel Biocare N1™

300892	Ø 3.5	3.5 x 9 mm
300893	Ø 3.5	3.5 x 11 mm
300894	Ø 3.5	3.5 x 13 mm
300895	Ø 4.0	4.0 x 7 mm
300896	Ø 4.0	4.0 x 9 mm
300897	Ø 4.0	4.0 x 11 mm
300898	Ø 4.0	4.0 x 13 mm
301179	Ø 4.0	4.0 x 15 mm



### Bone Mill Nobel Biocare N1™ TCC

	NP	RP
Ø 4.0	300909	
Ø 5.2	300910	300910



### Bone Mill Guide Nobel Biocare N1™ TCC


	NP	RP
Ø 4.0	300911	
Ø 5.2	300915	300916



**NOT INCLUDED IN THE SET**  
(available for separate purchase)

**NobelActive® PureSet**  
87294

**NobelParallel™ CC PureSet**  
87295



## NobelActive® PureSet 87294

PUR0200	NobelActive/NobelParallel CC PureSet Tray
34584	NobelActive Manual Torque Wrench Surgical*
32180	Surgical Driver
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37887	NobelActive Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

### Implant Drivers

	3.0	NP	RP	WP
28 mm	36773	36718	36720	37859
37 mm	36774	36719	36721	37860



## NobelParallel™ CC PureSet 87295

PUR0200	NobelActive/NobelParallel CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37792	NobelParallel Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

### Implant Drivers

	NP	RP	WP
28 mm	36718	36720	37859
37 mm	36719	36721	37860




INCLUDED IN THE SET

\* Article to be disassembled for cleaning

## NobelActive® PureSet 87294

## NobelParallel™ CC PureSet 87295



Depth Probe 7-18 mm  
37791



Precision Drill  
36118\*




Bone Mill (with Guide) CC  
3.0 Ø 4.0 37861A\*  
NP Ø 4.4 37863\*  
RP Ø 5.2 37864\*  
WP Ø 5.2 37866\*  
RP Ø 6.2 37867\*  
WP Ø 6.7 37869\*



Bone Mill Guide CC  
3.0 37862A\*  
NP 37865\*  
RP 37868\*  
WP 37870\*

- \* Article not included in this set
- \*\* Article to be disassembled for cleaning
- \*\*\* Screwdrivers available in other lengths


- A Specific to NobelActive
- P Specific to NobelParallel



Twist Drill  
Ø 1.5 33297A\*  
Ø 2 33297\*



Twist Step Drill  
Ø 2.4, 2.8 mm 32261\*  
Ø 2.8, 3.2 mm 34638\*  
Ø 3.6, 4.0 mm 32275\*  
Ø 3.8, 4.2 mm 32274\*  
Ø 4.2, 4.6 mm 34582\*  
Ø 4.2, 5.0 mm 37874\*



Implant Driver 3.0  
28 mm 36720  
37 mm 36774A



Implant Driver NP  
28 mm 36719  
37 mm 36719



Implant Driver RP  
28 mm 36720  
37 mm 36721




Implant Driver WP  
28 mm 37859  
37 mm 37860




Direction Indicator  
32112



Screwdriver  
manual Unigrip 28 mm  
29149\*\*\*



Shank Extension  
K/589B, 204



Surgical Driver  
32180\*



NobelActive Manual  
Torque Wrench Surgical  
34584A\*\*




NobelReplace Manual  
Torque Wrench Surgical  
28839P\*\*



Screw Tap<sup>A</sup>  
Ø 3.0 36816\*  
NP 3.5 36236\*  
RP 4.3 36237\*  
RP 5.0 36238\*  
WP 5.5 (7-10) 37871\*  
WP 5.5 (11.5-15) 37872\*



Cortical Drill<sup>P</sup>  
NP 3.75 38000\*  
RP 4.3 38001\*  
RP 5.0 38002\*  
WP 5.5 38003\*



Screw Tap<sup>P</sup>  
NP 3.75 (7-13) 37990\*  
NP 3.75 (7-18) 37991\*  
RP 4.3 (7-13) 37992\*  
RP 4.3 (7-18) 37993\*  
RP 5.0 (7-13) 37994\*  
RP 5.0 (7-18) 37995\*  
WP 5.5 (7-10) 37996\*  
WP 5.5 (7-15) 37997\*

PureSet™ A/P			
NobelActive®			
	Dense		
Ø 3.0	1	2	2.3
Ø 3.5	7 (0)	23 (0)	2.34
Ø 4.3	23 (0)	23.5	2.35 (0)
Ø 5.0	23.5	23.56	2.356 (7)
Ø 5.5	2.35 (0)	2.356 (7 (0))	2.356 (8)
	11.5-15		
(L) = Only widen cortex			

NobelParallel™ CC			
	Dense		
Ø 3.75	210	210 (0)	2.34 (0.5)
Ø 4.3	215 (0)	215 (0)	2.35 (0.5)
Ø 5.0	235 (0)	235.6 (0.5)	2.356 (0.5)
Ø 5.5	2.35 (7 (0))	2.356 (0.5)	2.356 (0.5)
(L) = Optional (C) = Cortical Drill    S = Screw Tap			

# NobelActive® PureSet

## 87294

# NobelParallel™ CC PureSet

## 87295

### Drills

Precision	36118			
Ø 1.5	31278 (7–15 mm) <sup>A</sup>			
Ø 2.0	32296 (7–10 mm)	32297 (7–15 mm)	32299 (10–18 mm)	
Ø 2.4/2.8	32260 (7–10 mm)	32261 (7–15 mm)	32262 (10–18 mm)	
Ø 2.8/3.2	37873 (7–10 mm)	34638 (7–15 mm)	34639 (10–18 mm)	
Ø 3.2/3.6	32263 (7–10 mm) <sup>P</sup>	32264 (7–15 mm)	32265 (10–18 mm)	
Ø 3.8/4.2	32275 (7–10 mm) <sup>P</sup>	32276 (7–15 mm)	32277 (10–18 mm)	
Ø 4.2/4.6	37874 (7–10 mm)	34582 (7–15 mm)	34583 (10–18 mm)	
Ø 4.2/5.0	37875 (7–10 mm)	37876 (7–15 mm)		



### Bone Mill Guides

3.0	NP	RP	WP
37862 <sup>A</sup>	37865	37868	37870



### Shank Extension

KI589B.204.



### Screw Taps NobelActive

	3.0	NP	RP	WP
Ø 3.0	36816			
Ø 3.5		36236		
Ø 4.3			36237	
Ø 5.0			36238	
Ø 5.5				37871 (7–10 mm) 37872 (11.5–15 mm)



### Screw Taps NobelParallel CC

	NP	RP	WP
Ø 3.75	37990 (7–13 mm) 37991 (7–18 mm)		
Ø 4.3		37992 (7–13 mm) 37993 (7–18 mm)	
Ø 5.0		37994 (7–13 mm) 37995 (7–18 mm)	
Ø 5.5			37996 (7–10 mm) 37997 (7–15 mm)



NOT INCLUDED IN THE SET  
(available for separate purchase)

### Cortical Drills NobelParallel CC

	NP	RP	WP
Ø 3.75	38000		
Ø 4.3		38001	
Ø 5.0		38002	
Ø 5.5			38003




### Bone Mills with Guide

	3.0	NP	RP	WP
Ø 4.0	37861 <sup>A</sup>			
Ø 4.4		37863		
Ø 5.2		37864	37866	
Ø 6.2			37867	
Ø 6.7				37869



A For NobelActive only  
P For NobelParallel CC only

# NobelReplace® CC PureSet 87296



# NobelReplace® CC PureSet 87296

INCLUDED IN THE SET

PUR0300	NobelReplace CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
2042	Irrigation Needle
37320	NobelReplace Radiographic Template
300567	NobelReplace CC PureSet Wall Chart

## Drills

	NP Ø 3.5	RP Ø 4.3	RP Ø 5.0
8 mm	32075	32076	32077
10 mm	29367	29370	29373
11.5 mm	36113	36114	36115
13 mm	29368	29371	29374
16 mm	29369	29372	29375



## Dense Bone Drills

	NP Ø 3.5	RP Ø 4.3	RP Ø 5.0
13 mm	29377	29380	29383
16 mm	29378	29381	29384



## Direction Indicators

NP	RP	WP
32255	32256	32257



## Screw Taps

NP	RP	WP
36717	32090	32091



## Implant Drivers


	NP	RP
28 mm	36718	36720
37 mm	36719	36721




\* Article to be disassembled for cleaning


# NobelReplace® CC PureSet 87296

  
Precision Drill  
36118\*

  
Drill with Tip  
Tapered 2 mm  
36117\*

  
Bone Mill (with Guide) CC  
NP Ø 4.4 37863\*  
NP Ø 5.2 37864\*  
RP Ø 5.2 37866\*  
RP Ø 6.2 37867\*


  
Bone Mill Guide CC  
NP 37865\*  
RP 37868\*



  
Direction Indicator  
Tapered  
NP 32255  
RP 32256  
WP 32257

  
Implant Driver NP  
28 mm 36718  
37 mm 36719

  
Implant Driver RP  
28 mm 36720  
37 mm 36721

  
Screw Tap Tapered  
NP 3.5 36717  
RP 4.3 32090  
WP 5.0 32091

\* Article not included in this set

\*\* Article to be disassembled for cleaning

\*\*\* Screwdrivers available in other lengths



# NobelReplace® CC PureSet 87296

NOT INCLUDED IN THE SET  
(available for separate purchase)

## Bone Mills with Guide

	NP	RP
Ø 4.4	37863	
Ø 5.2	37864	37866
Ø 6.2		37867



## Bone Mill Guides

NP	RP
37865	37868



## Drills

36117	Drill with Tip Tapered 2 mm	
36118	Precision Drill	

## Shank Extension

KI589B.204.




# NobelSpeedy® Groovy PureSet

87297

# Brånemark System® PureSet

87298



# NobelSpeedy® Groovy PureSet

## 87297

# Brånemark System® PureSet

## 87298

INCLUDED IN THE SET

PUR0900	NobelSpeedy Groovy/Brånemark System/ Replace Select TC PureSet Tray
DIB 097-0	Cover Screw Driver Brånemark System Hexagon
32110	Brånemark System Manual Torque Wrench Surgical*
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37628	NobelSpeedy Groovy Radiographic Template 7–25 mm**
33097	Brånemark System/NobelSpeedy Radiographic Template***
301381	NobelSpeedy Groovy/Brånemark System PureSet Wall Chart

### Implant Drivers

	NP	RP	WP
21 mm		29129	29134**
26 mm	29126**	29130	29135**
34 mm	29127**		



### Implant Driver Wrench


	NP	RP	WP
12 mm		29132	29136**
21 mm	29128**		



\* Article to be disassembled for cleaning  
 \*\* Article only included in NobelSpeedy Groovy PureSet  
 \*\*\* Article only included in Brånemark System PureSet

## NobelSpeedy® Groovy PureSet 87297

## Brånemark System® PureSet 87298



NobelSpeedy Depth Probe 18-25 mm  
37692\*\*

Depth Probe 7-18 mm Z-shaped  
37791



Precision Drill  
36118\*




Bone Mill (with Guide)  
Bmk Syst  
NP Ø 4.5 33392\*  
RP Ø 5.1 33393\*  
WP Ø 6.5 33495\*



Bone Mill Guide  
Bmk Syst  
NP 33496\*  
RP 33497\*  
WP 33498\*



Direction Indicator  
Ø 2/Ø 2.4\*\*  
32112



Screwdriver  
manual Unigrip 28 mm  
29149\*\*



Cover Screw Driver  
Bmk Syst Hexagon  
DIB 097-O



Brånemark System Manual  
Torque Wrench Surgical  
32110\*\*




Brånemark System  
Counterbore  
NP 32281\*  
RP 32283\*  
WP 32285\*




Shank Extension  
KI589B.204.




Screw Tap  
NP Ø 3.75 7-15 mm 37627\*  
RP Ø 3.75 7-13 mm 32290\*\*  
RP Ø 3.75 7-18 mm 32291\*\*  
RP Ø 4.7-13 mm 33072\*  
RP Ø 4.7-18 mm 33073\*  
WP Ø 5.7-13 mm 32292\*  
WP Ø 5.7-18 mm 32293\*  
WP Ø 6.7-13 mm 32294\*\*  
WP Ø 6.7-18 mm 32295\*\*




Drill  
7-10 mm  
Ø 2.0 (#1) 32266\*  
Ø 2.4 (#2) 32267\*  
Ø 2.8 (#3) 32268\*  
Ø 3.2 (#4) 32269\*  
Ø 3.6 (#5) 32270\*  
Ø 3.7/3.6 (#6) 32263\*  
Ø 3.8/4.2 (#7) 32275\*  
Ø 5.0 (#8) 32279\*\*




7-15 mm  
Ø 2.0 (#1) 32297\*  
Ø 2.4/2.8 (#2) 32261\*  
Ø 3.0 (#3) 32267\*  
Ø 3.2 (#4) 32270\*  
Ø 3.4 (#5) 32273\*  
Ø 3.6/4.2 (#7) 32274\*  
Ø 5.0 (#8) 32279\*\*



10-18 mm  
Ø 2.0 (#1) 32299\*  
Ø 2.4/2.8 (#2) 32262\*  
Ø 3.0 (#3) 32271\*  
Ø 3.2 (#4) 32274\*  
Ø 3.4 (#5) 32274\*  
Ø 3.7/3.6 (#6) 32265\*  
Ø 3.8/4.2 (#7) 32277\*  
Ø 5.0 (#8) 32280\*\*



18-25 mm  
Ø 2.0 (#1) 37614\*\*  
Ø 2.4/2.8 (#2) 37615\*\*  
Ø 3.0 (#3) 37616\*\*  
Ø 3.2 (#4) 37616\*\*  
Ø 3.4 (#5) 37617\*\*



Implant Driver Bmk Syst  
NP 26 mm 29126  
RP 26 mm 29126  
RP 34 mm 29129  
RP 26 mm 29130  
RP 34 mm 29131  
WP 21 mm 29134  
WP 26 mm 29135



Implant Driver Wrench  
Adap Bmk Syst  
NP 12 mm 29132  
RP 12 mm 29133  
RP 21 mm 29133  
WP 12 mm 29136

\* Article not included in this set  
\*\* Article to be disassembled for cleaning  
\*\*\* Screwdrivers available in other lengths

B Specific to Brånemark System  
S Specific to NobelSpeedy Groovy

# NobelSpeedy® Groovy PureSet 87297

# Brånemark System® PureSet 87298

**NOT INCLUDED IN THE SET**  
(available for separate purchase)

## Drills

Precision	36118				
Ø 2.0	32296 (7-10 mm)	32297 (7-15 mm)	32299 (10-18 mm)	37614 (18-25 mm) <sup>S</sup>	37619 (18-25 mm) LS <sup>S</sup>
Ø 2.4/2.8	32260 (7-10 mm)	32261 (7-15 mm)	32262 (10-18 mm)	37618 (18-25 mm) <sup>S</sup>	37620 (18-25 mm) LS <sup>S</sup>
Ø 3.0	32266 (7-10 mm)	32267 (7-15 mm)	32268 (10-18 mm)		
Ø 3.2	32269 (7-10 mm)	32270 (7-15 mm)	32271 (10-18 mm)	37616 (18-25 mm) <sup>S</sup>	
Ø 3.4	32272 (7-10 mm)	32273 (7-15 mm)	32274 (10-18 mm)	37617 (18-25 mm) <sup>S</sup>	
Ø 3.2/3.6	32263 (7-10 mm)	32264 (7-15 mm)	32265 (10-18 mm)		
Ø 3.8/4.2	32275 (7-10 mm)	32276 (7-15 mm)	32277 (10-18 mm)		
Ø 5.0	32278 (7-10 mm) <sup>S</sup>	32279 (7-15 mm) <sup>S</sup>	32280 (10-18 mm) <sup>S</sup>		



## Bone Mills (with Guide) Bmk Syst

	NP	RP	WP
Ø 4.5	33392		
Ø 5.1		33393	
Ø 6.5			33495



## Bone Mill Guides Bmk Syst

NP	RP	WP
33496	33497	33498



## Brånemark System® Counterbore

NP	RP	WP
32281	32283	32285




## Screw Taps

	NP	RP	WP
Ø 3.3	37627 (7-15 mm)		
Ø 3.75		32290 (7-13 mm) <sup>B</sup> 32291 (7-18 mm) <sup>B</sup>	
Ø 4.0		33072 (7-13 mm) 33073 (7-18 mm)	
Ø 5.0			32292 (7-13 mm) 32293 (7-18 mm)
Ø 6.0			32294 (7-13 mm) <sup>S</sup> 32295 (7-18 mm) <sup>S</sup>



## NobelSpeedy Depth Probe 18-25 mm

37692<sup>S</sup>



B Specific to Brånemark System  
S Specific to NobelSpeedy Groovy  
\* Article to be removed prior to placing lid on tray

# NobelSpeedy® Groovy PureSet

87297

# Brånemark System® PureSet

87298

**NOT INCLUDED IN THE SET**  
(available for separate purchase)

## Implant Drivers

	NP	RP	WP
21 mm			29134*
26 mm	29126*		29135*
34 mm	29127*	29131	



## Implant Driver Wrench

	NP	RP	WP
12 mm			29136*
21 mm	29128*	29133	




## Shank Extension

KI589B.204.



\* Article only included in NobelSpeedy Groovy PureSet

# Replace Select™ TC PureSet 87299



# Replace Select™ TC PureSet 87299

INCLUDED IN THE SET

PUR0900	NobelSpeedy Groovy/Brånemark System/ Replace Select TC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
29149	Screwdriver Manual Unigrip 28 mm
32112	Direction Indicator Ø 2/Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
300352	Implant Package Organizer
71580	Radiographic Template – Replace Select TC
301382	Replace Select™ TC PureSet Wall Chart

## Implant Drivers NobRpl

	NP	RP
28 mm	36124	36126




\* Article to be disassembled for cleaning



# Replace Select™ TC PureSet

87299



Depth Probe 7-18 mm Z-shaped  
37791



Precision Drill  
36118\*



Bone Mill with Guide NobRpl  
NP Ø 4.6 36830\*  
RP Ø 5.3 36831\*



- Drill**  
mm  
 Ø 2 (#1) 32296\*  
 Ø 2,4/2,8 (#2) 32260\*  
 Ø 3 (#3) 32266\*  
 Ø 3,2 (#4) 32269\*  
 Ø 3,4 (#5) 32272\*  
 7-15 mm  
 Ø 2 (#1) 32297\*  
 Ø 2,4/2,8 (#2) 32261\*  
 Ø 3 (#3) 32267\*  
 Ø 3,2 (#4) 32270\*  
 Ø 3,4 (#5) 32273\*  
 10-18 mm  
 Ø 2 (#1) 32299\*  
 Ø 2,4/2,8 (#2) 32262\*  
 Ø 3 (#3) 32268\*  
 Ø 3,2 (#4) 32271\*  
 Ø 3,4 (#5) 32274\*




Brånemark System® Manual  
Torque Wrench Surgical  
32110\*\*




Screwdriver  
manual Unigrip 28 mm  
29149\*\*\*




Shank Extension  
K1589B.204.



Direction Indicator  
Ø 2/Ø 2.4-2.8 mm  
32112



Implant Driver NobRpl  
NP Short 36124  
NP Long 36125\*  
RP Short 36126  
RP Long 36127\*



Screw Top  
NP Ø 3.3-7-15 mm 37627\*  
RP Ø 4.7-13 mm 33072\*  
RP Ø 4.7-18 mm 33073\*

\* Article not included in this set  
 \*\* Article to be disassembled for cleaning  
 \*\*\* Screwdrivers available in other lengths

# Replace Select™ TC PureSet 87299

**NOT INCLUDED IN THE SET**  
(available for separate purchase)

## Drills

Precision	36118		
Ø 2.0	32296 (7–10 mm)	32297 (7–15 mm)	32299 (10–18 mm)
Ø 2.4/2.8	32260 (7–10 mm)	32261 (7–15 mm)	32262 (10–18 mm)
Ø 3.0	32266 (7–10 mm)	32267 (7–15 mm)	32268 (10–18 mm)
Ø 3.2	32269 (7–10 mm)	32270 (7–15 mm)	32271 (10–18 mm)
Ø 3.4	32272 (7–10 mm)	32273 (7–15 mm)	32274 (10–18 mm)



## Screw Taps

	NP	RP
Ø 3.3	37627 (7–15 mm)	
Ø 4.0		33072 (7–13 mm)
		33073 (7–18 mm)



## Bone Mill with Guide NobelReplace

	NP	RP
Ø 4.6	36830	
Ø 5.3		36831



## Implant Drivers NobelReplace

	NP	RP
37 mm	36125	36127




## Shank Extension

K1589B.204.



# Guided SURGERY



# NobelActive® Guided PureSet 87305



# NobelActive® Guided PureSet 87305

INCLUDED IN THE SET

PUR0600	NobelActive Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
33065	Handpiece Connector
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
29167	Manual Torque Wrench Adapter Prosthetic*
34584	NobelActive Manual Torque Wrench Surgical
301165	NobelActive Guided PureSet Wall Chart

## Guided Drill Guides

	NP	RP
Ø 2.0	32814	32815
Ø 2.8	32817	32818
Ø 3.2	35882	32822
Ø 3.6		35883
Ø 4.2		35884



## Guided Implant Mounts

NP	RP
35887	35888



## Guided Template Abutments

NP	RP
35890	35891




## Guided Tissue Punches

NP	RP
37153	37154



\* Article to be disassembled for cleaning

# NobelActive® Guided PureSet 87305



\* Article not included in this set

\*\* Article not included in this set and also available in other lengths

\*\*\* Screwdrivers available in other lengths

# NobelActive® Guided PureSet 87305


NOT INCLUDED IN THE SET  
(available for separate purchase)

## Guided Twist Drill

Ø 1.5 x 20 mm	33066	
---------------	-------	---

## Guided Pilot Twist Drills

Ø 1.5 (10+) 10–15 mm	37397	
----------------------	-------	---

Ø 2.0 (10+) 7–13 mm	33107	
---------------------	-------	---

Ø 2.0 (10+) 7–18 mm	32746	
---------------------	-------	---

## Guided Screw Tap

	NP	RP	WP
--	----	----	----

Ø 3.5 mm	35878 (8.5–10 mm)		
----------	-------------------	--	--

Ø 4.3 mm		35879 (8.5–10 mm)	
----------	--	-------------------	--

Ø 5.0 mm		35880 (8.5–10 mm)	
----------	--	-------------------	--

Ø 5.5 mm			37946 (7–8.5 mm) 37947 (10–15 mm)
----------	--	--	--------------------------------------



## Guided Implant Mounts

RP	WP
----	----

35889	37948
-------	-------



## Guided Tissue Punches

RP	WP
----	----

37155	38060
-------	-------



## Handle for Guided Drill Guide

32813	
-------	--

## Guided Dense Bone Screw Tap

	NP	RP	WP
--	----	----	----

Ø 3.5 mm	38179 (8.5–18 mm)		
----------	-------------------	--	--

Ø 4.3 mm		38180 (8.5–18 mm)	
----------	--	-------------------	--

Ø 5.0 mm		38181 (8.5–18 mm)	
----------	--	-------------------	--

Ø 5.5 mm			38182 (7–8.5 mm) 38183 (10–15 mm)
----------	--	--	--------------------------------------



## Guided Drill Guides

	RP	WP
--	----	----

Ø 2.0 mm	32816	32816
----------	-------	-------

Ø 2.8 mm	38046	38046
----------	-------	-------

Ø 3.2 mm	35885	35885
----------	-------	-------

Ø 3.6 mm	32825	32825
----------	-------	-------

Ø 4.2 mm	35886	35886
----------	-------	-------

Ø 5.0 mm	32826	
----------	-------	--



## Guided Twist Step Drills

	7–13 mm	7–18 mm
--	---------	---------

Ø 2.4/2.8 (10+)	35839	35844
-----------------	-------	-------

Ø 2.8/3.2 (10+)	35840	35875
-----------------	-------	-------

Ø 3.2/3.6 (10+)	35841	35876
-----------------	-------	-------

Ø 3.8/4.2 (10+)	35842	35877
-----------------	-------	-------

Ø 4.2/4.6 (10+)	35843	38045
-----------------	-------	-------

Ø 4.2/5.0 (10+)	37934	37945
-----------------	-------	-------



## Guided Start Drill

Ø 2.0 (10+)	37152
-------------	-------



## Guided Template Abutments

RP	WP
----	----

35892	37949
-------	-------



## Guided Anchor Pins

Ø 1.5 mm Short Shaft	34761
----------------------	-------







# NobelParallel™ CC Guided PureSet 87306

INCLUDED IN THE SET

PUR0700	NobelParallel CC Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
28839	Manual Torque Wrench Surgical*
29167	Manual Torque Wrench Adapter Prosthetic
33065	Handpiece Connector
301166	NobelParallel CC Guided PureSet Wall Chart

## Guided Drill Guides

	NP	RP
Ø 2.0 mm	32814	32815
Ø 2.8 mm	32817	32818
Ø 3.2 mm	35882	
Ø 3.6 mm		35883



## Guided Tissue Punches

NP	RP
38059	37154



## Guided Implant Mounts

NP	RP
38065	38066




## Guided Template Abutments

NP	RP
38069	37158



\* Article to be disassembled for cleaning


# NobelParallel™ CC Guided PureSet 87306




Guided Twist Drill  
Ø 1.5 x 20 mm  
33066\*




Guided Anchor  
Pin Ø 1.5 mm  
30909



Guided Anchor Pin  
Ø 1.5 mm short shaft  
34761\*



Guided Tissue Punch  
NP 38059  
RP 37154  
WP/6.0 (S.0) 37155\*  
WP/6.0 (S.5) 38060\*



Guided Start Drill  
37152\*



Guided Twist Step Drill  
Ø 2.4/2.8 (#3) 35844\*\*  
Ø 2.8/3.2 (#4) 35854\*\*  
Ø 3.2/3.6 (#5) 35859\*\*  
Ø 3.6/4.0 (#6) 35864\*\*  
Ø 4.0/4.4 (#7) 38045\*\*  
Ø 4.4/5.0 (#8) 37945\*\*



Prosthetic Adapter  
29167



Surgical Torque  
Wrench Manual  
28839



Guided Implant  
Mount  
NP 38065  
RP Ø 4.3 38066  
RP Ø 5.0 38067\*  
WP Ø 5.5 38068\*



Guided Template  
Abutment  
NP 38069  
RP Ø 4.3 37158  
RP Ø 5.0 37159\*  
WP Ø 5.5 38070\*



Screwdriver  
Unigrip Machine  
29151\*\*\*



Screwdriver  
Unigrip Manual  
29149\*\*\*




Handpiece Connector  
33065



Guided Drill Guide NP  
#1: Article 32817  
#2: Article 32818  
#3: Article 32817  
#4: Article 32818  
#5: Article 32815  
#6: Article 32815  
#7: Article 32815  
#8: Article 32815



Guided Drill Guide RP 4.3  
#1: Article 32816  
#2: Article 32816  
#3: Article 32816  
#4: Article 32816  
#5: Article 32816  
#6: Article 32816  
#7: Article 32816  
#8: Article 32816



Guided Drill Guide RP 5.0  
#1: Article 32816  
#2: Article 32816  
#3: Article 32816  
#4: Article 32816  
#5: Article 32816  
#6: Article 32816  
#7: Article 32816  
#8: Article 32816




Guided Cortical Drill  
NP Ø 3.75 38061\*  
RP Ø 4.3 38062\*  
RP Ø 5.0 38063\*  
WP Ø 5.5 38064\*



Guided Screw Top  
7-10 mm  
NP Ø 3.75 38047\*  
RP Ø 4.3 38048\*  
RP Ø 5.0 38051\*  
WP Ø 5.5 38054\*



Guided Screw Top  
11.5-13 mm  
NP Ø 3.75 38049\*  
RP Ø 4.3 38052\*  
RP Ø 5.0 38054\*  
WP Ø 5.5 38057\*



Guided Screw Top  
15-18 mm  
NP Ø 3.75 38049\*  
RP Ø 4.3 38052\*  
RP Ø 5.0 38054\*  
WP Ø 5.5 38058\*

\* Article not included in this set




\*\* Article not included in this set and also available in other lengths

\*\*\* Screwdrivers available in other lengths

# NobelParallel™ CC Guided PureSet 87306

NOT INCLUDED IN THE SET  
(available for separate purchase)

## Guided Twist Drills

Ø 1.5 x 20 mm	33066	
Ø 2.0 (10+) 7-13 mm	33107	
Ø 2.0 (10+) 7-18 mm	32746	

## Guided Drill Guides

	RP	WP
Ø 2.0 mm	32816	32816
Ø 2.8 mm	38046	38046
Ø 3.6 mm	35885	35885
Ø 4.2 mm	32825	32825
Ø 4.6 mm		35886
Ø 5.0 mm		32826



## Guided Implant Mounts

RP	WP
38067	38068



## Guided Template Abutments

RP	WP
37159	38070



## Guided Screw Tap

	NP	RP	WP
Ø 3.75	38047 (7-10 mm)		
	38048 (11.5-13 mm)		
	38049 (15-18 mm)		
Ø 4.3		38050 (7-10 mm)	
		38051 (11.5-13 mm)	
		38052 (15-18 mm)	
Ø 5.0		38053 (7-10 mm)	
		38054 (11.5-13 mm)	
		38055 (15-18 mm)	
Ø 5.5			38056 (7-10 mm)
			38057 (11.5-13 mm)
			38058 (11.5-15 mm)



## Handle for Guided Drill Guide

32813	
-------	--

## Guided Tissue Punches

RP	WP
37155	38060



## Guided Start Drill

Ø 2.0 (10+)	37152
-------------	-------



## Guided Anchor Pins

Ø 1.5 mm Short Shaft	34761
----------------------	-------



## Guided Cortical Drills

Ø 3.75	38061
Ø 4.3	38062
Ø 5.0	38063
Ø 5.5	38064




## Guided Twist Step Drills

	7-13 mm	7-18 mm
Ø 2.4/2.8 (10+)	35839	35844
Ø 2.8/3.2 (10+)	35840	35875
Ø 3.2/3.6 (10+)	35841	35876
Ø 3.8/4.2 (10+)	35842	35877
Ø 4.2/4.6 (10+)	35843	38045
Ø 4.2/5.0 (10+)	37934	37945



# NobelReplace® CC Guided PureSet 87307





# NobelReplace® CC Guided PureSet 87307

INCLUDED IN THE SET

PUR0800	NobelReplace CC Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
33065	Handpiece connector
28839	Manual Torque Wrench Surgical*
29167	Manual Torque Wrench Adapter Prosthetic
301167	NobelReplace CC Guided PureSet Wall Chart

## Guided Drill Guides

	NP	RP	
Ø 2.0 mm	32814	32815	
RP to NP		33018	

## Guided Dense Bone Drills Tapered

	NP	RP
13 mm	32844	32847
16 mm	32845	32848



## Guided Drills Tapered

	NP	RP
(+) 8 mm	32827	32831
(+) 10 mm	32828	32832
(+) 11.5 mm	36119	36120
(+) 13 mm	32829	32833
(+) 16 mm	32830	32834



## Guided Implant Mounts

NP	RP
37149	37150



## Guided Template Abutments

NP	RP
37157	37158



## Guided Tissue Punches

NP	RP
37153	37154



## Guided Screw Taps Tapered

NP	RP
37171	32858



\* Article to be disassembled for cleaning



# NobelReplace® CC Guided PureSet 87307

NOT INCLUDED IN THE SET  
(available for separate purchase)

## Guided Twist Drills

Ø 1.5 x 20 mm 33066



Ø 2.0 x (10+) 8–16 mm 36123



## Guided Drill Guides

RP

Ø 2.0 mm 32816



WP to NP 33019



WP to RP 33020



## Guided Drills Tapered

WP

(+) 8 mm 32835

(+) 10 mm 32836

(+) 11.5 mm 36121

(+) 13 mm 32837

(+) 16 mm 32838



## Handle for Guided Drill Guide

32813



## Irrigation Needle

2042



## Guided Template Abutments

RP

37159



## Guided Implant Mount

RP

37151



## Guided Tissue Punch

RP

37155



## Guided Screw Tap Tapered

WP

32860



## Guided Dense Bone Drills Tapered

WP

13 mm

32850

16 mm

32851



## Guided Cortical Drills

Ø 3.5 37167

Ø 4.3 37168

Ø 5.0 37169



## Guided Start Drill

Ø 2.0 (10+) 37152




## Guided Anchor Pins

Ø 1.5 mm Short Shaft 34761




# Prosthetics





# Prosthetic PureSet



# Prosthetic PureSet Basic 87301

PUR0500 Prosthetic PureSet™ Tray

301232 Wall Chart Prosthetic PureSet Basic

## Screwdriver Machine UniGrip

20 mm 29151



30 mm 29153



## Screwdriver Machine Omnigrip

20 mm 37379



30 mm 37381



## Screwdriver Machine Multi-unit

21 mm 29158



## Handle for Machine Instruments

29161



## Manual Torque Wrench Prosthetic



29165



Handle for machine instruments  
29161




Torque Wrench Prosthetic  
29165



Omnigrip machine\*  
20 mm 37379  
30 mm 37381



Unigrip machine\*  
20 mm 29151  
30 mm 29153



Multi-unit machine  
21 mm 29158

\* Screwdrivers available in other lengths

# Prosthetic PureSet Nobel Biocare N1™ 87302

PUR0500 Nobel Biocare N1™ PureSet™ Tray

301076 Wall Chart Prosthetic PureSet

## Screwdriver Machine UniGrip

30 mm 29153



## Screwdriver Machine Omnigrip Mini

28 mm 300853



## Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



## Screwdriver Machine Multi-unit

21 mm 29158



## Handle for Machine Instruments

29161



## Manual Torque Wrench Prosthetic Nobel Biocare N1™


301082



Handle for machine instruments  
29161



Torque Wrench Prosthetic NB N1  
301082



UniGrip machine\*  
30 mm 29153

Machine Omnigrip Mini  
28 mm 300853

Machine NB N1 Base  
28 mm 300903

Multi-unit machine  
21 mm 29158

\* Screwdrivers available in other lengths

# Prosthetic PureSet 87353

PUR0500 Prosthetic PureSet™ Tray

301076 Wall Chart Prosthetic PureSet

## Screwdriver Machine UniGrip

30 mm 29153



## Screwdriver Machine Omnigrip Mini

28 mm 300853



## Screwdriver Machine Omnigrip

30 mm 37381



## Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



## Screwdriver Machine Multi-unit

21 mm 29158



\* Screwdrivers available in other lengths

## Handle for Machine Instruments

29161




## Manual Torque Wrench Prosthetic

29165



## Manual Torque Wrench Prosthetic Nobel Biocare N1™

301082



## Additional prosthetic instruments

NOT INCLUDED IN THE SET  
(available for separate purchase)

### Screwdrivers Omnigrip Mini

300852	Screwdriver Machine Omnigrip™ Mini 20 mm
300854	Screwdriver Machine Omnigrip™ Mini 36 mm
300855	Screwdriver Manual Omnigrip™ Mini 20 mm
300856	Screwdriver Manual Omnigrip™ Mini 36 mm

### Screwdrivers Nobel Biocare N1™ Base

300904	Screwdriver Machine Nobel Biocare N1™ Base 36 mm
301135	Screwdriver Manual Nobel Biocare N1™ Base 28 mm
301136	Screwdriver Manual Nobel Biocare N1™ Base 36 mm

### Screwdrivers Unigrip

29148	Screwdriver Manual UniGrip 20 mm
29149	Screwdriver Manual UniGrip 28 mm
29150	Screwdriver Manual UniGrip 36 mm
29151	Screwdriver Machine UniGrip 20 mm
29152	Screwdriver Machine UniGrip 25 mm
29154	Screwdriver Machine UniGrip 35 mm

### Screwdrivers Multi-unit

29156	Screwdriver Manual Multi-unit 25 mm
29157	Screwdriver Manual Multi-unit Brånemark System WP 25 mm
29159	Screwdriver Machine Multi-unit Brånemark System WP 20 mm

### Screwdrivers Omnigrip

37376	Omnigrip Screwdriver Manual 20 mm
37377	Omnigrip Screwdriver Manual 28 mm
37378	Omnigrip Screwdriver Manual 36 mm
37379	Omnigrip Screwdriver Machine 20 mm
37380	Omnigrip Screwdriver Machine 25 mm
37382	Omnigrip Screwdriver Machine 35 mm

### Screwdrivers On1 Base

38752	On1 Screwdriver Machine 21 mm
38753	On1 Screwdriver Manual 25 mm

### Screwdrivers NobelPearl

300661	NobelPearl Screwdriver RP/WP 19 mm
300678	NobelPearl Screwdriver RP/WP 23 mm
300681	NobelPearl Screwdriver RP/WP 28 mm

### Screwdriver Cover Screw Brånemark System

DIB 097-0	Cover Screw Driver Brånemark System Hexagon
-----------	---


### Screwdrivers Locator System by Zest Dental Solutions and Ball Abutments

REF08008	.050" (1.25 mm) Hex Latch Driver Torque Wrench
REF08913	Locator Latch Type Torque Wrench Driver 23 mm
REF08914	Locator Latch Type Torque Wrench Driver 29 mm
29026	Screwdriver Manual Ball Abutment 22 mm
29027	Screwdriver Machine Ball Abutment 24 mm

### Screwdrivers Scan Body by Elos

C13485	Elos Accurate Screwdriver Manual Short 0.9 Hex
C13048	Elos Accurate Screwdriver Manual Long 0.9 Hex

# Accessories



# PureSet Tray accessories

## Spare parts

---

PUR0201 NobelActive/NobelParallel CC PureSet Plate

---

PUR0301 NobelReplace CC PureSet Plate

---

PUR0401 Nobel Biocare N1™ PureSet Plate

---

PUR0501 Prosthetic PureSet Tray Plate

---

PUR0601 NobelActive Guided PureSet Plate

---

PUR0701 NobelParallel CC Guided PureSet Plate


---

PUR0801 NobelReplace CC Guided PureSet Plate

---

PUR0901 NobelSpeedy Groovy/Brånemark System/  
Replace Select TC PureSet Plate

---



---

## ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](https://nobelbiocare.com/store)

---

## ORDER BY PHONE

Call our customer service team, or contact your sales representative.

---

## LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

[nobelbiocare.com/warranty](https://nobelbiocare.com/warranty)

[nobelbiocare.com/pureset](https://nobelbiocare.com/pureset)



88118 GB 2207 Printed in the EU © Nobel Biocare Services AG, 2022. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.