



All-on-4[®]

TREATMENT concept

with NobelParallel™
Conical Connection TiUltra™
and NobelActive® TiUltra™

nobelbiocare.com/all-on-4

Product overview



Implants

NobelParallel™ Conical Connection TiUltra™

Platform	Implant Ø	Length*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
NP 3.5	3.75 mm	300295	300296	300297	300298	300299	300300	300301
RP 4.3	4.3 mm	300302	300303	300304	300305	300306	300307	300308
RP 5.0	5.0 mm	300309	300310	300311	300312	300313	300314	300315
WP 5.5	5.5 mm	300316	300317	300318	300319	300320	300321	




NobelActive® TiUltra™

Platform	Implant Ø	Length*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
NP 3.5	3.5 mm		300244	300245	300246	300247	300248	300249
RP 4.3	4.3 mm		300250	300251	300252	300253	300254	300255
RP 5.0	5.0 mm		300256	300257	300258	300259	300260	300261
WP 5.5	5.5 mm	300262	300263	300264	300265	300266	300267	




* Please note that actual implant length is 0.5 mm shorter than indicated in name

PureSet™ Trays



300403 NobelParallel CC PureSet

All-in-one: Store NobelParallel CC and NobelActive surgical instruments in one single tray. Includes surgical instruments to perform implant surgery with NobelParallel Conical Connection implants. Drills, screw taps and bone mills must be ordered separately.



300402 NobelActive PureSet

All-in-one: Store NobelActive and NobelParallel CC surgical instruments in one single tray. Includes surgical instruments to perform implant surgery with NobelActive implants. Drills, screw taps and bone mills must be ordered separately.

Surgical components

Drills and screw taps

Precision Drill

36118



Guide Drill

35426



Twist Drills

	7-10 mm	7-15 mm	10-18 mm
Ø 2.0 mm	32296	32297	32299



Twist Step Drills

	7-10 mm	7-15 mm	10-18 mm
Ø 2.4/2.8 mm	32260	32261	32262
Ø 2.8/3.2 mm	37873	34638	34639
Ø 3.2/3.6 mm	32263	32264	32265
Ø 3.8/4.2 mm	32275	32276	32277
Ø 4.2/4.6 mm	37874	34582	34583
Ø 4.2/5.0 mm	37875	37876	



Cortical Drills NobelParallel CC

Ø 3.75 mm	38000
Ø 4.3 mm	38001
Ø 5.0 mm	38002
Ø 5.5 mm	38003



Screw Taps NobelParallel CC

	7-10 mm	7-13 mm	7-15 mm	7-18 mm
Ø 3.75 mm		37990		37991
Ø 4.3 mm		37992		37993
Ø 5.0 mm		37994		37995
Ø 5.5 mm	37996		37997	



Screw Taps NobelActive

Ø 3.5 mm	36236
Ø 4.3 mm	36237
Ø 5.0 mm	36238
Ø 5.5 7.0-10.0 mm	37871
Ø 5.5 11.5-15.0 mm	37872



Drill Kit 7-15 mm

35425

Guide Drill	35426
Twist Drill 2.0 mm	32297
Twist Step Drill 2.4/2.8 mm	32261

Surgical components for all conical connection implants

Implant Drivers

	28 mm	37 mm
NP	36718	36719
RP	36720	36721
WP	37859	37860



All-on-4 Guide

	32068
--	-------



Bone Mills and Guides

		Bone Mill with Guide	Bone Mill Guide
NP	Ø 4.4 mm	37863	37865
NP	Ø 5.2 mm	37864	37865
RP	Ø 5.2 mm	37866	37868
RP	Ø 6.2 mm	37867	37868
WP	Ø 6.7 mm	37869	37870



Prosthetic components

Multi-unit Abutments Xeal

Multi-unit Abutment Xeal¹

Height	1.5 mm	2.5 mm	3.5 mm	4.5 mm
NP	300171	300174	300177	
RP	300172	300175	300178	300180
WP	300173	300176	300179	



17° Multi-unit Abutment Xeal¹

Height	2.5 mm	3.5 mm
NP	300181	300184
RP	300182	300185
WP	300183	300186



30° Multi-unit Abutment Xeal¹

Height	3.5 mm	4.5 mm
NP	300187	300189
RP	300188	300190
WP		



Healing Cap Multi-unit⁴ (2/pkg)

Height	4.1 mm	5.5 mm
Ø 5.0	300162	300163
Ø 6.0	300164	300165
Wide	300166	300167



Multi-unit Aligning Instrument



300161

Temporary Snap Coping Multi-unit^{2,3}

38915



Temporary Coping Multi-unit^{2,4}

29046



Universal Base Multi-unit Abutment

Height 0.9 mm

301238



Consistent color coding

Prosthetic components are color coded for accurate and fast identification.

NP

Narrow
Platform

RP

Regular
Platform

WP

Wide
Platform

- 1 Clinical screw and pre-mounted holder included.
- 2 Prosthetic screw included.
- 3 Only compatible with Multi-unit Abutment Xeal.
- 4 Compatible with all Multi-unit Abutments.

Impression copings abutment level

Impression copings Abutment level⁴

	8 mm	11 mm
Closed Tray	38924	
Open Tray		29089 ⁵



Guide Pins Multi-unit⁴

	1/pkg	5/pkg
10 mm	31154	29102
20 mm	31155	29103



Intraoral scanning – Elos Position Locators

	Kit	
	IO 2C-A	IO 2C Kit



Elos Accurate IO Driver

	Long	Short
	C13048	C13485



Protection analogs and screw access guidance

Apical Drill⁴

	38853
--	-------



Drill Guide Multi-unit⁴

	38917
--	-------



Protection Analog Multi-unit⁴

	5/pkg
	29123



Abutment replicas

Abutment Replica Multi-unit⁴

	1/pkg	5/pkg
	38918	38919



Multi-unit Replica for IOS Model – Elos Analog

PMA-MUA45-1



Elos Analog Insertion tools

Elos Accurate Analog Insertion Pin	PMA-AIP-2
Elos Accurate Analog Insertion Screw	PMA-AIS-5
Elos Accurate Analog Pliers	PMA-AIT-1



⁴ Compatible with all Multi-unit Abutments.

⁵ 15 mm guide pin included.

Clinical and laboratory screws

Screw Multi-unit Angled Abutment⁴

NP	36892
RP/WP	37893



Laboratory Screw Multi-unit Angled Abutment⁴

NP	37896
RP/WP	37897



Prosthetic Screw Multi-unit Abutment⁴

NP	29285
RP/WP	29285



Laboratory Prosthetic Screw Multi-unit Abutment⁴

NP	29287 (5/pkg)
RP/WP	29287 (5/pkg)



Screwdrivers

Screwdriver Machine Multi-unit⁴

21 mm	
29158	



Screwdriver Manual Multi-unit⁴

25 mm	
29156	



Screwdriver Machine Unigrip

20 mm	25 mm	30 mm	35 mm
29151	29152	29153	29154



Screwdriver Manual Unigrip

20 mm	28 mm	36 mm
29148	29149	29150



301233 Prosthetic PureSet Basic

Unigrip Screwdriver Machine

20mm	29151
30mm	29153



Omnigrip Screwdriver Machine

20mm	37379
30mm	37381




Multi-unit Screwdriver Machine

21mm	29158
------	-------



Manual Torque Wrench Prosthetic

29165



⁴ Compatible with all Multi-unit Abutments.

**Provide optimized
esthetics**

NobelProcera® fixed restorations

Choose from a variety of precision-manufactured CAD/CAM suprastructures in titanium and zirconia and combine these with acrylic, composite or porcelain. Other options include individualized crowns on an implant bridge for maximum esthetics.

Customized implant bridges



Titanium implant bridge framework



Zirconia implant bridge framework

Customized implant bars



Hybrid bar



Montreal bar with metallic lingual

ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

ORDER BY PHONE

Call our customer service team, or contact your sales representative.

LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

nobelbiocare.com/warranty

nobelbiocare.com/all-on-4



GMT 69945 GB 2007 © Nobel Biocare Services AG, 2020. All rights reserved. CHAIRSIDE, Color and Shape of Retention Inserts, LOCATOR, LOCATOR R-Tx and R-Tx, ZEST and Zest Dental Solutions are registered trademarks of ZEST IP Holdings, LLC. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.