

NobelProcera Scan & Design Services

Your time is precious.
Outsource scan,
design and production

For abutments, implant
crowns and implant bars

Do what you do best – we'll do the rest

NobelProcera Scan and Design Services

Too much work? Looking for an easy access to leading abutments? A case with implant platforms you rarely get? New to implant bars? Whatever the reason, NobelProcera Services give you more time to focus on the work that provides the greatest return for your lab.

Enjoy peace of mind

Providing the highest quality to your customers gives you peace of mind. Moreover, NobelProcera restorations activate lifetime warranty for implants, including implants not from Nobel Biocare.*

Save time and build your business

Save time by making your abutment orders easily online. The wide selection of supported implant platforms allows you to accept almost any case and to increase your case flow.

Invest in quality, not equipment

Avoid expensive issues associated with production equipment, maintenance or stock. Instead, use our state-of-the-art production and get products with consistently high quality and precision of fit.

Getting started has never been easier

We do the work for you, so you do not need to worry about CAD/CAM decisions and up-front investments. And our team of experts are there for you if questions arise.

"Many clients specifically ask for authentic manufacturer parts, so we are thankful to have both the design and milling completely done by Nobel Biocare."

Mr Joseph M. Apap, CDT, MDT, US



"There are many advantages to outsourcing scan and design. One very important issue is that you don't have the cost of a scanner, scanner software and annual updates."

Mr Jerry Kaizer, B.S., CDT, US



All NobelProcera prosthetics are produced at NobelProcera's state-of-the-art service center in Mahwah, NJ, US.



* For details see warranty program at nobelbiocare.com/warranty.

Unrivaled prosthetics at your service

Strong NobelProcera abutments

Save time with angulated screw channel solutions and work efficiently with full-contour zirconia.

NobelProcera Implant Crowns and ASC Abutments are available on Nobel Biocare conical connection. NobelProcera abutments are also available on other implant platforms.



Enjoy the flexibility you need

NobelProcera Implant Bars can be used for fixed and fixed-removable solutions at both implant and abutment level. Select any combination from more than 170 supported implant platforms. Choose from a wide variety of bar types, attachments and retentive elements.



NobelProcera Implant Bars deliver a precision of fit superior to that of conventional casts.

Get started

Fully digital order flow for abutments

Use our simple online order form to send in your case with scan file, receive design for approval in 24 hours, and receive precision-fit products for your customers.

Place your abutment order at
nobelbiocare.com/abutmentservice

Upload your scan file

Upload your scan file

Drop scan file to attach, or click to browse

Scanner type

3Shape

CONTINUE

General information

Tooth

02

Implant Brand


Nobel Biocare

Connection type

Conical Connect

Platform

RP



BACK

CONTINUE

Fill in the prescription

Define restoration design

Cement-retained Screw-retained Full-contour

Indicate ceramic thickness (mm)

0.5

ASC feature

ASC Feature

Indicate Screw channel position

Distal-Lingual

Indicate screw angulation (0-25°)

Angulated

BACK

CONTINUE

Order process for implant bars



1. You place the order

Fabricate the master model and send the case to us.

2. We scan and design the bar for you

Our highly trained technicians will scan and design the case according to your specifications.

3. You verify the design

An e-mail with the design suggestion is sent to you for verification. If changes are needed, simply send us your feedback and we'll redesign the case.

4. We produce the bar for you

After design approval, the bar is produced with our state-of-the-art machines, verified to guarantee passive fit and polish (if applicable) so you receive it ready to use.

5. You receive the case

We ship the bar to you together with your case materials, screws and attachments (if applicable).

nobelbiocare.com

 **Nobel
Biocare™**