

Contents

Freehand surgery

Nobel Biocare N1™ PureSet 87293

NobelActive® PureSet 87294

NobelParallel™ CC PureSet 87295

NobelReplace® CC PureSet 87296

NobelSpeedy® Groovy PureSet 87297

Brånemark System® PureSet 87298

Replace Select™ TC PureSet 87299

NobelZygoma™ PureSet 88521

Guided surgery

NobelActive® Guided PureSet 87305

NobelParallel™ CC Guided PureSet 87306

NobelReplace® CC Guided PureSet 87307

Prosthetics

Prosthetic PureSet Basic 87301

Prosthetic PureSet Nobel Biocare N1™ 87302

Prosthetic PureSet 87353

Accessories

PureSet Tray accessories

PureSet™

PureSet kits for freehand and guided surgery are available for NobelActive, NobelParallel CC, NobelReplace CC and Nobel Biocare N1.

The PureSet kits NobelSpeedy Groovy, Brånemark System, Replace Select TC and NobelZygoma are available for freehand surgery.

For even greater efficiency during freehand surgery, two of Nobel Biocare's leading implant systems – NobelActive and NobelParallel CC – can both be stored in a single tray, requiring fewer instruments, thanks to their similar drilling protocols.

Ease, speed and hygiene, all-in-one with PureSet

STRAIGHTFORWARD AND INTUITIVE LAYOUTS

Designed to reflect surgical workflows and make the components easy to identify

SECURE RETAINMENT OF INSTRUMENTS


Surgical instruments are contained and secured using special metal springs

UNCOMPROMISED HYGIENE

PureSet Trays are made of stainless steel¹ and designed for automated cleaning with few manual steps


SIGNIFICANT TIME SAVINGS

Significant time savings as tray and instruments remain together during cleaning and sterilization




¹ PureSet Trays are made of stainless steel except for the PureSet plates and grommets.

Freehand SURGERY



Nobel Biocare N1™ PureSet 87293



Nobel Biocare N1™ PureSet 87293

INCLUDED IN THE SET

PUR0400 Nobel Biocare N1 PureSet Tray

301074 Radiographic Template Nobel Biocare N1

301075 Wall Chart Nobel Biocare N1 PureSet

Implant Driver Nobel Biocare N1™ TCC NP Short

28 mm 300905



Implant Driver Nobel Biocare N1™ TCC NP Long

37 mm 300906



Implant Driver Nobel Biocare N1™ TCC RP Short

28 mm 300907



Implant Driver Nobel Biocare N1™ TCC RP Long

37 mm 300908



Direction Indicator Nobel Biocare N1™

300918




OsseoShaper™ Extension Nobel Biocare N1™

300923



Depth Probe Nobel Biocare N1™

300924



Manual Torque Wrench Surgical Nobel Biocare N1™

301083



Screwdriver Machine Omnigrip™ Mini

28 mm 300853



Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



Handle for Machine Instruments

29161




Nobel Biocare N1™ PureSet 87293



Handle for Machine Instruments
29161



Manual Torque Wrench Surgical NB N1
301083



Screwdriver Machine Nobel Biocare N1 Base
28 mm 300903



OsseoShaper Extension NB N1
300923




Direction Indicator NB N1
300918



Depth Probe NB N1
300924



Screwdriver Machine Omnigrip Mini
28 mm 300853



Implant Driver NB N1 TCC NP Short
28 mm 300905



Implant Driver NB N1 TCC NP Long
37 mm 300906



Implant Driver NB N1 TCC RP Short
28 mm 300907



Implant Driver NB N1 TCC RP Short
37 mm 300908

Nobel Biocare N1™ PureSet 87293

NOT INCLUDED IN THE SET
(available for separate purchase)

Implant site preparation instruments

301077		OsseoDirector Nobel Biocare N1 Ø 1.8–2.4 x 8–14 mm
301401		OsseoDirector Nobel Biocare N1 Ø 1.8–2.4 x 8–18 mm
301081		Guided Pilot Drill Nobel Biocare N1 Ø 2.0 mm
301402		Guided Pilot Drill Nobel Biocare N1 Ø 2.0 x 8–18 mm
301078	Ø 3.5	Twist Step Drill Ø 2.5/3.4 x 10–14 mm
301421	Ø 3.5	Twist Step Drill Ø 2.5/3.4 x 10–18 mm
301079	Ø 4.0	Twist Step Drill Ø 3.3/3.8 x 8–14 mm
301422	Ø 4.0	Twist Step Drill Ø 3.3/3.8 x 8–18 mm
301080	Ø 4.8	Twist Step Drill Ø 4.2/4.4 x 8–18 mm




OsseoShaper™ 2 Nobel Biocare N1™

300892	Ø 3.5	3.5 x 9 mm
300893	Ø 3.5	3.5 x 11 mm
300894	Ø 3.5	3.5 x 13 mm
301178	Ø 3.5	3.5 x 15 mm
300895	Ø 4.0	4.0 x 7 mm
300896	Ø 4.0	4.0 x 9 mm
300897	Ø 4.0	4.0 x 11 mm
300898	Ø 4.0	4.0 x 13 mm
301179	Ø 4.0	4.0 x 15 mm
300899	Ø 4.8	4.8 x 7 mm
300900	Ø 4.8	4.8 x 9 mm
300901	Ø 4.8	4.8 x 11 mm



Bone Mill Nobel Biocare N1™ TCC

	NP	RP
Ø 4.0	300909	
Ø 5.2	300910	300910



Bone Mill Guide Nobel Biocare N1™ TCC

	NP	RP
Ø 4.0	300911	
Ø 5.2	300915	300916



Nobel Biocare N1™ PureSet 87293

NOT INCLUDED IN THE SET
(available for separate purchase)

Screwdrivers

300852	Screwdriver Machine Omnigrip Mini 20 mm
300854	Screwdriver Machine Omnigrip Mini 36 mm
300855	Screwdriver Manual Omnigrip Mini 20 mm
300856	Screwdriver Manual Omnigrip Mini 36 mm
300904	Screwdriver Machine Nobel Biocare N1 Base 36 mm
301135	Screwdriver Manual Nobel Biocare N1 Base 28 mm
301136	Screwdriver Manual Nobel Biocare N1 Base 36 mm
29158	Screwdriver Machine Multi-unit 21 mm
29156	Screwdriver Manual Multi-unit 25 mm



Torque Wrench for Prosthetic Devices

301082	Manual Torque Wrench Prosthetic Nobel Biocare N1
--------	--



NOT INCLUDED IN THE SET
(available for separate purchase)

Abutment Retrieval Tool Nobel Biocare N1™ TCC

	Short	Long
NP	301533	301521
RP	301534	301522




Abutment screw retrieval instruments

Abutment Screw Retrieval Instrument CC 3.0/NP & TCC NP/RP	301297
Abutment Screw Retrieval Reverse Drill CC 3.0/NP & TCC NP/RP	301296
Rescue Drill Guide Nobel Biocare N1 TCC NP	300919
Rescue Drill Guide Nobel Biocare N1 TCC RP	300920
Screw Tap Repair Tool Nobel Biocare N1 TCC NP	301151
Screw Tap Repair Tool Nobel Biocare N1 TCC RP	301152
N1 Base Screw Removal Tool NP/RP	301702



NobelActive® PureSet
87294

NobelParallel™ CC PureSet
87295



NobelActive® PureSet 87294

PUR0200	NobelActive/NobelParallel CC PureSet Tray
34584	NobelActive Manual Torque Wrench Surgical*
32180	Surgical Driver
32112	Direction Indicator Ø 2/Ø 2.4 – 2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37887	NobelActive Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

Implant Drivers

	3.0	NP	RP	WP
28 mm	36773	36718	36720	37859
37 mm	36774	36719	36721	37860



NobelParallel™ CC PureSet 87295

PUR0200	NobelActive/NobelParallel CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
32112	Direction Indicator Ø 2/Ø 2.4 – 2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37792	NobelParallel Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

Implant Drivers

	NP	RP	WP
28 mm	36718	36720	37859
37 mm	36719	36721	37860



INCLUDED IN THE SET


* Article to be disassembled for cleaning

NobelActive® PureSet

87294

NobelParallel™ CC PureSet

87295



Depth Probe 7-18 mm
37791



Precision Drill
36118*





Bone Mill (with Guide) CC
3.0 Ø 4.0 37861A*
NP Ø 4.4 37863*
RP Ø 5.2 37864*
WP Ø 5.2 37866*
RP Ø 6.2 37867*
WP Ø 6.7 37869*



Bone Mill Guide CC
3.0 37862-A
NP 37865*
RP 37868*
WP 37870*

* Article not included in this set
** Article to be disassembled for cleaning
*** Screwdrivers available in other lengths

A Specific to NobelActive
P Specific to NobelParallel CC



NobelActive Manual
Torque Wrench Surgical
34584A***



NobelReplace Manual
Torque Wrench Surgical
28839P**



Screw Tap^A
Ø 3.0 36816*
NP 3.5 36236*
RP 4.3 36237*
RP 5.0 36238*
WP 5.5 (7-10) 37871*
WP 5.5 (11.5-15) 37872*



Cortical Drill^P
NP 3.75 38000*
RP 4.3 38001*
RP 5.0 38002*
WP 5.5 38003*



Screw Tap^P
NP 3.75 (7-13) 37990*
NP 3.75 (7-18) 37991*
RP 4.3 (7-13) 37992*
RP 4.3 (7-18) 37993*
RP 5.0 (7-13) 37994*
RP 5.0 (7-18) 37995*
WP 5.5 (7-10) 37996*
WP 5.5 (7-15) 37997*



Twist Drill
Ø 1.5 33297-A
Ø 2 33297*



Twist Step Drill
Ø 2.4, 2.8 mm 32261*
Ø 2.8, 3.2 mm 34658*
Ø 3.6, 4.0 mm 32274*
Ø 3.6, 4.0 mm 32274*
Ø 4.2, 4.6 mm 34582*
Ø 4.2, 5.0 mm 37874*




Implant Driver 3.0
28 mm 36770
37 mm 36774



Implant Driver NP
28 mm 36719



Implant Driver RP
28 mm 36720
37 mm 36721



Implant Driver WP
28 mm 37859
37 mm 37860

NobelActive® PureSet

87294

NobelParallel™ CC PureSet

87295

Drills

Precision	36118			
Ø 1.5	31278 (7-15 mm) ^A			
Ø 2.0	32296 (7-10 mm)	32297 (7-15 mm)	32299 (10-18 mm)	
Ø 2.4 / 2.8	32260 (7-10 mm)	32261 (7-15 mm)	32262 (10-18 mm)	
Ø 2.8 / 3.2	37873 (7-10 mm)	34638 (7-15 mm)	34639 (10-18 mm)	
Ø 3.2 / 3.6	32263 (7-10 mm) ^P	32264 (7-15 mm)	32265 (10-18 mm)	
Ø 3.8 / 4.2	32275 (7-10 mm) ^P	32276 (7-15 mm)	32277 (10-18 mm)	
Ø 4.2 / 4.6	37874 (7-10 mm)	34582 (7-15 mm)	34583 (10-18 mm)	
Ø 4.2 / 5.0	37875 (7-10 mm)	37876 (7-15 mm)		



Bone Mill Guides

3.0	NP	RP	WP
37862 ^A	37865	37868	37870



Shank Extension

KI589B.204.



Screw Taps NobelActive®

	3.0	NP	RP	WP
Ø 3.0	36816			
Ø 3.5		36236		
Ø 4.3			36237	
Ø 5.0			36238	
Ø 5.5				37871 (7-10 mm) 37872 (11.5-15 mm)



Screw Taps NobelParallel™ CC

	NP	RP	WP
Ø 3.75	37990 (7-13 mm) 37991 (7-18 mm)		
Ø 4.3		37992 (7-13 mm) 37993 (7-18 mm)	
Ø 5.0		37994 (7-13 mm) 37995 (7-18 mm)	
Ø 5.5			37996 (7-10 mm) 37997 (7-15 mm)



NOT INCLUDED IN THE SET
(available for separate purchase)

Cortical Drills NobelParallel™ CC

	NP	RP	WP
Ø 3.75	38000		
Ø 4.3		38001	
Ø 5.0		38002	
Ø 5.5			38003




Bone Mills with Guide

	3.0	NP	RP	WP
Ø 4.0	37861 ^A			
Ø 4.4		37863		
Ø 5.2		37864	37866	
Ø 6.2			37867	
Ø 6.7				37869



A Specific to NobelActive only
P Specific to NobelParallel CC only

NobelReplace® CC PureSet 87296



NobelReplace® CC PureSet 87296

INCLUDED IN THE SET

PUR0300	NobelReplace CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
2042	Irrigation Needle
37320	NobelReplace Radiographic Template
300567	NobelReplace CC PureSet Wall Chart

Drills

	NP Ø 3.5	RP Ø 4.3	RP Ø 5.0
8 mm	32075	32076	32077
10 mm	29367	29370	29373
11.5 mm	36113	36114	36115
13 mm	29368	29371	29374
16 mm	29369	29372	29375



Dense Bone Drills

	NP Ø 3.5	RP Ø 4.3	RP Ø 5.0
13 mm	29377	29380	29383
16 mm	29378	29381	29384



Direction Indicators

NP	RP	WP
32255	32256	32257



Screw Taps

NP	RP	WP
36717	32090	32091



Implant Drivers


	NP	RP
28 mm	36718	36720
37 mm	36719	36721




* Article to be disassembled for cleaning


NobelReplace® CC PureSet 87296

 Precision Drill
36118*

 Drill with Tip
Tapered 2 mm
36117*

 Bone Mill (with Guide) CC
NP Ø 4.4 37863*
NP Ø 5.2 37864*
RP Ø 5.2 37866*
RP Ø 6.2 37867*

 Bone Mill Guide CC
NP 37865*
RP 37868*



* Article not included in this set

** Article to be disassembled for cleaning

*** Screwdrivers available in other lengths

 Direction Indicator
Tapered
NP 32255
RP 32256
WP 32257


 Implant Driver NP
28 mm 36718
37 mm 36719

 Implant Driver RP
28 mm 36720
37 mm 36721

 Screw Tap Tapered
NP 3.5 36717
RP 4.3 32090
WP 5.0 32091

 Drill Tapered NP 3.5
8 mm 32075
10 mm 29267
11.5 mm 36113
13 mm 29270
16 mm 29269

 Drill Tapered RP 4.3
8 mm 32076
10 mm 29270
11.5 mm 36114
13 mm 29271
16 mm 29272


 Drill Tapered WP 5
8 mm 32077
10 mm 29273
11.5 mm 36115
13 mm 29274
16 mm 29275




NobelReplace Manual
Torque Wrench Surgical
28839**

 Screwdriver
manual Unigrip 28 mm
29149***

 Shank Extension
K1589B.204.

 Dense Bone Drill
Tapered 13 mm
NP 3.5 29377
RP 4.3 29380
WP 5.0 29383

 Dense Bone Drill
Tapered 16 mm
NP 3.5 29378
RP 4.3 29381
WP 5.0 29384

NobelReplace® CC PureSet 87296

NOT INCLUDED IN THE SET
(available for separate purchase)

Bone Mills with Guide

	NP	RP
Ø 4.4	37863	
Ø 5.2	37864	37866
Ø 6.2		37867



Bone Mill Guides

NP	RP
37865	37868



Drills

36117	Drill with Tip Tapered 2 mm	
36118	Precision Drill	

Shank Extension


KI589B.204.

NobelSpeedy® Groovy PureSet

87297

Brånemark System® PureSet

87298



NobelSpeedy® Groovy PureSet

87297

Brånemark System® PureSet

87298

INCLUDED IN THE SET

PUR0900	NobelSpeedy Groovy/Brånemark System/ Replace Select TC PureSet Tray
DIB 097-0	Cover Screw Driver Brånemark System Hexagon
32110	Brånemark System Manual Torque Wrench Surgical*
32112	Direction Indicator Ø 2 / Ø 2.4–2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37628	NobelSpeedy Groovy Radiographic Template 7–25 mm**
33097	Brånemark System/NobelSpeedy Radiographic Template***
301381	NobelSpeedy Groovy/Brånemark System PureSet Wall Chart

Implant Drivers

	NP	RP	WP
21 mm		29129	29134**
26 mm	29126**	29130	29135**
34 mm	29127**		



Implant Driver Wrench


	NP	RP	WP
12 mm		29132	29136**
21 mm	29128**		



* Article to be disassembled for cleaning
 ** Article only included in NobelSpeedy Groovy PureSet
 *** Article only included in Brånemark System PureSet

NobelSpeedy® Groovy PureSet 87297

Brånemark System® PureSet 87298



NobelSpeedy Depth Probe 18-25 mm
37692**

Depth Probe 7-18 mm Z-shaped
37791



Precision Drill
36118*



Bone Mill (with Guide)
Bmk Syst
NP Ø 4.5 33392*
RP Ø 5.1 33393*
WP Ø 6.5 33495*



Bone Mill Guide
Bmk Syst
NP 33496*
RP 33497*
WP 33498*



Direction Indicator
Ø 2/Ø 2.4*
32112



Screwdriver
manual Unigrip 28 mm
29149**



Cover Screw Driver
Bmk Syst Hexagon
DIB 097-O




Brånemark System Manual
Torque Wrench Surgical
32110**

Brånemark System
Counterbore
NP 32281*
RP 32283*
WP 32285*

Shank Extension
KI589B.204.

Screw Tap
NP Ø 3.75-7.15 mm 37627*
RP Ø 3.75-7.18 mm 32290**
RP Ø 4.7-13 mm 33072*
RP Ø 4.7-18 mm 33073*
WP Ø 5.7-13 mm 32292*
WP Ø 5.7-18 mm 32293*
WP Ø 6.7-13 mm 32294**
WP Ø 6.7-18 mm 32295**




Drill
7-10 mm
Ø 2 (#1) 32266*
Ø 2.4 (#2) 32267*
Ø 2.8 (#3) 32268*
Ø 3.2 (#4) 32269*
Ø 3.6 (#5) 32270*
Ø 3.7/3.6 (#6) 32263*
Ø 3.8/4.2 (#7) 32275**
Ø 5 (#8) 32278**

7-15 mm
Ø 2 (#1) 32297*
Ø 2.4/2.8 (#2) 32261*
Ø 3 (#3) 32267*
Ø 3.2 (#4) 32270*
Ø 3.4 (#5) 32273*
Ø 3.6/4.2 (#7) 32274*
Ø 5 (#8) 32279**

10-18 mm
Ø 2 (#1) 32299*
Ø 2.4/2.8 (#2) 32262*
Ø 3 (#3) 32271*
Ø 3.4 (#4) 32274*
Ø 3.6/4.2 (#7) 32265*
Ø 5 (#8) 32280**

18-25 mm
Ø 2.0 (#1) 37614**
Ø 2.4 (#2) 37615**
Ø 2.8 (#3) 37616**
Ø 3.2 (#4) 37617**
Ø 3.6 (#5) 37618**



Implant Driver Bmk Syst
NP 26 mm 29126
RP 26 mm 29127
RP 34 mm 29129
RP 26 mm 29130
RP 34 mm 29131
WP 21 mm 29134
WP 26 mm 29135



Implant Driver Wrench
Adbp Bmk Syst
NP 12 mm 29132
RP 12 mm 29133
RP 21 mm 29136
WP 12 mm 29136

* Article not included in this set
** Article to be disassembled for cleaning
*** Screwdrivers available in other lengths

B Specific to Brånemark System
S Specific to NobelSpeedy Groovy

NobelSpeedy® Groovy PureSet

87297

Brånemark System® PureSet

87298

NOT INCLUDED IN THE SET
(available for separate purchase)

Drills

Precision	36118				
Ø 2.0	32296 (7-10 mm)	32297 (7-15 mm)	32299 (10-18 mm)	37614 (18-25 mm) ^S	37619 (18-25 mm) LS ^S
Ø 2.4/2.8	32260 (7-10 mm)	32261 (7-15 mm)	32262 (10-18 mm)	37618 (18-25 mm) ^S	37620 (18-25 mm) LS* ^S
Ø 3.0	32266 (7-10 mm)	32267 (7-15 mm)	32268 (10-18 mm)		
Ø 3.2	32269 (7-10 mm)	32270 (7-15 mm)	32271 (10-18 mm)	37616 (18-25 mm) ^S	
Ø 3.4	32272 (7-10 mm)	32273 (7-15 mm)	32274 (10-18 mm)	37617 (18-25 mm) ^S	
Ø 3.2/3.6	32263 (7-10 mm)	32264 (7-15 mm)	32265 (10-18 mm)		
Ø 3.8/4.2	32275 (7-10 mm)	32276 (7-15 mm)	32277 (10-18 mm)		
Ø 5.0	32278 (7-10 mm) ^S	32279 (7-15 mm) ^S	32280 (10-18 mm) ^S		



Bone Mills (with Guide) Bmk Syst

	NP	RP	WP
Ø 4.5	33392		
Ø 5.1		33393	
Ø 6.5			33495



Bone Mill Guides Brånemark System®

NP	RP	WP
33496	33497	33498




Brånemark System® Counterbore

NP	RP	WP
32281	32283	32285



NobelSpeedy® Depth Probe 18-25 mm

37692^S



B Specific to Brånemark System
S Specific to NobelSpeedy Groovy
* Article to be removed prior to placing lid on tray

NobelSpeedy® Groovy PureSet

87297

Brånemark System® PureSet

87298

NOT INCLUDED IN THE SET
(available for separate purchase)

Implant Drivers

	NP	RP	WP
21 mm			29134*
26 mm	29126*		29135*
34 mm	29127*	29131	



Implant Driver Wrench

	NP	RP	WP
12 mm			29136*
21 mm	29128*	29133	




Shank Extension

KI589B.204.



* Article only included in NobelSpeedy Groovy PureSet

Replace Select™ TC PureSet 87299



Replace Select™ TC PureSet 87299

INCLUDED IN THE SET

PUR0900	NobelSpeedy Groovy/Brånemark System/ Replace Select TC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical*
29149	Screwdriver Manual Unigrip 28 mm
32112	Direction Indicator Ø 2 / Ø 2.4 – 2.8 mm x 4
37791	Depth Probe 7–18 mm Z-shaped
300352	Implant Package Organizer
71580	Radiographic Template – Replace Select TC
301382	Replace Select TC PureSet Wall Chart


Implant Drivers NobelReplace®

	NP	RP
28 mm	36124	36126



* Article to be disassembled for cleaning

Replace Select™ TC PureSet 87299



Depth Probe 7-18 mm Z-shaped
37791



Precision Drill
36118*



Bone Mill with Guide NobRpl
NP Ø 4.6 36830*
RP Ø 5.3 36831*



Drill bits
 Ø 2 (H1) 32296*
 Ø 2,4/2,8 (H2) 32260*
 Ø 3 (H3) 32266*
 Ø 3,2 (H4) 32269*
 Ø 3,4 (H5) 32272*
7-15 mm
 Ø 2 (H1) 32297*
 Ø 2,4/2,8 (H2) 32261*
 Ø 3 (H3) 32267*
 Ø 3,2 (H4) 32270*
 Ø 3,4 (H5) 32273*
10-18 mm
 Ø 2 (H1) 32299*
 Ø 2,4/2,8 (H2) 32262*
 Ø 3 (H3) 32268*
 Ø 3,2 (H4) 32271*
 Ø 3,4 (H5) 32274*




Brånemark System Manual
Torque Wrench Surgical
32110**



Screwdriver
manual Unigrp 28 mm
29149***



Shank Extension
K1589B.204.



Direction Indicator
Ø 2/Ø 2.4-2.8 mm
32112



Implant Driver NobRpl
NP Short 36124
NP Long 36125*
RP Short 36126
RP Long 36127*



Screw Tap
NP Ø 3.3 37627*
RP Ø 4.7-13 mm 33072*
RP Ø 4.7-18 mm 33073*

* Article not included in this set

** Article to be disassembled for cleaning

*** Screwdrivers available in other lengths

Replace Select™ TC PureSet 87299

NOT INCLUDED IN THE SET
(available for separate purchase)

Drills

Precision	36118		
Ø 2.0	32296 (7–10 mm)	32297 (7–15 mm)	32299 (10–18 mm)
Ø 2.4 / 2.8	32260 (7–10 mm)	32261 (7–15 mm)	32262 (10–18 mm)
Ø 3.0	32266 (7–10 mm)	32267 (7–15 mm)	32268 (10–18 mm)
Ø 3.2	32269 (7–10 mm)	32270 (7–15 mm)	32271 (10–18 mm)
Ø 3.4	32272 (7–10 mm)	32273 (7–15 mm)	32274 (10–18 mm)



Bone Mill with Guide NobelReplace®

	NP	RP
Ø 4.6	36830	
Ø 5.3		36831



Implant Drivers NobelReplace®

	NP	RP
37 mm	36125	36127



Screw Taps

	NP	RP
Ø 3.3	37627 (7–15 mm)	
Ø 4.0		33072 (7–13 mm)
		33073 (7–18 mm)



Shank Extension

K1589B.204. 

NobelZygoma™ PureSet 88521



NobelZygoma™ PureSet 88521

INCLUDED IN THE SET


General

PUR1000	NobelZygoma PureSet Tray
301711	NobelZygoma PureSet Wallchart
37786	Zygoma Handle
29165	Manual Torque Wrench Prosthetic*



Depth Indicator

37789	Zygoma Depth Indicator Straight
37790	Zygoma Depth Indicator Angled




Drivers

29081	Connection to Handpiece
29149	Screwdriver Manual Unigrip 28 mm
29152	Screwdriver Machine Unigrip 25 mm
DIB 097-0	Cover Screwdriver Brånemark System Hexagon
29158	Screwdriver Machine Multi-unit 21 mm
29129	Implant Driver Brånemark System RP 21 mm
29132	Impl Driver Wrench Adap Brånemark System RP 12 mm









Drill Guard

37787	Zygoma Drill Guard
37788	Zygoma Drill Guard Short



* Article to be disassembled for cleaning

NobelZygoma™ PureSet 88521


















- * Article not included in this set
- ** Article to be disassembled for cleaning
- *** Screwdrivers available in other lengths

NobelZygoma™ PureSet 88521

NOT INCLUDED IN THE SET
(available for separate purchase)


Drills

DIA 578-0	Brånemark System Zygoma Round Bur	
32630	Brånemark System Zygoma Pilot Drill 3.5 mm	
32791	Brånemark System Zygoma Pilot Drill 3.5 mm Short	
32628	Brånemark System Zygoma Twist Drill 2.9 mm	
37766	NobelZygoma 0° Twist Drill 2.9 mm	
32629	Brånemark System Zygoma Twist Drill 2.9 mm Short	
37767	NobelZygoma 0° Twist Drill 2.9 mm Short	
32631	Brånemark System Zygoma Twist Drill 3.5 mm	
37768	NobelZygoma 0° Twist Drill 3.5 mm	
32632	Brånemark System Zygoma Twist Drill 3.5 mm Short	
37769	NobelZygoma 0° Twist Drill 3.5 mm Short	
37770	NobelZygoma 0° Twist Drill 4.0 mm	
37771	NobelZygoma 0° Twist Drill 4.0 mm Short	
37772	NobelZygoma 0° Twist Drill 4.4 mm	
37773	NobelZygoma 0° Twist Drill 4.4 mm Short	

Guided

SURGERY

NobelActive® Guided PureSet 87305



NobelActive® Guided PureSet 87305

INCLUDED IN THE SET

PUR0600	NobelActive Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
33065	Handpiece Connector
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
29167	Manual Torque Wrench Adapter Prosthetic*
34584	NobelActive Manual Torque Wrench Surgical
301165	NobelActive Guided PureSet Wall Chart

Guided Drill Guides

	NP	RP
Ø 2.0	32814	32815
Ø 2.8	32817	32818
Ø 3.2	35882	32822
Ø 3.6		35883
Ø 4.2		35884



Guided Implant Mounts

NP	RP
35887	35888



Guided Template Abutments

NP	RP
35890	35891




Guided Tissue Punches

NP	RP
37153	37154



* Article to be disassembled for cleaning

NobelActive® Guided PureSet 87305



* Article not included in this set

** Article not included in this set and also available in other lengths

*** Screwdrivers available in other lengths




NobelActive® Guided PureSet 87305

NOT INCLUDED IN THE SET
(available for separate purchase)

Guided Twist Drill

Ø 1.5 x 20 mm	33066	
---------------	-------	---

Guided Pilot Twist Drills

Ø 1.5 (10+) 10-15 mm	37397	
Ø 2.0 (10+) 7-13 mm	33107	
Ø 2.0 (10+) 7-18 mm	32746	

Guided Screw Tap

	NP	RP	WP
Ø 3.5 mm	35878 (8.5-10 mm)		
Ø 4.3 mm		35879 (8.5-10 mm)	
Ø 5.0 mm		35880 (8.5-10 mm)	
Ø 5.5 mm			37946 (7-8.5 mm) 37947 (10-15 mm)



Guided Implant Mounts

RP	WP
35889	37948



Guided Tissue Punches

RP	WP
37155	38060



Handle for Guided Drill Guide

32813



Guided Dense Bone Screw Tap

	NP	RP	WP
Ø 3.5 mm	38179 (8.5-18 mm)		
Ø 4.3 mm		38180 (8.5-18 mm)	
Ø 5.0 mm		38181 (8.5-18 mm)	
Ø 5.5 mm			38182 (7-8.5 mm) 38183 (10-15 mm)



Guided Drill Guides

	RP	WP
Ø 2.0 mm	32816	32816
Ø 2.8 mm	38046	38046
Ø 3.2 mm	35885	35885
Ø 3.6 mm	32825	32825
Ø 4.2 mm	35886	35886
Ø 5.0 mm	32826	



Guided Twist Step Drills

	7-13 mm	7-18 mm
Ø 2.4/2.8 (10+)	35839	35844
Ø 2.8/3.2 (10+)	35840	35875
Ø 3.2/3.6 (10+)	35841	35876
Ø 3.8/4.2 (10+)	35842	35877
Ø 4.2/4.6 (10+)	35843	38045
Ø 4.2/5.0 (10+)	37934	37945



Guided Start Drill

Ø 2.0 (10+)	37152
-------------	-------



Guided Template Abutments

RP	WP
35892	37949




Guided Anchor Pins

Ø 1.5 mm Short Shaft	34761
----------------------	-------



NobelParallel™ CC Guided PureSet 87306



NobelParallel™ CC Guided PureSet 87306

INCLUDED IN THE SET

PUR0700	NobelParallel CC Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
28839	Manual Torque Wrench Surgical*
29167	Manual Torque Wrench Adapter Prosthetic
33065	Handpiece Connector
301166	NobelParallel CC Guided PureSet Wall Chart

Guided Drill Guides

	NP	RP
Ø 2.0 mm	32814	32815
Ø 2.8 mm	32817	32818
Ø 3.2 mm	35882	
Ø 3.6 mm		35883



Guided Tissue Punches

NP	RP
38059	37154



Guided Implant Mounts

NP	RP
38065	38066




Guided Template Abutments

NP	RP
38069	37158



* Article to be disassembled for cleaning

NobelParallel™ CC Guided PureSet 87306



* Article not included in this set




** Article not included in this set and also available in other lengths

*** Screwdrivers available in other lengths

NobelParallel™ CC Guided PureSet 87306

NOT INCLUDED IN THE SET
(available for separate purchase)

Guided Twist Drills

Ø 1.5 x 20 mm	33066	
Ø 2.0 (10+) 7-13 mm	33107	
Ø 2.0 (10+) 7-18 mm	32746	

Guided Drill Guides

	RP	WP
Ø 2.0 mm	32816	32816
Ø 2.8 mm	38046	38046
Ø 3.6 mm	35885	35885
Ø 4.2 mm	32825	32825
Ø 4.6 mm		35886
Ø 5.0 mm		32826



Guided Implant Mounts

RP	WP
38067	38068



Guided Template Abutments

RP	WP
37159	38070



Guided Screw Tap

	NP	RP	WP
Ø 3.75	38047 (7-10 mm) 38048 (11.5-13 mm) 38049 (15-18 mm)		
Ø 4.3		38050 (7-10 mm) 38051 (11.5-13 mm) 38052 (15-18 mm)	
Ø 5.0		38053 (7-10 mm) 38054 (11.5-13 mm) 38055 (15-18 mm)	
Ø 5.5			38056 (7-10 mm) 38057 (11.5-13 mm) 38058 (11.5-15 mm)



Handle for Guided Drill Guide

32813	
-------	--

Guided Tissue Punches

RP	WP
37155	38060



Guided Start Drill

Ø 2.0 (10+)	37152
-------------	-------



Guided Anchor Pins

Ø 1.5 mm Short Shaft	34761
----------------------	-------



Guided Cortical Drills

Ø 3.75	38061
Ø 4.3	38062
Ø 5.0	38063
Ø 5.5	38064




Guided Twist Step Drills

	7-13 mm	7-18 mm
Ø 2.4 / 2.8 (10+)	35839	35844
Ø 2.8 / 3.2 (10+)	35840	35875
Ø 3.2 / 3.6 (10+)	35841	35876
Ø 3.8 / 4.2 (10+)	35842	35877
Ø 4.2 / 4.6 (10+)	35843	38045
Ø 4.2 / 5.0 (10+)	37934	37945



NobelReplace® CC Guided PureSet 87307




NobelReplace® CC Guided PureSet 87307

INCLUDED IN THE SET


PUR0800	NobelReplace CC Guided PureSet Tray
30909	Guided Anchor Pin Ø 1.5 mm
29149	Screwdriver Manual Unigrip 28 mm
29151	Screwdriver Machine Unigrip 20 mm
33065	Handpiece connector
28839	Manual Torque Wrench Surgical*
29167	Manual Torque Wrench Adapter Prosthetic
301167	NobelReplace CC Guided PureSet Wall Chart

Guided Drill Guides

	NP	RP	
Ø 2.0 mm	32814	32815	
RP to NP		33018	


Guided Dense Bone Drills Tapered

	NP	RP
13 mm	32844	32847
16 mm	32845	32848



Guided Drills Tapered

	NP	RP
(+) 8 mm	32827	32831
(+) 10 mm	32828	32832
(+) 11.5 mm	36119	36120
(+) 13 mm	32829	32833
(+) 16 mm	32830	32834



Guided Implant Mounts

NP	RP
37149	37150



Guided Template Abutments

NP	RP
37157	37158



Guided Tissue Punches

NP	RP
37153	37154




Guided Screw Taps Tapered

NP	RP
37171	32858



* Article to be disassembled for cleaning

NobelReplace® CC Guided PureSet 87307



* Article not included in this set

** Article not included in this set and also available in other lengths

*** Screwdrivers available in other lengths

NobelReplace® CC Guided PureSet 87307

NOT INCLUDED IN THE SET
(available for separate purchase)

Guided Twist Drills

Ø 1.5 x 20 mm 33066



Ø 2.0 x (10+) 8-16 mm 36123



Guided Drill Guides

RP

Ø 2.0 mm 32816



WP to NP 33019



WP to RP 33020



Guided Drills Tapered

WP


(+) 8 mm 32835

(+) 10 mm 32836

(+) 11.5 mm 36121

(+) 13 mm 32837

(+) 16 mm 32838



Handle for Guided Drill Guide

32813



Irrigation Needle


2042



Guided Template Abutments

RP


37159



Guided Implant Mount

RP

37151



Guided Tissue Punch

RP

37155



Guided Screw Tap Tapered

WP

32860




Guided Dense Bone Drills Tapered

WP

13 mm 32850

16 mm 32851



Guided Cortical Drills

Ø 3.5 37167

Ø 4.3 37168

Ø 5.0 37169



Guided Start Drill

Ø 2.0 (10+) 37152




Guided Anchor Pins

Ø 1.5 mm Short Shaft 34761



Prosthetics



Prosthetic PureSet Basic 87301

PUR0500 Prosthetic PureSet Tray

301232 Wall Chart Prosthetic PureSet Basic

Screwdriver Machine Unigrip™

20 mm 29151



30 mm 29153



Screwdriver Machine Omnigrip™

20 mm 37379



30 mm 37381



Screwdriver Machine Multi-unit

21 mm 29158



Handle for Machine Instruments

29161



Manual Torque Wrench Prosthetic



29165



Handle for Machine Instruments
29161



Torque Wrench Prosthetic
29165



Omnigrip machine*
20 mm 37379
30 mm 37381



Unigrip machine*
20 mm 29151
30 mm 29153



Multi-unit machine
21 mm 29158

* Screwdrivers available in other lengths

Prosthetic PureSet Nobel Biocare N1™ 87302

PUR0500 Nobel Biocare N1 PureSet Tray

301076 Wall Chart Prosthetic PureSet

Screwdriver Machine Unigrip™

30 mm 29153



Screwdriver Machine Omnigrip™ Mini

28 mm 300853



Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



Screwdriver Machine Multi-unit


21 mm 29158



* Screwdrivers available in other lengths

Handle for Machine Instruments

29161



Manual Torque Wrench Prosthetic Nobel Biocare N1™


301082



Handle for Machine Instruments
29161



Torque Wrench Prosthetic NB N1
301082



Unigrip machine*
30 mm 29153

Machine Omnigrip Mini
28 mm 300853

Machine NB N1 Base
28 mm 300903

Multi-unit machine
21 mm 29158

Prosthetic PureSet 87353

PUR0500 Prosthetic PureSet Tray

301076 Wall Chart Prosthetic PureSet

Screwdriver Machine Unigrip™

30 mm 29153



Screwdriver Machine Omnigrip™ Mini

28 mm 300853



Screwdriver Machine Omnigrip™

30 mm 37381



Screwdriver Machine Nobel Biocare N1™ Base

28 mm 300903



Screwdriver Machine Multi-unit

21 mm 29158



* Screwdrivers available in other lengths

Handle for Machine Instruments

29161



Manual Torque Wrench Prosthetic

29165



Manual Torque Wrench Prosthetic Nobel Biocare N1™


301082



Handle for Machine Instruments
29161



Torque Wrench Prosthetic
29165



Unigrip machine*
30 mm 29153

Machine Omnigrip Mini
28 mm 300853

Omnigrip machine*
30 mm 37381

Machine NB N1 Base
28 mm 300903

Multi-unit machine
21 mm 29158

Torque Wrench Prosthetic: NB N1
301082

Additional prosthetic instruments

NOT INCLUDED IN THE SET
(available for separate purchase)

Screwdrivers Omnigrip™ Mini

300852	Screwdriver Machine Omnigrip Mini 20 mm
300854	Screwdriver Machine Omnigrip Mini 36 mm
300855	Screwdriver Manual Omnigrip Mini 20 mm
300856	Screwdriver Manual Omnigrip Mini 36 mm

Screwdrivers Nobel Biocare N1™ Base

300904	Screwdriver Machine Nobel Biocare N1 Base 36 mm
301135	Screwdriver Manual Nobel Biocare N1 Base 28 mm
301136	Screwdriver Manual Nobel Biocare N1 Base 36 mm

Screwdrivers Unigrip™

29148	Screwdriver Manual Unigrip 20 mm
29149	Screwdriver Manual Unigrip 28 mm
29150	Screwdriver Manual Unigrip 36 mm
29151	Screwdriver Machine Unigrip 20 mm
29152	Screwdriver Machine Unigrip 25 mm
29154	Screwdriver Machine Unigrip 35 mm

Screwdrivers Multi-unit

29156	Screwdriver Manual Multi-unit 25 mm
29157	Screwdriver Manual Multi-unit Brånemark System WP 25 mm
29159	Screwdriver Machine Multi-unit Brånemark System WP 20 mm

Screwdrivers Omnigrip™

37376	Omnigrip Screwdriver Manual 20 mm
37377	Omnigrip Screwdriver Manual 28 mm
37378	Omnigrip Screwdriver Manual 36 mm
37379	Omnigrip Screwdriver Machine 20 mm
37380	Omnigrip Screwdriver Machine 25 mm
37382	Omnigrip Screwdriver Machine 35 mm

Screwdrivers On1™ Base

38752	On1 Screwdriver Machine 21 mm
38753	On1 Screwdriver Manual 25 mm

Screwdrivers NobelPearl™

300661	NobelPearl Screwdriver RP/WP 19 mm
300678	NobelPearl Screwdriver RP/WP 23 mm
300681	NobelPearl Screwdriver RP/WP 28 mm

Screwdriver Cover Screw Brånemark System®

DIB 097-0	Cover Screw Driver Brånemark System Hexagon
-----------	---


Screwdrivers Locator® System by Zest Dental Solutions and Ball Abutments

REF08008	.050" (1.25 mm) Hex Latch Driver Torque Wrench
REF08913	Locator Latch Type Torque Wrench Driver 23 mm
REF08914	Locator Latch Type Torque Wrench Driver 29 mm
29026	Screwdriver Manual Ball Abutment 22 mm
29027	Screwdriver Machine Ball Abutment 24 mm

Screwdrivers Scan Body by Elos Accurate™

C13485	Elos Accurate Screwdriver Manual Short 0.9 Hex
C13048	Elos Accurate Screwdriver Manual Long 0.9 Hex

Accessories



PureSet Tray accessories

Plates

PUR0201 NobelActive/NobelParallel CC PureSet Plate

PUR0301 NobelReplace CC PureSet Plate

PUR0401 Nobel Biocare N1 PureSet Plate

PUR0501 Prosthetic PureSet Tray Plate

PUR0601 NobelActive Guided PureSet Plate

PUR0701 NobelParallel CC Guided PureSet Plate

PUR0801 NobelReplace CC Guided PureSet Plate

PUR0901 NobelSpeedy Groovy/Brånemark System/
Replace Select TC PureSet Plate

PUR1001 NobelZygoma PureSet Plate



ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

store.nobelbiocare.com

ORDER BY PHONE

Contact your sales representative, or call our customer service team.

800 322 5001

LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

nobelbiocare.com/warranty

nobelbiocare.com/pureset



MKT-5423 Rev 01 (10/23, GMT 87124) © Nobel Biocare USA, LLC, 2023. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.