



All-on-4®

TREATMENT concept

with NobelParallel™
Conical Connection TiUltra™
and NobelActive® TiUltra™

nobelbiocare.com/all-on-4

Product overview



Implants

NobelParallel Conical Connection TiUltra

Platform	Implant Ø	Length*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
NP 3.5	3.75 mm	300295	300296	300297	300298	300299	300300	300301
RP 4.3	4.3 mm	300302	300303	300304	300305	300306	300307	300308
RP 5.0	5.0 mm	300309	300310	300311	300312	300313	300314	300315
WP 5.5	5.5 mm	300316	300317	300318	300319	300320	300321	



NobelActive TiUltra

Platform	Implant Ø	Length*						
		7 mm	8.5 mm	10 mm	11.5 mm	13 mm	15 mm	18 mm
NP 3.5	3.5 mm		300244	300245	300246	300247	300248	300249
RP 4.3	4.3 mm		300250	300251	300252	300253	300254	300255
RP 5.0	5.0 mm		300256	300257	300258	300259	300260	300261
WP 5.5	5.5 mm	300262	300263	300264	300265	300266	300267	



* Please note that actual implant length is 0.5 mm shorter than indicated in name


PureSet™ Trays

NobelActive PureSet 87294

PUR0200	NobelActive/NobelParallel CC PureSet Tray
34584	NobelActive Manual Torque Wrench Surgical
32180	Surgical Driver
32112	Direction Indicator Ø 2 / Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip™ 28 mm
300352	Implant Package Organizer
37887	NobelActive Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

Implant Drivers

	28 mm	37 mm
NP	36718	36719
RP	36720	36721
WP	37859	37860



All-in-one: Store NobelActive and NobelParallel CC surgical instruments in one single tray.
Includes surgical instruments to perform implant surgery with NobelActive implants.
Drills, screw taps and bone mills must be ordered separately.


PureSet Trays

NobelParallel CC PureSet 87295

PUR0200	NobelActive/NobelParallel CC PureSet Tray
28839	NobelReplace Manual Torque Wrench Surgical
32112	Direction Indicator Ø 2 / Ø 2.4–2.8 mm × 4
37791	Depth Probe 7–18 mm Z-shaped
29149	Screwdriver Manual Unigrip 28 mm
300352	Implant Package Organizer
37792	NobelParallel Radiographic Template
300781	NobelActive/NobelParallel CC PureSet Tray Wall Chart

Implant Drivers

	28 mm	37 mm
NP	36718	36719
RP	36720	36721
WP	37859	37860



All-in-one: Store NobelParallel CC and NobelActive surgical instruments in one single tray. Includes surgical instruments to perform implant surgery with NobelParallel Conical Connection implants. Drills, screw taps and bone mills must be ordered separately.

Surgical components

Drills and screw taps

Precision Drill

36118



Guide Drill

35426



Twist Drills

	7-10 mm	7-15 mm	10-18 mm
Ø 2.0 mm	32296	32297	32299



Twist Step Drills

	7-10 mm	7-15 mm	10-18 mm
Ø 2.4/2.8 mm	32260	32261	32262
Ø 2.8/3.2 mm	37873	34638	34639
Ø 3.2/3.6 mm	32263	32264	32265
Ø 3.8/4.2 mm	32275	32276	32277
Ø 4.2/4.6 mm	37874	34582	34583
Ø 4.2/5.0 mm	37875	37876	



Cortical Drills NobelParallel CC

Ø 3.75 mm	38000
Ø 4.3 mm	38001
Ø 5.0 mm	38002
Ø 5.5 mm	38003



Screw Taps NobelParallel CC

	7-10 mm	7-13 mm	7-15 mm	7-18 mm
Ø 3.75 mm		37990		37991
Ø 4.3 mm		37992		37993
Ø 5.0 mm		37994		37995
Ø 5.5 mm	37996		37997	



Screw Taps NobelActive

Ø 3.5 mm	36236
Ø 4.3 mm	36237
Ø 5.0 mm	36238
Ø 5.5 7.0-10.0 mm	37871
Ø 5.5 11.5-15.0 mm	37872



Drill Kit 7-15 mm

35425

Guide Drill	35426
Twist Drill 2.0 mm	32297
Twist Step Drill 2.4/2.8 mm	32261

Surgical components for all conical connection implants

All-on-4 Guide

32068



Bone Mills and Guides

		Bone Mill with Guide	Bone Mill Guide
NP	Ø 4.4 mm	37863	37865
NP	Ø 5.2 mm	37864	37865
RP	Ø 5.2 mm	37866	37868
RP	Ø 6.2 mm	37867	37868
WP	Ø 6.7 mm	37869	37870



Prosthetic components

Multi-unit Abutments Xeal™

Multi-unit Abutment Xeal¹

Height	1.5 mm	2.5 mm	3.5 mm	4.5 mm
NP	300171	300174	300177	
RP	300172	300175	300178	300180
WP	300173	300176	300179	



17° Multi-unit Abutment Xeal¹

Height	2.5 mm	3.5 mm
NP	300181	300184
RP	300182	300185
WP	300183	300186



30° Multi-unit Abutment Xeal¹

Height	3.5 mm	4.5 mm
NP	300187	300189
RP	300188	300190
WP		



Temporary Snap Coping Multi-unit^{2,3}

38915



Temporary Coping Multi-unit^{2,4}

29046



Universal Base Multi-unit Abutment²

Height 0.9 mm	301238
---------------	--------

Height 0.3 mm	301365
---------------	--------



Healing Cap Multi-unit⁴ (2/pkg)

Height	4.1 mm	5.5 mm
Ø 5.0	300162	300163
Ø 6.0	300164	300165
Wide	300166	300167



Multi-unit Aligning Instrument

300161



Consistent color coding

Prosthetic components are color coded for accurate and fast identification.

NP

Narrow
Platform

RP

Regular
Platform

WP

Wide
Platform

¹ Clinical screw and pre-mounted holder included.

² Prosthetic screw included.

³ Only compatible with Multi-unit Abutment Xeal.

⁴ Compatible with all Multi-unit Abutments.

Impression copings abutment level

Impression copings Abutment level⁴

	8 mm	11 mm
Closed Tray	38924	
Open Tray		29089 ⁵



Guide Pins Multi-unit⁴

	1/pkg	5/pkg
10 mm	31154	29102
20 mm	31155	29103



Intraoral and desktop scanning – Elos Accurate™ Position Locators

	Kit
	IO 2C-A IO 2C Kit



Elos Accurate IO Driver

	Long	Short	Mini
	C13048	C13485	C18512



Protection analogs and screw access guidance

Apical Drill⁴

	38853
--	-------



Drill Guide Multi-unit⁴

	38917
--	-------



Protection Analog Multi-unit⁴

	5/pkg
	29123



Abutment replicas

Abutment Replica Multi-unit⁴

	1/pkg	5/pkg
	38918	38919



Multi-unit Replica for IOS Model – Elos Accurate Analog

	PMA-MUA45-1
--	-------------



Elos Accurate Analog Insertion tools

Elos Accurate Analog Insertion Pin	PMA-AIP-2
Elos Accurate Analog Insertion Screw	PMA-AIS-5
Elos Accurate Analog Pliers	PMA-AIT-1



⁴ Compatible with all Multi-unit Abutments.

⁵ 15 mm guide pin included.

Clinical and laboratory screws

Screw Multi-unit Angled Abutment

NP	36892
----	-------

RP/WP	37893
-------	-------



Prosthetic Screw Multi-unit Abutment

NP/RP/WP	29285
----------	-------



Prosthetic Screw Multi-unit Abutment – Omnigrip™ Mini

NP/RP/WP	301203
----------	--------



Laboratory Screw Multi-unit Angled Abutment

NP	37896
----	-------

RP/WP	37897
-------	-------



Laboratory Prosthetic Screw Multi-unit Abutment

NP/RP/WP	29287 (5/pkg)
----------	---------------



Laboratory Prosthetic Screw Multi-unit Abutment – Omnigrip Mini

NP/RP/WP	301201
----------	--------



Screwdrivers

Screwdriver Machine Multi-unit

21 mm
29158



Screwdriver Machine Unigrip

20 mm	25 mm	30 mm	35 mm
29151	29152	29153	29154



Screwdriver Machine Omnigrip Mini⁶

20 mm	28 mm	36 mm
300852	300853	300854



Screwdriver Manual Multi-unit

25 mm
29156



Screwdriver Manual Unigrip

20 mm	28 mm	36 mm
29148	29149	29150



Screwdriver Manual Omnigrip Mini⁶

20 mm	36 mm
300855	300856



⁶ Compatible only with Zirconia bridge.

Prosthetic PureSet Basic

87301

PUR0500	Prosthetic PureSet Tray
301232	Wall Chart Prosthetic PureSet Basic

Unigrip Screwdriver Machine

20 mm	29151
30 mm	29153



Multi-unit Screwdriver Machine

21 mm	29158
-------	-------



Handle for Machine Instruments

	29161
---	-------


Omnigrip Screwdriver Machine

20 mm	37379
30 mm	37381



Manual Torque Wrench Prosthetic

	29165
---	-------



**Provide optimized
esthetics**

NobelProcera® fixed restorations

Choose from a variety of precision-manufactured CAD/CAM suprastructures in titanium and zirconia and combine these with acrylic, composite or porcelain. Other options include individualized crowns on an implant bridge for maximum esthetics.

Customized implant bridges



Titanium implant bridge framework



Zirconia implant bridge framework



Zirconia implant bridge

Customized implant bars



Hybrid bar



Montreal bar with metallic lingual

ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

ORDER BY PHONE

Call our customer service team, or contact your sales representative.

LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

nobelbiocare.com/warranty

nobelbiocare.com/all-on-4



88398 A GB 2310 Printed in the EU © Nobel Biocare Services AG, 2023. All rights reserved. Distributed by: Nobel Biocare. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. Universal Base Multi-Unit Abutments, Elos Accurate Position Locators, Elos Accurate Analog, Elos Accurate IO drivers and tools are manufactured by Elos Medtech Pinol A/S, Engvej 33, 3330 Gørlose, Denmark and distributed by Nobel Biocare Services AG.