


The logo for X-Guide, featuring a stylized 'X' made of three parallel lines in green and blue, followed by the word 'GUIDE' in a bold, blue, sans-serif font.

**X-GUIDE**

The logo for Nobel Biocare, consisting of a large, stylized 'N' in a dark blue color, followed by the words 'Nobel Biocare' in a smaller, dark blue, sans-serif font.

**Nobel  
Biocare™**

A close-up photograph of a surgical tracker system mounted on a clear, anatomical model of a human head. The tracker consists of a metal arm with a checkered pattern on its side, attached to the model via several screws. The background is a solid teal color.

Perform  
**edentulous**  
CASES WITH  
**X-Guide®**

**EDX system benefits**

**Tracker arm designed to clear the lip and provide retraction**

**Secure and stable attachment**

**Patient tracker remains out of the way during surgery**

**Navigate a wide range of cases including terminal dentition, full-arch, and dual full-arch cases**


**Special bone-screw design – flapped or flapless**

**Different arm designs based on handedness of the user**

# Elevate your X-Guide game

## Full-arch workflow with X-Mark™

1. Scan and plan
2. Attach patient tracker using edentulous hardware
3. Virtually register the patient with X-Mark
4. Perform surgery



[nobelbiocare.com/x-guide](https://nobelbiocare.com/x-guide)



GMT 73545 GB 2102 © Nobel Biocare Services AG, 2021. All rights reserved. Distributed by: Nobel Biocare. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Please contact the local Nobel Biocare sales office for current product assortment and availability. Disclaimer: For prescription use only. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. X-Guide is manufactured by X-Nav Technologies, LLC and distributed by Nobel Biocare.