



Design with



3shape▶

exocad

NobelProcera®
Scan and
Design Services




NobelProcera®

Product overview

nobelbiocare.com/lab

Implant crowns and abutments




	Implant abutments					Implant crown	
Description/features	Titanium abutment	Titanium ASC Abutment – Angulated screw channel 0-25°	Zirconia abutment	Zirconia abutment – Mechanically retained adapter	Zirconia abutment – Mechanically retained adapter – Angulated screw channel 0-25°	Full-contour/cut-back zirconia implant crown Full-contour/cut-back zirconia implant crown – Angulated screw channel 0-25°	Full-contour/cut-back zirconia implant crown – Mechanically retained adapter – Angulated screw channel 0-25°
Tooth replacement options	Single tooth	Single tooth	Single tooth	Single tooth	Single tooth	Single tooth	Single tooth
Available shades	Titanium 	Anodized Titanium or  Titanium 	Zirconia  White Light Medium Intense			Full-contour zirconia  OM2 A1 A2 A3 A3.5 B1 B2 C1 C2 D2	Full-contour zirconia  S0 S1 S2 S3 S4 S5 S6 S7
Platform							
Internal conical connection	3.0	●					
	NP, RP, WP	●	●		●		●
Internal tri-channel connection	NP, RP, WP, 6.0	●	●		●		
External hex connection	NP, RP, WP	●		●			
NobelBiocare N1™ Base	NP, RP					●*	

* N1 base sold separately

All Nobel Biocare restorations are co-packed with clinical screws.

Implant bridges



Description/Features		Zirconia implant bridge Full-contour, cut-back and/or framework – ASC 0-25 degrees for CC implant Level and Multi-unit Abutments – Mechanically retained adapter for CC implant-level full-contour/cut-back zirconia implant bridge	Titanium implant bridge Framework
Tooth replacement options		2-14 units	2-14 units
Available shades		Full-contour zirconia or esthetic zirconia 	Titanium 
Platform			
Internal conical connection	3.0		
	NP, RP, WP	●	●
Internal tri-channel connection	NP, RP, WP, 6.0	●	●
External hex connection	NP, RP, WP	●	●
Abutment connection Multi-unit Abutment	NP, RP, WP	●	●

Also available for other implant systems

On Multi-unit Abutment level: Biomet 3i®, Astra Tech Implant System™, BioHorizons®, Neodent®, Straumann®, Zimmer®

Implant bars

Available only from DTX Studio™ Lab software or services via DTX Studio™ Go



Description/Features		Titanium fixed implant bars	Titanium fixed-removable implant bars
Tooth replacement options		2-14 units	2-14 units
Available shades		Titanium 	Titanium 
Platform			
Internal conical connection	3.0		
	NP, RP, WP	●	●
Internal tri-channel connection	NP, RP, WP, 6.0	●	●
External hex connection	NP, RP, WP	●	●
Abutment connection Multi-unit Abutment	NP, RP, WP	●	●

Example of other implant systems

On implant level (use Scan and Design services): Straumann® Bone Level, Straumann® Standard/Standard Plus, Astra Tech Implant System™, Biomet 3i® Certain Internal, Biomet 3i® External Hex, Zimmer® Hex, Camlog® K-Series, Ankylos®

On Multi-unit Abutment level: Biomet 3i®, Astra Tech Implant System™, BioHorizons®, Neodent®, Straumann®, Zimmer®

Cemented crowns & bridges



Tooth replacement options	1 unit
Available shades	Full-contour zirconia
	OM2 A1 A2 A3 A3.5
	B1 B2 C1 C2 D2



Tooth replacement options	2–5 units
Available shades	Full-contour zirconia
	OM2 A1 A2 A3 A3.5
	B1 B2 C1 C2 D2

Single-unit screws

Internal conical connection



Titanium abutment



Titanium ASC Abutment



Zirconia abutment




Full-contour/cut-back
zirconia implant crown

Platform		Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab
Internal conical connection		 37890	 36805						
NP	3.0								
	NP	 37891	 37894	 301721	 301722	 37367	 37374	 37367	 37374
	RP	 37892	 37895 (5/pkg)	 301723	 301724	 37606	 37607	 37606	 37607
WP	WP	 37892	 37895 (5/pkg)	 301723	 301724	 37606	 37607	 37606	 37607



Screwdriver Manual Unigrip™



Screwdriver Manual Omnigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Single-unit screws

Internal tri-channel connection



Titanium abutment



Titanium ASC Abutment




Customized zirconia abutment

Platform		Clinical	Lab	Clinical	Lab	Clinical	Lab
Internal tri-channel connection	NP	 36818	 31170	 301725	 301726	 28837	
	RP	 29475	 29293 (5/pkg)	 301727	 301728	 28816	
	WP	 29475	 29293 (5/pkg)	 301727	 301728	 28816	
	6.0	 29475	 29293 (5/pkg)	 301727	 301728	 28816	

All Nobel Biocare restorations are co-packed with clinical screws.



Screwdriver Manual Unigrip™



Screwdriver Manual Omnigrip™

Single-unit screws

External hex connection




Titanium abutment



Zirconia abutment

Platform		Clinical	Lab	Clinical	Lab
External hex connection	NP	 29282	 37894	 31171	 37894
	RP	 29283	 38417	 28815	 37895 (5/pkg)
	WP	 29284	 31169	 28844	



Screwdriver Manual Unigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Single-unit screws


Nobel Biocare N1™ Base Connection



Implant abutment

Implant crown

Platform		Clinical	Lab	Clinical	Lab
Nobel Biocare N1™ Base Connection	NP	 301031	 301035	 301031	 301035
	RP	 301032	 301036	 301032	 301036



Screwdriver Manual Omnigrip™ Mini

All Nobel Biocare restorations are co-packed with clinical screws.

Implant bridge screws

Internal conical connection



Zirconia implant bridge
Full-contour, cut-back and/or framework



Titanium fixed-removable implant bars




Titanium fixed implant bars



Titanium implant bridge
Framework

Platform		Clinical	Lab	Clinical	Lab
Internal conical connection	3.0				
	NP	 37367	 37374	 37891	 37894
	RP	 37606	 37607	 37892	 37895 (5/pkg)
	WP	 37606	 37607	 37892	 37895 (5/pkg)

All Nobel Biocare restorations are co-packed with clinical screws.



Screwdriver Manual Omnigrip™



Screwdriver Manual Unigrip™

Implant bridge screws

Internal tri-channel connection



Zirconia implant bridge
Full-contour, cut-back and/or framework



**Titanium fixed-removable
implant bars**



Titanium fixed implant bars



Titanium implant bridge
Framework

Platform		Clinical	Lab	Clinical	Lab
Internal tri-channel connection	NP	 28837		 36818	 31170
	RP	 28816		 29475	 38419 29293 (5/pkg)
	WP	 28816		 29475	 38419 29293 (5/pkg)
	6.0	 28816		 29475	 38419 29293 (5/pkg)

All Nobel Biocare restorations are co-packed with clinical screws.



Screwdriver Manual Unigrip™

Implant bridge screws

External hex and Multi-unit
Abutment connection



Zirconia implant bridge
Full-contour, cut-back and/or framework



**Titanium fixed-removable
implant bars**




Titanium fixed implant bars



Titanium implant bridge Framework

Platform		Clinical	Lab	Clinical	Lab
External hex connection	NP	 31171	 37894	 29282	 37894
	RP	 28815	 37895 (5/pkg)	 29283	 38417
	WP	 28844		 29284	 31169
Multi-unit Abutment connection					
	NP, RP, WP	 301203	 301201	 29285	 29287 (5/pkg)
Multi-unit Abutment External Hex WP	WP	 301200	 301202	 29286	 31163

All Nobel Biocare restorations are
co-packed with clinical screws.



Screwdriver Manual Unigrip™



Screwdriver Manual Omnigrip™ Mini

NobelProcera® shade selection

Full-contour zirconia and esthetic zirconia

Shade	VITA classic shade
 OM2	Bleach Shades
 A1	A1
 A2	A2
 A3	A3
 A3.5	A3.5 / A4
 B1	B1
 B2	B2 / B3 / B4
 C1	C1
 C2	C2 / C3 / C4
 D2	D2 / D3 / D4

	Multi layer	Yttria content
Full-contour zirconia	No	3y
Esthetic zirconia	Yes	4–5y

Original zirconia

Nobel Biocare shade	VITA classic shade
 White	Bleach shades
 Light	A1 / B1
 Medium	A2 / A3 / B2 / C1 / C2 / D2
 Intense	A3.5 / A4 / C3 / C4 / D3 / D4

Full-contour zirconia (for single units)

Nobel Biocare shade	VITA classic shade
 S0	OM1 / OM2 / OM3
 S1	Bleach with yellow tinge
 S2	A1 / B1 / C1 / D2 (incisal)
 S3	A2 / B2
 S4	A3 / A3.5 / B3 / C2 / D3 / D4 / D2 (body)
 S5	C3
 S6	C4 / A4 / 5M2
 S7	B4

Scan accessories



For single-unit implant restorations	For multiple-unit implant restorations	For models produced with physical impression	For models produced with intra-oral scan (IOS) data* Elos	For models produced with intra-oral scan (IOS) data* Nobel Biocare	
Elos IO/ Desktop Scan Body	Elos IO/ Desktop Scan Body	Plaster model implant analog	IOS model implant analog	IOS analog insertion tool	IOS model implant analog IOS analog insertion tool

Nobel Biocare implant connections

Platform								
Internal conical connection	3.0	IO 2B-D SA		36791	PMA-NBA30-1	PMA-AIT-1 PMA-AIS-5 PMA-AIP-2	38188	38200
	NP	IO 2B-A SA	IO 2B-A	36697	PMA-NBA35-1		38189	38201
	RP	IO 2B-B SA	IO 2B-B	36698	PMA-NBA43-1		38190	38203
	WP	IO 2B-C SA	IO 2B-C	37879	PMA-NBA60-1		38191	38204
Internal tri-channel connection	NP	IO 2A-A	IO 2A-A	29498	PMA-NBR35-1			
	RP	IO 2A-B	IO 2A-B	29500	PMA-NBR43-1			
	WP	IO 2A-C	IO 2A-C	29502	PMA-NBR50-1			
	6.0	IO 2A-D	IO 2A-D	29995	PMA-NBR60-1			
External hex connection	NP	IO 6A-A	IO 6A-A	31158	PMA-BRA35-1			
	RP	IO 6A-B	IO 6A-B	31159	PMA-BRA41-1			
	WP	IO 6A-C	IO 6A-C	31160	PMA-BRA51-1			
Abutment connection	NP, RP		IO 2C-A	38918	PMA-MUA45-1			
Multi-unit Abutment	WP		IO 2C-B	31162	PMA-MUA60-1			

* The IOS workflow supports single and implant bridge restorations up to 6 units.

Scan accessories



For single-unit implant restorations	For single-unit implant restorations	For models produced with physical impression	For models produced with intra-oral scan (IOS) data*	
Desktop Scan Body	IO Scan Body	Plaster model implant analog	IOS model implant analog	IOS analog insertion tool

Nobel Biocare implant connections

Platform						
N1 Base™ Connection	NP	301039	300991	301021	301024	300966
	RP	301040	300992	301022	301025	

For intraoral and desktop scanning

Abutment Connection - Multi-unit Abutment

Connection	Multi-Unit	
Platform	NP / RP / WP ¹	Brånemark WP ²
Nobel Biocare Multi-unit PoLo 8 mm	301945	301946
		
Nobel Biocare Multi-unit PoLo 11 mm	301985	301986
		
Nobel Biocare Multi-unit PoLo Link 10 mm	301947	
		
Nobel Biocare Multi-unit PoLo Link 15 mm	301948	
		
Nobel Biocare Multi-unit PoLo Link 20 mm	301948	
Screwdriver for Multi-unit PoLo	Omnigrip Mini	
Replacement Screws (5/pkg)	301983	301984

For desktop scanning

Abutment Connection - Multi-unit Abutment


Connection	Multi-Unit	
Platform	NP / RP / WP ¹	Brånemark WP ²
Multi-unit (NP/RP/WP)	301950	301951
Screwdriver for ReverseScan Multi-unit	Depending on where to be used	



¹ Complatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP external hex (Brånemark)

² Complatible ONLY with Multi-unit Abutment WP external hex (Brånemark)

Prosthetic PureSet Basic





Prosthetic PureSet™ Basic

87301*

Prosthetic Pure Set Tray	PUR0500
Wall Chart Prosthetic Pure Set Basic	301232
Unigrip Screwdriver Machine 20 mm	29151
Unigrip Screwdriver Machine 30 mm	29153
Omnigrip Screwdriver Machine 20 mm	37379
Omnigrip Screwdriver Machine 20 mm	37381
Omnigrip Mini Screwdriver Machine 28 mm	300853
Multi-unit Screwdriver Machine	29158
Handle for Machine Instrument	29161
Manual Torque Wrench Prosthetic	29165

* The content of the set may vary by market. Please contact the local Nobel Biocare sales office for current product assortment and availability. Articles can also be purchased individually.

Screwdrivers

	Machine Omnigrip™ 30 mm	Machine Unigrip™ 30 mm	Machine Omnigrip™ Mini 28 mm
Clinic	 37381	 29153	 300853
	Manual Omnigrip™ 28 mm	Manual Unigrip™ 28 mm	Manual Omnigrip™ Mini 36mm
Lab	 37377	 29149	 300856

Additional screwdrivers

Omnigrip™ Screwdriver Machine 20 mm	37379
Omnigrip™ Screwdriver Machine 25 mm	37380
Omnigrip™ Screwdriver Machine 30 mm	37381
Omnigrip™ Screwdriver Machine 35 mm	37382
Omnigrip™ Screwdriver Manual 20 mm	37376
Omnigrip™ Screwdriver Manual 28 mm	37377
Omnigrip™ Screwdriver Manual 36 mm	37378

Screwdriver Machine Unigrip™ 20 mm	29151
Screwdriver Machine Unigrip™ 25 mm	29152
Screwdriver Machine Unigrip™ 30 mm	29153
Screwdriver Machine Unigrip™ 35 mm	29154
Screwdriver Manual Unigrip™ 20 mm	29148
Screwdriver Manual Unigrip™ 28 mm	29149
Screwdriver Manual Unigrip™ 36 mm	29150

Screwdriver Machine Omnigrip™ Mini 20 mm	300852
Screwdriver Machine Omnigrip™ Mini 28 mm	300853
Screwdriver Machine Omnigrip™ Mini 36 mm	300854
Screwdriver Manual Omnigrip™ Mini 20 mm	300855
Screwdriver Manual Omnigrip™ Mini 36 mm	300856



3Shape NA System Screwdriver (NobelActive)	80610223
3Shape NR System Screwdriver (Replace Select)	80610230
3Shape NB System Screwdriver (Brånemark)	80610217

Elos Accurate™ IO Driver Mini	C18512
Elos Accurate™ IO Driver Short	C13485
Elos Accurate™ IO Driver Long	C13048

Screwdriver Manual Multi-unit 25mm	29156
------------------------------------	-------

START NOW!

Register at nobelbiocare.com/openaccess to download DME files and gain access to a reliable partner with over 30 years of experience.

nobelbiocare.com/openaccess

nobelbiocare.com/lab



GMT 94338 GB 2501 Printed in the EU © Nobel Biocare Services AG, 2025. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. IO scan bodies, PMA analogs and Elos Accurate™ tools are manufactured by Elos Medtech Pinol A/S, Engvej 33, 3330 Gørøse, Denmark and distributed by Nobel Biocare Services AG.