

# Paths to NobelProcera® restorations



## Desktop scanners

| Vendor           | Model   | Supported NobelProcera restorations                     | Software or program used  | Scan locator used                         |
|------------------|---|---|---|---|
| Nobel Biocare    | <b>LS 3</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| Nobel Biocare    | <b>2G</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*                            | Nobel Biocare                             |
| 3Shape           | <b>D640</b>   | Single-unit cement & single screw-retained restorations | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | 3Shape,<br>Nobel Biocare,<br>Elos Medtech |
|                  | <b>D500, D600, D710/720/750, D800/810/850, D900/L, D1000/2000, E1, E2, E3</b> | All restorations  | DTX Studio Lab  | Nobel Biocare,<br>Elos Medtech            |
|                  |   | All restorations, excluding implant bars                | NobelProcera Services*<br>3Shape Dental System™                     | 3Shape,<br>Nobel Biocare,<br>Elos Medtech |
| Amann Girschbach | <b>Map 200</b>  | Single-unit cement & single screw-retained restorations | DTX Studio Lab**<br>NobelProcera Services*<br>3Shape Dental System™ | Nobel Biocare,<br>Elos Medtech            |
|                  | <b>Map 400/600</b>  | All restorations  | DTX Studio Lab**<br>NobelProcera Services*<br>3Shape Dental System™ | Nobel Biocare,<br>Elos Medtech            |
| Dentalwings      | <b>7 series</b>   | Single-unit cement & single screw-retained restorations | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| DOF              | <b>Freedom HD</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| Medit            | <b>Identica T300/T500, t500</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| Sirona           | <b>inEos XS</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| SmartOptics      | <b>Activity 888</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |
| Zirkonzahn       | <b>S600</b>   | All restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™   | Nobel Biocare,<br>Elos Medtech            |

\* Today, NobelProcera Services supports implant bars on over 170 implant platforms and different type of abutments and crowns. Visit [nobelbiocare.com/services](http://nobelbiocare.com/services) to find out the most recent selection of restorations supported through NobelProcera Services.

\*\* Supported DTX Studio Lab 1.10 release

## Intra-oral scanners



| Vendor     | Model                                 | Supported NobelProcera restorations  | Software or program used  | Scan locator used       |
|------------|---------------------------------------|--|---|-------------------------|
| 3M         | <b>True Definition</b>                | Single-unit cement & single screw-retained restorations, 2–6 unit cement bridges                                     | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
| 3Shape     | <b>Trios Mono/Color<br/>Trios 2/3</b> | Single-unit cement, 2–6 unit cement bridges.<br>Single screw-retained restorations, 2–3 unit screw-retained bridges  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | 3Shape,<br>Elos Medtech |
| Carestream | <b>3500</b>                           | Single-unit cement & single screw-retained restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
|            | <b>3600</b>                           | Single-unit cement, 2–6 unit cement bridges.<br>Single screw-retained restorations, 2–3 unit screw-retained bridges  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
| iTero      | <b>I/II Generation</b>                | Single-unit cement & single screw-retained restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
|            | <b>Elements</b>                       | Single-unit cement, 2–6 unit cement bridges.<br>Single screw-retained restorations, 2–3 unit screw-retained bridges  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
| Medit      | <b>i500</b>                           | Single-unit cement, 2–14 unit cement bridges.<br>Single screw-retained restorations, 2–3 unit screw-retained bridges | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
| Sirona     | <b>CEREC Omnicam</b>                  | Single-unit cement & single screw-retained restorations, 2–6 unit cement bridges                                     | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |
|            | <b>Bluecam,<br/>Apollo DI</b>         | Single-unit cement & single screw-retained restorations  | DTX Studio Lab<br>NobelProcera Services*<br>3Shape Dental System™ | Elos Medtech            |

\* Today, NobelProcera Services supports implant bars on over 170 implant platforms and different type of abutments and crowns. Visit [nobelbiocare.com/services](http://nobelbiocare.com/services) to find out the most recent selection of restorations supported through NobelProcera Services.