



NobelProcera® OpenAccess

Product overview

Design with 3Shape®.
Produce with NobelProcera®.

nobelbiocare.com/openaccess



WE FOLLOW NO
ONE.

Implant crowns and abutments





		Implant abutments			Implant crown
Description/features		Customized titanium abutment	Customized zirconia abutment	Customized zirconia abutment – Mechanically retained adapter	Customized zirconia abutment – Mechanically retained adapter – Angulated screw channel 0-25°
Tooth replacement options		Single tooth	Single tooth	Single tooth	Single tooth
Available shades		Titanium 	Zirconia 	Zirconia 	Zirconia
Platform					
Internal conical connection	3.0	●			
	NP, RP	●	●		●
	WP	●			●
Internal tri-channel connection	NP, RP, WP, 6.0	●		●	
External hex connection	NP, RP, WP	●	●		

All Nobel Biocare restorations are co-packed with clinical screws.

Implant bridges



Description/Features	Zirconia implant bridge Full-contour, cut-back and/or framework - ASC 0-25 degrees for CC implant Level and Multi-unit Abutments - Mechanically retained adapter for CC implant-level full-contour/cut-back HTML zirconia implant bridge	Titanium implant bridge Framework
Tooth replacement options	2-14 units	2-14 units
Available shades	Full-contour zirconia 	Titanium 
Platform		
Internal conical connection	3.0 NP, RP WP	● ● ●
Internal tri-channel connection	NP, RP, WP, 6.0	●
External hex connection	NP, RP, WP	●
Abutment connection Multi-unit Abutment		●


Cemented crowns & bridges



Tooth replacement options 1 unit

Available shades


Full-contour zirconia



Tooth replacement options 2-5 units

Available shades

Full-contour zirconia



Single-unit screws

Internal conical connection



Customized titanium abutment



Customized zirconia abutment



























Customized zirconia abutment



Customized zirconia abutment



Full-contour/cut-back zirconia implant crown

Platform		Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab
Internal conical connection	3.0	 37890	 36805								
	NP	 37891	 37894	 37891	 37894			 37367	 37374	 37367	 37374
	RP	 37892	 37895 (5/pkg)	 37892	 37895 (5/pkg)			 37606	 37607	 37606	 37607
	WP	 37892	 37895 (5/pkg)					 37606	 37607	 37606	 37607



Screwdriver Manual Unigrip™



Screwdriver Manual Omnigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Single-unit screws

Internal tri-channel connection



Customized titanium abutment



Customized zirconia abutment



Customized zirconia abutment




Customized zirconia abutment



Full-contour/cut-back zirconia implant crown

Platform		Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab
Internal tri-channel connection											
NP											
		36818	31170			28837					
	RP										
			29475	29293 (5/pkg)			28816				
WP											
		29475	29293 (5/pkg)			28816					
6.0											
		29475	29293 (5/pkg)			28816					



Screwdriver Manual Unigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Single-unit screws

External hex connection



Customized titanium abutment



Customized zirconia abutment














Customized zirconia abutment



Customized zirconia abutment



Full-contour/cut-back zirconia implant crown

Platform		Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab	Clinical	Lab
External hex connection	NP	 29282	 37894	 31171	 37894						
	RP	 29283	 38417	 28815	 37895 (5/pkg)						
	WP	 29284	 31169	 28844							



Screwdriver Manual Unigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Implant bridge screws













Internal conical connection




Zirconia implant bridge
Full-contour, cut-back and/or framework



Titanium implant bridge
Framework

Platform		Clinical	Lab	Clinical	Lab
Internal conical connection	3.0				
	NP	 37367	 37374	 37891	 37894
	RP	 37606	 37607	 37892	 37895 (5/pkg)
WP	 37606	 37607	 37892	 37895 (5/pkg)	



Screwdriver Manual Omnigrip™



Screwdriver Manual Unigrip™

All Nobel Biocare restorations are co-packed with clinical screws.

Implant bridge screws













Internal tri-channel connection



Zirconia implant bridge
Full-contour, cut-back and/or framework



Titanium implant bridge
Framework

Platform		Clinical	Lab	Clinical	Lab
Internal tri-channel connection	NP	 28837		 36818	 31170
	RP	 28816		 29475	 38419 29293 (5/pkg)
	WP	 28816		 29475	 38419 29293 (5/pkg)
	6.0	 28816		 29475	 38419 29293 (5/pkg)

All Nobel Biocare restorations are co-packed with clinical screws.



Screwdriver Manual Unigrip™

Implant bridge screws



















External hex and Multi-unit
Abutment connection




Zirconia implant bridge
Full-contour, cut-back and/or
framework




Titanium implant bridge
Framework

Platform		Clinical	Lab	Clinical	Lab
External hex connection	NP	 31171	 37894	 29282	 37894
	RP	 28815	 37895 (5/pkg)	 29283	 38417
	WP	 28844		 29284	 31169
Multi-unit Abutment connection	NP, RP, WP	 301203	 301201		 29287 (5/pkg)
	Multi-unit Abutment External Hex WP	 301200	 301202	 29286	 31163



Screwdriver Manual Unigrip™

All Nobel Biocare restorations are
co-packed with clinical screws.



Screwdriver Manual Omnigrip™ Mini

Scan accessories



For single-unit implant restorations	For multiple-unit implant restorations	For models produced with physical impression	For models produced with intra-oral scan (IOS) data**	
IO Scan Body/ Desktop Scan Body	IO Scan Body/ Desktop Scan Body	Plaster model implant analog	IOS model implant analog	IOS analog insertion tool







Nobel Biocare implant connections

Platform						
Internal conical connection	3.0	IO 2B-D SA	IO 2B-D	36791	38188	38200
	NP	IO 2B-A SA	IO 2B-A	36697	38189	38201
	RP	IO 2B-B SA	IO 2B-B	36698	38190	38203
	WP	IO 2B-C SA	IO 2B-C	37879	38191	38204
Internal tri-channel connection	NP	IO 2A-A	IO 2A-A	29498	MA-NBR35-1	AT-18-1
	RP	IO 2A-B	IO 2A-B	29500	MA-NBR43-1	AT-20-1
	WP	IO 2A-C	IO 2A-C	29502	MA-NBR50-1	AT-20-1
	6.0	IO 2A-D	IO 2A-D	29995	MA-NBR60-1	AT-20-1
External hex connection	NP	IO 6A-A	IO 6A-A	31158	MA-BRA35-1*	AT-16-1
	RP	IO 6A-B	IO 6A-B	31159	MA-BRA41-1*	AT-20-1
	WP	IO 6A-C	IO 6A-C	31160	MA-BRA51-1*	AT-25-1*
Abutment connection	NP, RP		IO 2C-A	38918		
Multi-unit Abutment	WP		IO 2C-B	38918		

* Needs to be ordered from Elos on elos.net.dynamicweb-cms.com. All other articles can be purchased on store.nobelbiocare.com

** The IOS workflow supports single and implant bridge restorations up to 3 units.

Screwdrivers

	Machine Omnigrip™ 30 mm	Machine Unigrip™ 30 mm	Machine Omnigrip™ Mini 28 mm
Clinic	 37381	 29153	 300853
	Manual Omnigrip™ 28 mm	Manual Unigrip™ 28 mm	Manual Omnigrip™ Mini 36mm
Lab	 37377	 29149	 300856

Additional screwdrivers

Omnigrip™ Screwdriver Machine 20 mm	37379
Omnigrip™ Screwdriver Machine 25 mm	37380
Omnigrip™ Screwdriver Machine 30 mm	37381
Omnigrip™ Screwdriver Machine 35 mm	37382
Omnigrip™ Screwdriver Manual 20 mm	37376
Omnigrip™ Screwdriver Manual 28 mm	37377
Omnigrip™ Screwdriver Manual 36 mm	37378

Screwdriver Machine Unigrip™ 20 mm	29151
Screwdriver Machine Unigrip™ 25 mm	29152
Screwdriver Machine Unigrip™ 30 mm	29153
Screwdriver Machine Unigrip™ 35 mm	29154
Screwdriver Manual Unigrip™ 20 mm	29148
Screwdriver Manual Unigrip™ 28 mm	29149
Screwdriver Manual Unigrip™ 36 mm	29150

Screwdriver Machine Omnigrip™ Mini 20 mm	300852
Screwdriver Machine Omnigrip™ Mini 28 mm	300853
Screwdriver Machine Omnigrip™ Mini 36 mm	300854
Screwdriver Manual Omnigrip™ Mini 20 mm	300855
Screwdriver Manual Omnigrip™ Mini 36 mm	300856



Elos Accurate Driver Long 0.9 Hex	C13048
Elos Accurate IO Driver Short 0.9 Hex	C13485
Screwdriver Manual Multi-unit 25mm	29156

START NOW!

Register at nobelbiocare.com/openaccess to download DME files and gain access to a reliable partner with over 30 years of experience.

nobelbiocare.com/openaccess

nobelbiocare.com/openaccess



GMT 73183 GB 2102 © Nobel Biocare Services AG, 2021. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.