



# Rescue instrumentation

Solutions for  
EVERY SITUATION

Product overview





## ABUTMENT RETRIEVAL

Removal of an abutment held by a tight conical connection is simple with the Abutment Release Pin.

For titanium and zirconia restoration, including the NobelProcera® Angulated Screw Channel (ASC) abutment and the Full-contour Zirconia (FCZ) implant crown, removal is simplified with a range of dedicated abutment retrieval instruments.



## ABUTMENT / CLINICAL SCREW REMOVAL

In the event that an abutment screw breaks, there's a safe solution, even when part of the screw remains in the screw channel.

While the Abutment Screw Remover will resolve the majority of such instances efficiently, Reverse Drills and Abutment Retrieval Instruments enable retrieval in more complex cases.



## IMPLANT RETRIEVAL

Managing a retrieval case where the implant driver no longer engages?

The Implant Retrieval Instrument is designed precisely for this scenario, regardless of whether it's due to implant fracture, a damaged connection or strong osseointegration. For more difficult cases, Trepine Drills are available.

# Implant retrieval instruments

## Implant Retrieval Instruments

		Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
22 mm	37471			NP, RP	NP, RP
31 mm	37472			NP, RP	NP, RP
31 mm	37473		3.0		
22 mm	37474		NP		WP
22 mm	37475		RP	WP	
22 mm	37476			6.0	
22 mm	37927		WP		
31 mm	301294	NP	3.0		
22 mm	301295	RP	RP	WP	



## Implant Rescue Collar

Ø		Tri-channel
3.5 mm	37478	NP
4.3 mm	37479	RP



## Trephine Drills

Ø		Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
3.2/4.0 mm	37477		3.0		
3.8/4.6 mm	37928	NP	NP	NP	NP
4.4/5.2 mm	37929	RP	RP	RP, NobelReplace Platform Shift	NP RP
5.2/6.2 mm	37930		RP	WP, NobelReplace Platform Shift	RP WP
5.6/6.6 mm	37931		WP		
6.2/7.0 mm	37932			6.0, NobelReplace Platform Shift	WP WP



Handle for Implant Rescue Collar & Guides

37480



# Abutment screw retrieval instruments

## Abutment Screw Remover

	Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
3.0		37505		
NP	37505, 37506, 37507	37506	37506	37506
RP	37505, 37506, 37507	37507	37507	37507
WP, 6.0		37507	37507	37507



## Nobel Biocare N1™ Base Screw Removal Tool

	Trioval Conical Connection
NP, RP	301702



## Screw Tap Repair

		Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
M1.4	37497		3.0		
M1.6	37498		NP		NP
M1.8	37499			NP	
M2	37517		RP, WP	RP, WP, 6.0	RP
M2.5	37500				WP
	301151	NP			
	301152	RP			



# Abutment screw retrieval instruments

## Abutment Screw Retrieval Reverse Drill

	Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
3.0		37501, 301296		
NP	301296	37501, 301296	37501	37501
RP	301296	37502	37502	37502
WP, 6.0		37502	37502	37502



## Abutment Screw Retrieval Instrument

	Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
3.0		37503, 301297		
NP	301297	37503, 301297	37503	37503
RP	301297	37504	37504	37504
WP, 6.0		37504	37504	37504



## Rescue Drill Guide

	Trioval Conical Connection	Conical Connection	Tri-channel	External Hex
3.0		37485		
NP	300919	37486	37488	37482
RP	300920	37487	37489	37483
WP		37933	37490	37484
6.0			37491	



Handle for Implant Rescue Collar & Guides

37480



Handle for Machine Instruments

29161



# Abutment retrieval instruments

## Abutment Retrieval Instrument Zirconia

	Conical Connection
NP	37512
RP, WP	37882



## Abutment Retrieval Instrument Titanium

	Trioval Conical Connection	Conical Connection
NP	301521 (long), 301533 (short)	36247
RP	301522 (long), 301534 (short)	37881
WP		37881



## Abutment Release Pin

	Conical Connection
3.0	37509
NP	37510
RP, WP	37511



### Consistent color coding

Prosthetic components are color coded for accurate and fast identification.

<b>3.0</b>	<b>NP</b>	<b>RP</b>	<b>WP</b>
Platform 3.0	Narrow Platform	Regular Platform	Wide Platform

# Sets

## Implant retrieval

### 87309 Implant Retrieval Set

Set includes:		QTY
301343	Implant Retrieval Set Organizer	1
37471	Impl Retrieval Instrument Ext Hex & Tri-Channel NP/RP 22 mm	1
37472	Impl Retrieval Instrument Ext Hex & Tri-Channel NP/RP 31 mm	1
37473	Implant Retrieval Instrument CC 3.0 31 mm	1
37474	Implant Retrieval Instrument CC NP & Ext Hex WP 22 mm	1
37475	Implant Retrieval Instrument CC RP & Tri-Channel WP 22 mm	1
37476	Implant Retrieval Instrument Tri-Channel 6.0 22 mm	1
37927	Implant Retrieval Instrument CC WP 22 mm	1
37477	Trephine Drill 3.2/4.0 mm	1
37928	Trephine Drill 3.8/4.6 mm	1
37929	Trephine Drill 4.4/5.2 mm	1
37930	Trephine Drill 5.2/6.2 mm	1
37931	Trephine Drill 5.6/6.6 mm	1
37932	Trephine Drill 6.2/7.0 mm	1
37478	Implant Rescue Collar Tri-Channel Ø 3.5	1
37479	Implant Rescue Collar Tri-Channel Ø 4.3	1
37480	Handle for Implant Rescue Collar & Drill Guides	1



**Note:** Please note that the content of the sets may vary by market and that the set article number refers to an internal ordering number and is not visible on the delivery note or invoice. The articles in the sets can also be purchased individually.



# Sets

## Abutment screw retrieval

### 87308 Abutment Screw Retrieval Set

Set includes:	QTY
301344 Abutment Screw Retrieval Set Organizer	1
37482 Rescue Drill Guide External Hex NP	1
37483 Rescue Drill Guide External Hex RP	1
37484 Rescue Drill Guide External Hex WP	1
37485 Rescue Drill Guide Conical Connection 3.0	1
37486 Rescue Drill Guide Conical Connection NP	1
37487 Rescue Drill Guide Conical Connection RP	1
37933 Rescue Drill Guide Conical Connection WP	1
37488 Rescue Drill Guide Tri-Channel NP	1
37489 Rescue Drill Guide Tri-Channel RP	1
37490 Rescue Drill Guide Tri-Channel WP	1
37491 Rescue Drill Guide Tri-Channel 6.0	1
37501 Abutment Screw Retrieval Reverse Drill 3.0/NP	2
37502 Abutment Screw Retrieval Reverse Drill RP/WP/6.0	2
37503 Abutment Screw Retrieval Instrument 3.0/NP	1
37504 Abutment Screw Retrieval Instrument RP/WP/6.0	1
37505 Abutment Screw Remover 3.0	1
37506 Abutment Screw Remover NP	1

### 87308 Abutment Screw Retrieval Set

Set includes:	QTY
37507 Abutment Screw Remover RP/WP/6.0	1
37497 Screw Tap Repair M1.4	1
37498 Screw Tap Repair M1.6	1
37499 Screw Tap Repair M1.8	1
37517 Screw Tap Repair M2.0	1
37500 Screw Tap Repair M2.5	1
37480 Handle for Implant Rescue Collar & Drill Guides	1
29161 Handle for Machine Instruments	1



**Note:** Please note that the content of the sets may vary by market and that the set article number refers to an internal ordering number and is not visible on the delivery note or invoice. The articles in the sets can also be purchased individually.

# Sets

## Abutment retrieval

### 87338 Abutment Retrieval Set

Set includes:		QTY
37509	Abutment Release Pin CC 3.0	1
37510	Abutment Release Pin CC NP	1
37511	Abutment Release Pin CC RP/WP	1
37512	Abutment Retrieval Instrument Zirconia CC NP	1
37882	Abutment Retrieval Instrument Conical Connection RP/WP Zirconia	1
36247	Abutment Retrieval Tool Conical Connection NP Titanium	1
37881	Abutment Retrieval Instrument Conical Connection RP/WP Titanium	1

**Note:** Please note that the content of the sets may vary by market and that the set article number refers to an internal ordering number and is not visible on the delivery note or invoice. The articles in the sets can also be purchased individually.



---

## ORDER ONLINE

Order our complete range of implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](https://nobelbiocare.com/store)

[nobelbiocare.com/rescue](https://nobelbiocare.com/rescue)



88565 GB 2303 Printed in the EU © Nobel Biocare Services AG, 2023. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.