


NobelReplace® Conical Connection Tapered Drill Usage Chart



This chart is designed to help you track the number of times each tapered drill and screw tap is used.

Tapered drills and screw taps are reusable and should be replaced after 20 uses, or when cutting efficiency declines. Worn-out and damaged drills need to be discarded and replaced with new sharp drills. The over-extended use may cause bone overheating and lead to implant failure. Tapered drills are to be cooled internally via irrigation and require specific cleaning procedures prior to sterilization. See current cleaning and sterilization guidelines for details: nobelbiocare.com/sterilization

Account number: _____

Implant Ø		Article No.	Size	Please tick for each drill used per procedure																			
3.5mm		32075	3.5 × 8 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29367	3.5 × 10 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		36113	3.5 × 11.5 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29368	3.5 × 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29369	3.5 × 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29377	Dense Bone 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29378	Dense Bone 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		36717	Screw Tap	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
4.3mm		32076	4.3 × 8 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29370	4.3 × 10 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		36114	4.3 × 11.5 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29371	4.3 × 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29372	4.3 × 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29380	Dense Bone 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29381	Dense Bone 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
			32090	Screw Tap	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
5.0mm		32077	5.0 × 8 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29373	5.0 × 10 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		36115	5.0 × 11.5 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29374	5.0 × 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29375	5.0 × 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29383	Dense Bone 13 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
		29384	Dense Bone 16 mm	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
			32091	Screw Tap	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

Note: Tapered drills and screw taps are used for NobelReplace Conical Connection TiUltra, NobelReplace Conical Connection, NobelReplace Conical Connection PMC, Replace Select Tapered, Replace Select Tapered PMC, NobelReplace Tapered Groovy, NobelReplace Platform Shift and NobelDirect. They are made from surgical stainless steel and coated with an amorphous diamond coating, which gives them their black color.

nobelbiocare.com/en-us