

Product overview



Rescue instrumentation

Solutions
FOR EVERY
SITUATION

Abutment / clinical screw removal

In the event that an abutment screw breaks, there's a safe solution, even when part of the screw remains in the screw channel.

While the Abutment Screw Remover will resolve the majority of such instances efficiently, Reverse Drills and Abutment Retrieval Instruments enable retrieval in more complex cases.



Implant retrieval

Managing a retrieval case where the implant driver no longer engages?

The Implant Retrieval Instrument is designed precisely for this scenario, regardless whether it's due to implant fracture, a damaged connection or strong osseointegration. For more difficult cases, Trephine Drills are available.

Abutment retrieval

Removal of an abutment held by a tight conical connection is simple with the Abutment Release Pin.

For titanium and zirconia restoration, including the NobelProcera ASC Abutment and the FCZ Implant Crown, removal is simplified with a range of dedicated abutment retrieval instruments.

Surgical components

Implant retrieval

Implant Retrieval Instruments

		Conical Connection	Tri-channel	External Hex
22 mm	37471		NP, RP	NP, RP
31 mm	37472		NP, RP	NP, RP
31 mm	37473	3.0		
22 mm	37474	NP		WP
22 mm	37475	RP	WP	
22 mm	37476		6.0	
22 mm	37927	WP		

STERILE R



Implant Rescue Collar

Ø		Tri-channel
3.5 mm	37478	NP
4.3 mm	37479	RP



Abutment screw retrieval

Abutment Screw Remover

	Conical Connection	Tri-channel	External Hex
3.0	37505		
NP	37506	37506	37506
RP, WP, 6.0	37507	37507	37507

STERILE R



Rescue Drill Guide

	Conical Connection	Tri-channel	External Hex
3.0	37485		
NP	37486	37488	37482
RP	37487	37489	37483
WP	37933	37490	37484
6.0		37491	



Abutment Screw Retrieval Instrument

	Conical Connection	Tri-channel	External Hex
3.0, NP	37503	37503	37503
RP, WP, 6.0	37504	37504	37504

STERILE R



Abutment retrieval

Abutment Retrieval Instrument

		Conical Connection
NP	Zr	37512
RP, WP	Zr	37882
NP	Ti	36247
RP, WP	Ti	37881



Trephine Drills

Ø		Conical Connection	Tri-channel	External Hex
3.2/4.0 mm	37477	3.0		
3.8/4.6 mm	37928	NP	NP	NP
4.4/5.2 mm	37929	RP	RP, NobelReplace Platform Shift NP	RP
5.2/6.2 mm	37930	RP	WP, NobelReplace Platform Shift RP	WP
5.6/6.6 mm	37931	WP		
6.2/7.0 mm	37932		6.0, NobelReplace Platform Shift WP	WP

STERILE R



Handle for Implant Rescue Collar & Guides

37480



Screw Tap Repair

	Conical Connection	Tri-channel	External Hex
M1.4	37497	3.0	
M1.6	37498	NP	NP
M1.8	37499	NP	
M2	37517	RP, WP	RP, WP, 6.0
M2.5	37500		WP

STERILE R



Abutment Screw Retrieval Reverse Drill

	Conical Connection	Tri-channel	External Hex
3.0, NP	37501	37501	37501
RP, WP, 6.0	37502	37502	37502

STERILE R



Handle for Machine Instruments

29161



Abutment Release Pin

	Conical Connection
3.0	37509
NP	37510
RP, WP	37511



ORDER ONLINE

Order our complete range of implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

nobelbiocare.com/rescue

