

Now available

NobelProcera® Esthetic Zirconia Implant Bridge

Nobel Biocare has over 18 years of experience producing cement-free zirconia full-arch restorations. The NobelProcera Esthetic Zirconia Implant Bridge provides enhanced esthetics while maintaining the proven strength.

Up to
12-year
follow-up

98.2%
Prosthetic
survival rate¹

NobelProcera® restorations backed by real science

January 1, 2025²



132 Publications



>5,500 Patients



>7,500 Restorations



Made for
esthetics



Cement-free



Extensive
flexibility

Manufacturing excellence

NobelProcera guarantees a 100% inspection to ensure a passive fit against your digital design. Verification of dimensional tolerances as low as $\pm 5\mu\text{m}$.

Warranty you can rely on

As an original prosthetic, the product enables life-time warranty of the implant and comes with a 10-year prosthetic warranty.



100% verification of passive fit as part of the manufacturing process


10 YEAR
PRODUCT
WARRANTY
NobelProcera®

Product information

- Multi-layered zirconia material
- Available in 10 VITA shades
- Between 4 and 5Y (4–5% mol yttria depending on layer and shade)
- Available to design as full contour, cut-back, thimble, and framework with or without soft tissue
- Restore on Multi-unit Abutment or direct to Nobel Biocare implant connections
- Angulated screw channel for Multi-unit Abutment and conical connection platforms
- Recommended for anterior and posterior indications from 2 to 14 units



Case courtesy of Dr. Francesco Mintrone



Michiel Wouters
Dental Technician
The Netherlands

"My first experience with this material is just great. The translucency combined with the shade matching is among the best materials I have ever used, which allows me to spend less time to finalize my cases."



Dr. Markus Blatz
DMD, PhD
USA

"This material sets new standards in optical and physical properties, accuracy, and stability. Combined with industrial precision milling, it's a winning solution for demanding cases."



Sean Han
Ceramist & Full Arch Specialist
USA

"Without layering, I can craft highly esthetic restorations that capture light, depth, and natural beauty – all with precision and efficiency. It's not just about saving time; it's about pushing the boundaries of dentistry."



Watch Sean finalize a case in under two hours



Get the latest updates on NobelProcera



Clinical study Prof. Pozzi et al.



Discover more about our All-on-4® digital workflow

1. Pozzi A, Arcuri L, Fabbri G, Singer G, Londono J. Long-term survival and success of zirconia screw-retained implant-supported prostheses for up to 12 years: A retrospective multicenter study. J Prosthet Dent. 2021 Jun 26;S0022-3913(21)00265-1. Online ahead of print. | 2. Data on file

MKT-5954 Rev 00 (02/25, GMT 95242) © Nobel Biocare USA, LLC, 2025. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or law in your jurisdiction restricts this device to sale by or on the order of a dentist or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings, and precautions.