



# Universal Base

Product overview



[nobelbiocare.com/universal-base](https://nobelbiocare.com/universal-base)

# Universal Base – Engaging

Single unit restoration

35  
Uni



## Universal Base engaging



Connection	CC			External hex			Tri-channel		
	NP	RP	WP	NP	RP	WP	NP	RP	WP
Platform									
Article nr. collar 1.5 mm	38213	38214	38215	38227	38228	38229	38171	38172	38173
Article nr. collar 3.0 mm	38216	38217	38218	38230	38231	38232	38174	38225	38226
Post height	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm
Diameter	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm	Ø 4.775 mm	Ø 4.775 mm	Ø 6.515 mm
Screw	37891	37892	37892	29282	29283	29284	36818	29475	29475
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37894	37895*	37895*	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A SA	IO 2B-B SA	IO 2B-C SA	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B SA Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300457	300458	300459	300460	300461	300462	300452	300453	300454
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

\* 5/pkg

Universal Base Engaging: Burn-out coping and clinical screw included.

# Universal Base – Non-Engaging

Multiple unit restoration – Conical design

35  
Uni



## Universal Base non-engaging – Conical design



Connection	CC			External hex			Tri-channel		
	NP	RP	WP	NP	RP	WP	NP	RP	WP
Platform									
Article nr.	301235	301236	301237	301240	301241	301242	301243	301244	301245
Abutment collar height	1.8 mm	1.0 mm	1.0 mm	1.5 mm	1.0 mm	0.8 mm	1.0 mm	0.6 mm	0.6 mm
Post height	4.1 mm	4.1 mm	4.1 mm	4.1mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm	4.1 mm
Diameter	Ø 5.0 mm	Ø 5.0 mm	Ø 5.7 mm	Ø 5.0 mm	Ø 5.0 mm	Ø 6.0 mm	Ø 5.0 mm	Ø 5.0 mm	Ø 5.7 mm
Screw	37891	37892	37892	29282	29283	29284	36818	29475	29475
Screwdriver	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37894	37895*	37895*	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A	IO 2B-B	IO 2B-C	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300470	300471	300472	300463	300464	300465	300466	300467	300468
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

\* 5/pkg

Universal Base non-engaging: clinical screw included.

# Universal Base – Non-Engaging

Multiple unit restoration – Straight design

35  
Omni

35  
Uni



## Universal Base non-engaging – Straight design



Connection	CC			External hex			Tri-channel		
	NP	RP	WP	NP	RP	WP	NP	RP	WP
Platform									
Article nr.	301362	301363	301364	301367	301368	301369	301370	301371	301372
Abutment collar height	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm	1.5 mm
Post height	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm	4.0 mm
Diameter	Ø 3.6 mm	Ø 4.2 mm	Ø 5.1 mm	Ø 3.9 mm	Ø 4.4 mm	Ø 5.4 mm	Ø 3.7 mm	Ø 4.5 mm	Ø 5.1 mm
Screw	37367	37606	37606	29282	29283	29284	36818	29475	29475
Screwdriver	OmniGrip	OmniGrip	OmniGrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip	Unigrip
Lab Screw	37374	37607	37607	37894	37895*	31169	31170	29293	29293
Elos Accurate™ Position Locator	IO 2B-A	IO 2B-B	IO 2B-C	IO 6A-A	IO 6A-B	IO 6A-C	IO 2A-A	IO 2A-B	IO 2A-C
Elos Accurate™ Position Locator Kit	IO 2B Kit			IO 6A Kit			IO 2A Kit		
Replica for IOS - Elos Accurate™ Analog	PMA-NBA35-1	PMA-NBA43-1	PMA-NBA60-1	PMA-BRA35-1	PMA-BRA41-1	PMA-BRA51-1	PMA-NBR35-1	PMA-NBR43-1	PMA-NBR50-1
Desktop scanbody DTX	300470	300471	300472	300463	300464	300465	300466	300467	300468
Plaster model implant analog	36697	36698	37879	31158	31159	31160	29498	29500	29502

\* 5/pkg

Universal Base non-engaging: clinical screw included.

# Universal Base – Multi unit

Multiple unit and full arch restoration – Conical design

15  
Uni



## Universal Base Multi unit – Conical design

Connection	Multi-Unit	
Platform	NP / RP <sup>1</sup>	WP <sup>2</sup>
Article nr.	301238	301239
Abutment collar height	0.9 mm	1.2 mm
Post height	4.1 mm	4.1 mm
Diameter	Ø 5.0 mm	Ø 5.7 mm
Screw	29285	29286
Screwdriver	Unigrip	Unigrip
Lab Screw	29287*	31163
Elos Accurate™ Position Locator	IO 2C-A	IO 2C-B
Elos Accurate™ Position Locator Kit	IO 2C Kit	
Replica for IOS - Elos Accurate™ Analog	PMA-MUA45-1	PMA-MUA60-1
Desktop scanbody DTX	300473	300474
Plaster model implant analog	38918	31162

\* 5/pkg

<sup>1</sup> Compatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP ex-hex (Brånemark)

<sup>2</sup> Compatible ONLY with Multi-unit Abutments WP ex-hex (Brånemark)

Universal Base non-engaging: clinical screw included.

# Universal Base – Multi unit

Multiple unit and full arch restoration – Straight design

15  
O-Mini



## Universal Base Multi unit – Straight design

Connection	Multi-Unit	
Platform	NP / RP <sup>1</sup>	WP <sup>2</sup>
Article nr.	301365	301366
Abutment collar height	0.3 mm	0.3 mm
Post height	4.5 mm	4.5 mm
Diameter	Ø 5.0 mm	Ø 5.9 mm
Screw	301203	301200
Screwdriver	OmniGrip Mini	OmniGrip Mini
Lab Screw	301201	301202
Elos Accurate™ Position Locator	IO 2C-A	IO 2C-B
Elos Accurate™ Position Locator Kit	IO 2C Kit	
Replica for IOS - Elos Accurate™ Analog	PMA-MUA45-1	PMA-MUA60-1
Desktop scanbody DTX	300473	300474
Plaster model implant analog	38918	31162

<sup>1</sup> Compatible with all Multi-unit Abutments EXCEPT Multi-unit Abutment WP ex-hex (Brånemark)

<sup>2</sup> Compatible ONLY with Multi-unit Abutments WP ex-hex (Brånemark)

Universal Base non-engaging: clinical screw included.

**On1™ Universal Base**

**Universal Abutment**

**Nobel Biocare N1™**

# On1™ Universal Base

Single unit restoration – Engaging

35  
Uni



## On1 Universal Abutment – Engaging

Connection	On1 Base		
	NP	RP	WP
Platform			
Article nr. collar 0.3 mm	38707	38709	38711
Article nr. collar 1.25 mm	38708	38710	38712
Post height	3.98 mm	3.98 mm	3.98 mm
Diameter	Ø 4.775 mm	Ø 4.775 mm	Ø 4.775 mm
Screw	38750	38750	38750
Screwdriver	Unigrip	Unigrip	Unigrip
Lab Screw	38751	38751	38751
IOS scanbody & Healing cap (4.5 mm / 5.0 mm / 6.0 mm)	38735 / - / 38734	38737 / - / 38736	38739 / 38738 / -
IOS / Desktop analog	38740	38741	38742
Insertion tool	38201	38203	38204
Desktop scanbody DTX	300437	300437	300437
Plaster model implant analog	38731	38732	38733



Delivered with the corresponding On1 Prosthetic screw and On1 Burn-out coping.



# On1™ Universal Base

Multiple unit restoration – Non-Engaging

35  
Uni



## On1 Universal Abutment – Non-engaging

Connection	On1 Base		
	NP	RP	WP
Platform			
Article nr. collar 0.3 mm	300361	300363	300365
Article nr. collar 1.25 mm	300362	300364	300366
Post height	3.98 mm	3.98 mm	3.98 mm
Diameter	Ø 4.775 mm	Ø 5.3 mm	Ø 6.515 mm
Screw	38750	38750	38750
Screwdriver	Unigrip	Unigrip	Unigrip
Lab Screw	38751	38751	38751
IOS scanbody & Healing cap (4.5 mm / 5.0 mm / 6.0 mm)	38735 / - / 38734	38737 / - / 38736	38739 / 38738 / -
IOS / Desktop analog	38740	38741	38742
Insertion tool	38201	38203	38204
Desktop scanbody DTX	300437	300437	300437
Plaster model implant analog	38731	38732	38733



Delivered with the corresponding On1 Prosthetic screw and On1 Burn-out coping.

# Universal Abutment Nobel Biocare N1™ Base

Single unit restoration – Engaging – Base level

20  
O-Mini



## Universal Abutment N1 Base – Engaging

Connection	N1 Base	
	NP	RP
Platform		
Article nr.	301006	301007
Post height	4.0 mm	4.0 mm
Screw	301029	301030
Screwdriver	OmniGrip Mini	OmniGrip Mini
Lab Screw	301033	301034
IOS scanbody	301039	301040
IOS scanbody & Healing cap	300991	300992
IOS analog	301024	301025
IOS insertion tool	300966	300966
Desktop scanbody DTX	301039	301040
Plaster model implant analog	301021	301022



Delivered with the corresponding N1 Prosthetic screw.  
ASC on N1 Bases with 20° angulation.

# Universal Abutment Nobel Biocare N1™ Base

Bridge – Non-Engaging – Base level

20  
O-Mini



Universal Abutment N1 Base – Bridge

Connection	N1 Base	
	NP	RP
Platform		
Article nr.	301009	301010
Post height	4.0 mm	4.0 mm
Screw	301029	301030
Screwdriver	OmniGrip Mini	OmniGrip Mini
Lab Screw	301033	301034
IOS scanbody	301039	301040
IOS scanbody & Healing cap	300991	300992
IOS analog	301024	301025
IOS insertion tool	300966	300966
Desktop scanbody DTX	301039	301040
Plaster model implant analog	301021	301022



Delivered with the corresponding N1 Prosthetic screw.  
ASC on N1 Bases. Max. 6 units and up to 20° divergence between implants.

# Universal Abutment Nobel Biocare N1™

Single unit restoration – Engaging – Implant level

20  
O-Mini



## Universal Abutment N1 Base – Engaging

Connection	TCC	
	NP	RP
Platform		
Article nr. collar 1.5 mm	300946	300948
Article nr. collar 3.0 mm	300947	300949
Post height	4.05 mm	4.05 mm
Screw	300968	300969
Screwdriver	OmniGrip Mini	OmniGrip Mini
Lab Screw	300972	300973
IOS scanbody	300976	300977
IOS analog	300964	300965
IOS insertion tool	300966	300966
Desktop scanbody DTX	300976	300977
Plaster model implant analog	300962	300963



Delivered with clinical screw.

# Screwdrivers & Tools

# Screwdrivers

## OmniGrip™


### OmniGrip™ Screwdriver Machine

20 mm	37379
25 mm	37380
30 mm	37381
35 mm	37382



### OmniGrip™ Screwdriver Manual

20 mm	37376
28 mm	37377
36 mm	37378



## OmniGrip™ Mini

### OmniGrip™ Mini Screwdriver Machine

20 mm	300852
28 mm	300853
36 mm	300854



### OmniGrip™ Mini Screwdriver Manual

20 mm	300855
36 mm	300856



## UniGrip™

### UniGrip™ Screwdriver Machine

20 mm	29151
25 mm	29152
30 mm	29153
35 mm	29154



### UniGrip™ Screwdriver Manual

20 mm	29148
28 mm	29149
36 mm	29150



# Tools

## Elos Accurate™ Analog insertion tool and IO driver

Elos Accurate™ Analog Insertion Pin	PMA-AIP-2
Elos Accurate™ Analog Insertion Screw	PMA-AIS-5
Elos Accurate™ Analog Pliers	PMA-AIT-1
Elos Accurate™ IO driver – mini	C18512
Elos Accurate™ IO driver – short	C13485
Elos Accurate™ IO driver – long	C13048

## Compatible softwares



3shape

exocad

dental wings

## Other tools – Manual Torque Wrench

Manual Torque Wrench Prosthetic	29165
Manual Torque Wrench Adapter	29167



Manual Torque Wrench Prosthetic Nobel Biocare N1™	301082
---	--------



For availability of compatible software and download of local production DME's, please go to <https://www.nobelbiocare.com/local-production>

---

## ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

[nobelbiocare.com/store](https://nobelbiocare.com/store)

---

## ORDER BY PHONE

Call our customer service team, or contact your sales representative.

---

## LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

[nobelbiocare.com/warranty](https://nobelbiocare.com/warranty)



[nobelbiocare.com/universal-base](https://nobelbiocare.com/universal-base)



87967 GB 2203 Printed in the EU © Nobel Biocare Services AG, 2022. All rights reserved. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. Universal Base Non-engaging Abutments, Universal Base Multi-Unit Abutments, Elos Accurate™ Position Locators, Elos Accurate™ Position Locator Kits, Replica for IOS Elos Accurate™ Analog and Elos Accurate™ Analog insertion tool and IO driver are manufactured by Elos Medtech Pinol A/S, Engvej 33, 3330 Gørløse, Denmark and distributed by Nobel Biocare.