

NEW

Screw-retained restorations for any tooth

with the new
NobelProcera® ASC
Titanium Abutment





NobelProcera® ASC Titanium Abutment

Avoid intraoral cementation

Easier handling with Omnigrip®

- Torque at an angle from 0–25°

Easier access in the posterior

- Angulate the screw-channel for easier access
- Avoid intraoral cementation

Improved esthetics

- Angulate the screw-channel out of the esthetic zone
- The golden hue presents a natural appearance in the transmucosal zone

Increased surface area for cementation

- New retention grooves

Save time by producing the crown and abutment in parallel

- Download your NobelProcera® production core-file from your customer account at nobelbiocare.com

Available for



Start now

Watch the how-to video and download our NobelProcera® libraries directly in 3Shape® or exocad

nobelbiocare.com/openaccess

89111 GB 2402 Printed in the EU © Nobel Biocare Services AG, 2024. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. 3Shape® and 3Shape Dental System™ are trademarks of 3Shape A/S. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.

