

Paths to NobelProcera® restorations



Desktop scanners

Vendor	Model	Supported NobelProcera restorations	Software or program used	Scan locator used
Nobel Biocare	LS 3	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
Nobel Biocare	2G	All restorations	DTX Studio Lab NobelProcera Services*	Nobel Biocare
3Shape	D640	Single-unit cement & single screw-retained restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	3Shape, Nobel Biocare, Elos Medtech
	D500, D600, D710/720/750, D800/810/850, D900/L, D1000/2000, E1, E2, E3	All restorations	DTX Studio Lab	Nobel Biocare, Elos Medtech
		All restorations, excluding implant bars	NobelProcera Services* 3Shape Dental System™	3Shape, Nobel Biocare, Elos Medtech
Amann Girschbach	Map 200	Single-unit cement & single screw-retained restorations	DTX Studio Lab** NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
	Map 400/600	All restorations	DTX Studio Lab** NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
Dentalwings	7 series	Single-unit cement & single screw-retained restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
DOF	Freedom HD	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
Medit	Identica T300/T500, t500	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
Sirona	inEos XS	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
SmartOptics	Activity 888	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech
Zirkonzahn	S600	All restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Nobel Biocare, Elos Medtech

* Today, NobelProcera Services supports implant bars on over 170 implant platforms and different type of abutments and crowns. Visit nobelbiocare.com/services to find out the most recent selection of restorations supported through NobelProcera Services.

** Supported DTX Studio Lab 1.10 release

Intra-oral scanners



Vendor	Model	Supported NobelProcera restorations	Software or program used	Scan locator used
3M	True Definition	Single-unit cement & single screw-retained restorations, 2–6 unit cement bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
3Shape	Trios Mono/Color Trios 2/3	Single-unit cement, 2–6 unit cement bridges. Single screw-retained restorations, 2–3 unit screw-retained bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	3Shape, Elos Medtech
Carestream	3500	Single-unit cement & single screw-retained restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
	3600	Single-unit cement, 2–6 unit cement bridges. Single screw-retained restorations, 2–3 unit screw-retained bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
iTero	I/II Generation	Single-unit cement & single screw-retained restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
	Elements	Single-unit cement, 2–6 unit cement bridges. Single screw-retained restorations, 2–3 unit screw-retained bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
Medit	i500	Single-unit cement, 2–14 unit cement bridges. Single screw-retained restorations, 2–3 unit screw-retained bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
Sirona	CEREC Omnicam	Single-unit cement & single screw-retained restorations, 2–6 unit cement bridges	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech
	Bluecam, Apollo DI	Single-unit cement & single screw-retained restorations	DTX Studio Lab NobelProcera Services* 3Shape Dental System™	Elos Medtech

* Today, NobelProcera Services supports implant bars on over 170 implant platforms and different type of abutments and crowns. Visit nobelbiocare.com/services to find out the most recent selection of restorations supported through NobelProcera Services.