


Product overview



On1™ concept

Preserving
SOFT-TISSUE
attachment



On1™ concept

Product overview

Simplify all restorative procedures

The On1 Base moves the prosthetic platform from bone level to tissue level. This makes attachment of restorative components and impression taking much easier as the connection is clearly visible.

Gain peace of mind


The On1 Base has a unique prosthetic connection, ensuring that only precision-engineered Nobel Biocare restorations are used.

Preserve soft tissue

The soft tissue is left undisturbed as the On1 Base remains in position during the entire restorative workflow and throughout the lifetime of the restoration.

Maintain surgical flexibility


The On1 concept can be used with any Nobel Biocare conical connection implant system, each designed for high primary stability and built-in platform shifting. These include NobelActive, NobelParallel and NobelReplace Conical Connection.



One base, many restorative options

Surgical components


On1 Base incl. handle	NP	RP	WP
1.75 mm	38689	38691	38693
2.5 mm	38690	38692	38694



Delivered with the corresponding On1 Clinical Screw

STERILE R

On1 Healing Cap	NP	RP	WP
1.5 mm	38695	38697	38699
2.5 mm	38696	38698	38700



STERILE R

On1 IOS Healing Cap incl. handle	NP	RP	WP
4 mm			38739
4.5 mm	38735	38737	
5 mm			38738
6 mm	38734	38736	




STERILE R

On1 Clinical Screw	NP	RP	WP
1.75 mm	38744	38746	38748
2.5 mm	38745	38747	38749



STERILE R

On1 Base Screwdriver	NP	RP	WP
Machine 21 mm		38752	
Manual 25 mm		38753	



NON STERILE

Lab components

On1 Base Replica	NP	RP	WP
	38731	38732	38733



On1 IOS Replica	NP	RP	WP
	38740	38741	38742



NON STERILE

On1 Prosthetic Lab Screw	NP	RP	WP
	38751	38751	38751




NON STERILE

On1 IOS Replica Insertion Tool	NP	RP	WP
	38201	38201	38201




NON STERILE



Easy handling


The On1 Base, IOS Healing Cap and Temporary Abutment are all delivered with a pre-mounted handle for easier placement.

Restorative components



On1 Temporary Abutment incl. handle		NP	RP	WP
	Engaging	38701	38703	38705

Delivered with the corresponding On1 Prosthetic Screw

STERILE


On1 Impression Coping		NP	RP	WP
	Closed Tray	38725	38726	38727
	Open Tray	38728	38729	38730

NON
STERILE


On1 Final Abutment		NP	RP	WP
	Universal Abutment Engaging			
	0.3 mm	301121	301123	301125
	1.25 mm	301122	301124	301126
	Esthetic Abutment Titanium			
	0.3 mm	38713	38714	38715

Delivered with the corresponding On1 Prosthetic Screw


NON
STERILE

On1 Prosthetic Screw		NP	RP	WP
		38750	38750	38750



STERILE

Manual Torque Wrench Prosthetic	
	29165

NON
STERILE

Manual Torque Wrench Adapter	
	29167

NON
STERILE

Unigrip Screwdriver					
	Manual	20 mm	28 mm	36 mm	
		29148	29149	29150	
	Machine	20 mm	25 mm	30 mm	35 mm
		29151	29152	29153	29154

NON
STERILE



Peace of mind

The warranty covers all Nobel Biocare implants including prefabricated prosthetic components. For further information visit nobelbiocare.com/warranty

ORDER ONLINE

Order our complete range of implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

ORDER BY PHONE

Call our customer service team or contact your sales representative.

nobelbiocare.com/on1

