Procera™

Solution overview for labs



PROCERA™ YOUR TRUSTED PARTNER FOR GROWTH AND EXCELLENCE

40+
years of CAD/CAM production

60+
countries served worldwide

5M+
restorations produced over decades

Case courtesy of Dr. Édouard Cristofari and Luis Saraiva



1983

Procera is born: Dr. Matts Andersson creates a breakthrough method of high-precision repeatable manufacturing of customized dental crowns.

A LEGACY THAT RESONATES



"Procera was born to deliver the first individualized CAD/CAM crown, following the Brånemark spirit of bringing biocompatibility to the restoration level."

Dr. Matts Andersson
Founder of Procera

We built the foundation.
You carry it forward.

Discover how Procera was created

A PARTNER FOR

GROWTH

Capture demand. Scale smart. Focus on what matters most, with minimal investment.

Scale across more indications

and deliver Procera quality across implant systems, strengthening your competitive positioning.

Benefit from flexible

open, and validated workflows tailored to your digital footprint.

Invest in quality, not equipment

by relying on authentic Procera solutions in a worry-free workflow that lets you focus on what you do best.



"In my lab we want the most esthetic workflows possible. The Procera ASC abutment allows us that flexibility to correct angulations and have non-cement retained restorations that adhere to the emergence profile predictably and consistently giving the dentists and patients the best possible esthetic and soft tissue outcomes."

Rob Dulieu (Dental technician)
RWD Dental Image - Australia

PROCERA SOLUTIONS



Procera implant crowns & abutments

On Nobel Biocare and non-Nobel Biocare implants



Procera implant bridges

Zirconia or titanium implant bridges on Nobel Biocare and non-Nobel Biocare implants



Procera bars

Fixed and fixed-removable implant bars on Nobel Biocare and non-Nobel Biocare implants



Open and validated workflows

Access to precision-manufactured Procera restorations through validated workflows.

What's available for me


Procera™
100%
Inspection for every restoration through industrial quality control, verifying dimensional tolerances as low as ±5 microns.
10
year Procera™ product warranty
LIFE-TIME
warranty on Nobel Biocare implants
What if predictability became something you chose, not something you chased?

A PARTNER FOR

EXCELLENCE

Excellence through predictable quality, esthetic mastery, and Procera expertise.

Standardize for predictability

using industrialized quality control to ensure consistent passive fit, helping reduce rework and inefficiencies.

Express your expertise

relying on carefully selected materials and cement-free ASC solutions that bring esthetic excellence to life.

Build on Procera's expert support

and proven processes to deliver consistent, high-quality outcomes with confidence.





"Procera is the only brand that offers angulated screw-channel zirconia bridges, cement-free, with a warranty to back it up. A warranty is only possible when you are fully confident in the product you deliver."

Luc and Patrick Rutten (Dental technicians)
Dental Team – Belgium

A PARTNER FOR

TRUST

Proven science. Proven quality.
Trusted long-term partner.

Build credibility

through science-backed solutions rooted in the legacy of the CAD/CAM pioneer, driving confident adoption of higher-value cases.

Protect your reputation

with proven workflows and clinically validated Procera solutions, helping you to be seen as a long-term partner.

Stay connected

through education and commercial programs that elevate and expand your professional network.




"Over the years, I've learned that trust is built through consistency. When clinicians ask why my cases are reliable long term, the answer is simple: they're based on clinically validated Procera solutions. That scientific foundation allows me to stand behind my work with confidence and be seen as a long-term partner, not just a lab."


Diogo Firmino (Dental technician)
Diogofirminodental - Spain



Connecting expertise and opportunities through education programs

The Procera difference





ONE CASE.
TWO OUTCOMES.



"Local milling can result in rough surfaces, limited precision, and inconsistent detail. Even with optimized protocols, milling and sintering may cause dimensional deformation, leading to clinical misfit."

"In full-arch restorations on Multi-unit Abutments, even minimal vertical discrepancies at the framework–Multi-unit Abutment interface can create micro gaps, that lead to plaque retention and risks of peri-implant soft-tissue inflammation."

Dr. Francesco Mintrone
Modena - Italy

***Imagine** what Procera could change for you.*

What's available for me

Procera, where precision meets beauty







A clinician's view

Dr. Édouard Cristofari
Paris - France

"In my daily practice, particularly in full-arch restorations, I rely on solutions that are reliable and reproducible. The precision of Procera milling ensures a consistently passive fit, which is critical for long-term mechanical stability and clinical confidence.

Beyond mechanical performance, the NobelProcera® Esthetic Zirconia material also enables restorations that are truly biomimetic, achieving natural, lifelike esthetics further enhanced by the ASC solution.

Finally, the direct screw-retained Multi-unit interface, without cement, is a key advantage, reducing biological risks associated with excess cement or Ti-base debonding".



A technician's view

Luis Saraiva
Clarity Lab - Portugal

"Our role is to put our passion and techniques, on our work to deliver consistency, esthetic, and long-term performance.

High-quality, highly reliable milling and quality-control help me to reduce adjustments, standardize our workflow, and stay focused on what I do best.

This is why Procera has become a trusted part of my routine, with the right protocol between dentist and dental technician I can rely on solid foundations to fully express my personal touch."

Procera™



Ready for
ceramic layering



Initial build up



Second build up



Final esthetic result

Case courtesy of
Dr. Édouard Cristofari
and Luis Saraiva

PROCERA RESTORATIONS

Customized implant crowns and abutments

Material choice: Zirconia and titanium with a true patented angulated screw channel, 100% cement-free in zirconia

Design flexibility with full-contour, partial cut-back, and customized options

Compatible with Nobel Biocare and non-Nobel Biocare implant interfaces.



Procera zirconia abutments are also available as full-contour or partial cut-back for veneering

10 year Procera™ product warranty

LIFE-TIME warranty on Nobel Biocare implants

Case courtesy of Dr. Giacomo Fabbri

More information

More science

Product options







100% cement free

Abutment

		Titanium abutment	Titanium ASC abutment	Zirconia abutment	Zirconia abutment	Zirconia abutment
Features			Angulated screw channel 0–25° Core-file			Angulated screw channel 0–25° 100% cement free
Material options		Titanium grade 5	Titanium grade 5 Anodized titanium	Legacy zirconia		
Platform	Conical connection	3.0, NP, RP, WP	NP, RP, WP	–	–	NP, RP, WP
	Tri-channel	NP, RP, WP, 6.0	NP, RP, WP, 6.0	–	NP, RP, WP, 6.0	–
	External hex	NP, RP, WP	–	NP, RP, WP	–	–
	Nobel Biocare N1™ Base*	–	–	–	–	–




100% cement free

100% cement free

Implant crown

		Implant crown	Implant crown
Features		Full-contour or cut-back zirconia Angulated screw channel 0–25° 100% cement free	Full-contour or cut-back zirconia Angulated screw channel 0–25° 100% cement free
Material options		Full-contour zirconia	Full-contour zirconia
Platform	Conical connection	NP, RP, WP	–
	Tri-channel	–	–
	External hex	–	–
	Nobel Biocare N1™ Base*	–	NP, RP

* Nobel Biocare N1 Base sold separately

All Procera restorations are co-packed with clinical screws.

Also available for non-Nobel Biocare implant systems:

Implant system	Implant connection
3i®	Certain® internal hex
Astra Tech Implant System™	OsseoSpeed®
	EV
HiOssen	ET
Neodent®	Grand Morse
Straumann®	Bone Level®
	Standard (SynOcta®)
	BLX
Zimmer	Screw-Vent®




PROCERA IMPLANT BRIDGES

NobelProcera® Esthetic Zirconia Implant Bridge



Multi-layered zirconia material available in 10 VITA shades.

Angulated screw channel on Multi-unit Abutment and conical connection platforms.

Design flexibility with full-contour, partial cut-back or full cut-back for FP1 to FP3 indications

10 year Procera™ product warranty

LIFE-TIME warranty on Nobel Biocare implants

Case courtesy of Michiel Wouters

More information

More science

Product options





		Zirconia implant bridge	Titanium implant bridge
Features		Angulated screw channel 0–25° 100% cement free Design options: Full-contour, cut-back and/or framework Available for FP1, FP2, FP3	Design options: Framework Available for FP1, FP2, FP3
Material options		Esthetic zirconia	Titanium grade 2
Platform	Multi-unit Abutment	NP, RP, WP	NP, RP, WP
	Conical connection	NP, RP, WP	NP, RP, WP
	Tri-channel	NP, RP, WP, 6.0	NP, RP, WP, 6.0
	External hex	NP, RP, WP	NP, RP, WP

All Procera restorations are co-packed with clinical screws.

Also available for non-Nobel Biocare implant systems on Multi-unit Abutment level:

Implant system	Implant connection
3i®	Osseotite®
Astra Tech Implant System™	OsseoSpeed®
	EV
BioHorizons®	Internal
Dentium™	Superline™
Neodent®	Grand Morse™
Osstem®	TS
Straumann®	Bone Level
	BLX
Zimmer®	Screw-Vent®





PROCERA BARS

Fixed and fixed-removable implant bars



Comprehensive range of fixed and fixed-removable implant bar solutions

Compatible with Nobel Biocare and non-Nobel Biocare interfaces

10 year Procera™ product warranty

LIFE-TIME warranty on Nobel Biocare implants

More information

More science

Product options



Titanium fixed implant bars



Titanium fixed-removable implant bars



Now available in exocad

Split bar with exocad™


Material options		Titanium	Titanium	Titanium grade 5
Platform	Multi-unit Abutment	NP, RP, WP	NP, RP, WP	NP, RP, WP
	Conical connection	NP, RP, WP	NP, RP, WP	NP, RP, WP
	Tri-channel	NP, RP, WP, 6.0	NP, RP, WP, 6.0	–
	External hex	NP, RP, WP	NP, RP, WP	–

All Procera restorations are co-packed with clinical screws.

Also available for non-Nobel Biocare implant systems on Multi-unit Abutment level:

Implant system	Implant connection
3i®	Osseotite®
Astra Tech Implant System™	OsseoSpeed®
	EV
BioHorizons®	Internal
Dentium™	Superline™
Neodent®	Grand Morse™
Osstem®	TS
Straumann®	Bone Level
	BLX
Zimmer®	Screw-Vent®





PROCERA CROWNS AND BRIDGES

Cemented crowns and bridges

Full-contour zirconia available in 10 VITA shades.

Full-contour design reduces the risk of veneer chipping, remakes and stress.

Fully shaded zirconia maintains its color even after adjustments.

Case courtesy of
Dr. Francesco Mintrone

10 year Procera™ product warranty

More information

Product options

Now available in

esthetic zirconia



	Zirconia crown	Zirconia bridge
Indication	1 unit	2–5 units
Material options	Full-contour esthetic zirconia	Full-contour esthetic zirconia



Procera centralized manufacturing facility, Mahwah, New Jersey

Design your way

Stay in the CAD environment you know and trust. No software disruption – simply continue using your DTX Studio™ suite, exocad™ or 3Shape® software and scanners.


Produce with confidence

Send your design directly to Procera for precision manufacturing, and benefit from high quality restorations.

Stay open, stay simple

Keep full flexibility with your preferred design tools while accessing industrial-grade Procera quality.

Paths to Procera



You digitize



You design

DTX STUDIO™ exocad 3shape

Scan and Design Services by Procera™ COMPLETE



We produce



You finalize

Procera digital solutions

Download updated libraries

How to videos

Procera™

Your trusted partner for growth and excellence

More information



NOBELBIOCARE.COM/PROCERA-SOLUTIONS

96626 GB 2606 Printed in the EU © Nobel Biocare Services AG, 2026. All rights reserved. Disclaimer: Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily at scale. They are for illustration purposes only and may not be an exact representation of the product. Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a dentist or a physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.

Part of the **Nobel Biocare** family